



Aug 19, 2020

Direct ELISA for investigating the binding of recombinant or chemically-made Protein-LG to immunoglobulins.

DOI

dx.doi.org/10.17504/protocols.io.bjxvkpn6

Angel A Justiz-Vaillant¹

¹University of the West Indies St. Augustine

University of the West In...

angel.vaillant@sta.uwi.e...



Angel A Justiz-Vaillant

University of the West Indies St. Augustine

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bjxvkpn6>

Protocol Citation: Angel A Justiz-Vaillant 2020. Direct ELISA for investigating the binding of recombinant or chemically-made Protein-LG to immunoglobulins.. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bjxvkpn6>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: August 19, 2020

Last Modified: August 19, 2020

Protocol Integer ID: 40661

Keywords: mammalian immunoglobulin, fc regions of many mammalian immunoglobulin, immunoglobulin, many mammalian immunoglobulin, protein Ig, reactivity with immunoglobulin, binding bacterial protein, linked immunosorbent assay, bacterial protein, binding protein, immunosorbent assay, protein, direct elisa, splg, made protein, enzyme

Abstract

Protein LG (SpLG) is an immunoglobulin-binding protein that interacts with the Fab and Fc regions of many mammalian immunoglobulins [1].

References

1. Vaillant AJ, McFarlane-Anderson N, Wisdom B, Mohammed W, Vuma S, et al. (2013) Immunoglobulin-binding Bacterial Proteins (IBP) Conjugates and their Reactivity with Immunoglobulin in Enzyme-Linked Immunosorbent Assays (ELISA). J Anal Bioanal Tech 4: 175. doi:10.4172/2155-9872.1000175

Materials

MATERIALS



Nunc®; 96-Well Polystyrene Round Bottom Microwell Plates, V 96 well plate, Non-Treated, clear, without lid, Sterile **Thermo Fisher Catalog #260210**



Protein-L from P. Magnus



Streptococcal protein G by Sigma Aldrich

Troubleshooting

- 1 This ELISA is used to study the interaction of protein-LG (SpLG) with diverse immunoglobulins. or chemically
- 2 The 96 well microtitre plate is coated overnight at 4°C with 1 µg/µl per well of purified immunoglobulins or 50 µl of any animal sera in carbonate-bicarbonate buffer pH 9.6.
- 3 Then plate is treated with bovine serum albumin solution and washed 4X with PBS-Tween.
- 4 Then 50 µl of peroxidase-labeled-protein-LG conjugate diluted 1:3000 in PBS-non-fat milk is added to each well and incubated for 1.30h at RT. After that the plate is washed 4X with PBS-Tween.
- 5 Pipette 50 µl of 3,3',5,5' - tetramethylbenzidine (TMB; Sigma-Aldrich) to each well.
- 6 The reaction is stopped with 50 µl of 3M H₂SO₄ solution.
- 7 The plate is visually assessed for the development of colour and read in a microplate reader at 450 nm.
- 8 A cut-off point should be calculated as the mean of the optical density of negative controls x 2 SD.