


Mar 20, 2020

Version 1

# 🌐 Direct nuclear tagmentation and RNA-sequencing (DNTR-seq) V.1

DOI

[dx.doi.org/10.17504/protocols.io.bdthi6j6](https://dx.doi.org/10.17504/protocols.io.bdthi6j6)



Vasilios Zachariadis<sup>1</sup>, Huaitao Cheng<sup>1</sup>, Nathanael Andrews<sup>1</sup>, Martin Enge<sup>1</sup>

<sup>1</sup>Karolinska Institutet

Enge lab



Vasilios Zachariadis

Karolinska Institutet

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bdthi6j6>

External link: <https://www.biorxiv.org/content/10.1101/2020.03.04.976530v1.full>

**Protocol Citation:** Vasilios Zachariadis, Huaitao Cheng, Nathanael Andrews, Martin Enge 2020. Direct nuclear tagmentation and RNA-sequencing (DNTR-seq). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bdthi6j6>

**Manuscript citation:**

A highly scalable method for joint whole genome sequencing and gene expression profiling of single cells bioRxiv  
2020.03.04.976530; doi: <https://doi.org/10.1101/2020.03.04.976530>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**Created:** March 17, 2020

**Last Modified:** March 20, 2020

**Protocol Integer ID:** 34377

**Keywords:** single-cell, scRNA-seq, scWGS, direct nuclear tagmentation, whole genome, gene expression, rna, complex cell mixture, single cell, mrna, gene, question in complex cell mixture, central challenge in biology, genetic variation, genotype, cell, dntr,

## Abstract

Understanding how genetic variation alters gene expression - how genotype affects phenotype - is a central challenge in biology. To address this question in complex cell mixtures, we developed Direct Nuclear Tagmentation and RNA-sequencing (DNTR-seq), which enables whole genome and mRNA sequencing jointly in single cells.

## Guidelines

Oligonucleotides (all ordered from IDT using Standard desalting, except barcodes ordered in solution/plates)

Oligo-dT: AAGCAGTGGTATCAACGCAGAGTACTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTT(N1:34333300)  
(N2:25252525)

IS\_PCR: 5'-AAGCAGTGGTATCAACGCAGAGT-3'

TSO: 5'-AAGCAGTGGTATCAACGCAGAGTACATrGrG+G-3'

ME-A: 5'-TCGTCGGCAGCGTCAGATGTGTATAAGAGACAG-3'

ME-B: 5'-GTCTCGTGGGCTCGGAGATGTGTATAAGAGACAG-3'

ME-Rev: 5'-/5Phos/CTGTCTCTTATACACATCT-3'

Illumina-compatible barcodes used (Sxxx/Nxxx series, n=784) are available as a supplementary table in the manuscript.



## Materials

### MATERIALS

- HotStart ReadyMix (KAPA HiFi PCR kit) **Kapa Biosystems Catalog #KK2601**
- Proteinase K **Thermo Fisher Scientific Catalog #EO0491**
- Tween-20 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P-7949**
- psfTn5 **addgene Catalog #79107**
- 10% SDS solution **Teknova Catalog #S0287**
- SMARTScribe Reverse Transcriptase **Takara Bio Inc. Catalog #634888**
- Magnesium chloride solution for molecular biology (1.00 M) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M1028**
- Ice
- Triton X-100 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #93426**
- Microseal® 'F' Foil **Bio-Rad Laboratories Catalog #MSF-1001**
- dNTP Mix (dATP, dCTP, dGTP, and dTTP, each at 10mM) **Thermo Fisher Scientific Catalog #R0192**
- KAPA HiFi PCR kit with dNTPs **Fisher Scientific Catalog #NC0142652**
- Betaine 5M **Merck MilliporeSigma (Sigma-Aldrich) Catalog #B0300**
- Dry ice
- UltraPure™ DNase/RNase-Free Distilled Water **Thermo Fisher Catalog #10977035**
- ERCC RNA Spike-In Mix **Thermo Fisher Catalog #4456740**
- USB Dithiothreitol (DTT), 0.1M Solution **Thermo Fisher Catalog #707265ML**
- Sera-Mag Speed Beads **GE Healthcare Catalog #65152105050250**
- RNase Inhibitor **Takara Bio Inc. Catalog #2313A**
- Hard-Shell® 384-Well PCR Plates thin wall skirted **Bio-Rad Laboratories Catalog #HSP3801**

### Before start

Bleach clean environment - to avoid DNA contamination. And RNase away or similar to avoid degraded RNAs. Prepare solutions in a strictly pre-PCR environment. Keep plates and reagents on ice unless otherwise noted.



## Prepare lysis buffer plates for cell sorting

### 1 Prepare **lysis buffer mix**

NOTE: Reagents are prepared on ice, working quickly. ERCC is stored in single-use aliquots at -80 °C , thawed on ice and added last.

| Reagent                        | Reaction conc. | μL per reaction | 384w plate |
|--------------------------------|----------------|-----------------|------------|
| Nuclease free H <sub>2</sub> O | -              | 1.965           | 786        |
| RNase Inhibitor (40u/μL)       | 1 unit/μL      | 0.075           | 30         |
| ERCC (1:1 200 000)             | -              | 0.075           | 30         |
| Triton-X100 (10% solution)     | 0.2%           | 0.06            | 24         |
| dNTP (10mM each)               | 2.5mM/each     | 0.75            | 300        |
| Oligo-dT (100μM)               | 2.5μM          | 0.075           | 30         |
|                                |                |                 |            |
| To dispense                    |                | 3               | 1200       |

Add 3 μL lysis buffer mix to each well. Cover with appropriate lids. Spin down.  
Snap freeze on **dry ice**. Store until use at -80 °C

## Sort single-cells

### 2 Sort single cells into 3 μL lysis buffer mix

Immediately seal with appropriate seals (approved for -80C > 100C) and centrifuge at 2000 x g, 4°C, 00:05:00

Snap freeze on **dry ice**. Store until use at -80 °C

## Separation of nuclear and cytosolic fractions

### 3 Thaw plate on ice.

Centrifuge at 500 x g, 4°C, 00:05:00 .

Keep on ice.

### 4 Transfer 2 μL from each well of the sorted plate into an empty 384-well plate. Use a low flow rate (**2mm/s**) and an aspiration height of **0.9mm** above the bottom.

**Note**

**NOTE:** We use the Eppendorf EpMotion 5073m benchtop liquid handler. We have successfully used other solutions, including the Hamilton STARlet, a semi-manual Gilson Platemaster 96-well pipette, and even manual 8-channel pipettes.

- 5 Spin down and freeze nuclear fraction at -20 °C to aid complete lysis.

If proceeding with **cDNA protocol** → step 12.

If proceeding with **DNA protocol** (step 6): spin down and snap freeze cytosolic fraction on **dry ice** and store at -80 °C

**Note**

**NOTE:** We will typically proceed with cDNA synthesis, unless experimental design dictates otherwise, to avoid an additional freeze-thaw cycle for mRNAs in the cytosolic fraction.

## Single-cell genomic libraries

- 6 Using plate with nuclear fraction, with remaining volume 1µl/well.

**Proteinase K treatment**

1. Dilute Proteinase K (stock 20mg/ml) to 0.2mg/ml by 30mM Tris-HCl pH8.0
2. Add 2 µL diluted Proteinase K (0.2mg/ml) to each well. Makes 0.13mg/ml reaction concentration.
3. Incubate in thermocycler at:
  - 50 °C 01:00:00
  - 80 °C 00:30:00
  - 4 °C hold

7 **Tn5 digestion**

Tn5 is produced from psfTn5 (Addgene #79107), purified to ~3mg/ml and assembled with Illumina Tn5 adapters (see **oligos**) as in *Picelli et al, 2014*.



### Citation

Picelli S, Björklund AK, Reinius B, Sagasser S, Winberg G, Sandberg R (2014)  
. Tn5 transposase and tagmentation procedures for massively scaled sequencing projects..  
Genome research.

<https://doi.org/10.1101/gr.177881.114>

LINK

- 7.1 Prepare 2X Tn5 Buffer. Keep assembled Tn5 enzyme (Picelli *et al*, 2014) on ice block and add last.

| Reagent  | Reaction conc                                     | μL per reaction | 384 w plate |
|--|---|-----------------|-------------|
| 5X TAPS-PEG (50mM TAPS, 25mM MgCl <sub>2</sub> , 40% PEG-8000) | 10mM TAPS<br>5mM MgCl <sub>2</sub><br>8% PEG-8000 | 1.6             | 672         |
| psfTn5, loaded with 50μM MED S-A/B                             |   | 0.1             | 42          |
| Nuclease free H <sub>2</sub> O                                 |   | 3.3             | 1386        |
|  |   |                 |             |
| To disperse  |   | 5               | 2100        |

Add  5 µL per well. Vortex and spin down plate.

#### Note

**NOTE:** Buffer contains PEG, which is viscous. 5X TAPS-PEG buffer should be allowed to assume room temperature before dispensing to allow proper mixing.

7.2 Incubate in thermocycler:  55 °C  00:10:00

Remove immediately and stop reaction by adding  2 µL per well of 0.2% SDS.


Vortex, spin down and incubate  00:10:00 at  55 °C

## 8 PCR amplification and barcoding

1. Prepare PCR master-mix

| Reagent                        | Reaction conc. | µl per reaction | 384 w plate |
|--------------------------------|----------------|-----------------|-------------|
| Nuclease free H <sub>2</sub> O | -              | 3.2             | 1280        |
| KAPA HiFi Buffer (5X)          | 1X             | 4               | 1600        |
| dNTP (10mM/each)               | 0.3mM/each     | 0.6             | 240         |
| KAPA enzyme (1u/µl)            | 0.02u/µl       | 0.4             | 160         |
| Tween-20 (10%)                 | 0.1%           | 0.2             | 80          |
|                                |                |                 |             |
| To dispense                    |                | 8.4             | 3360        |


2. Dispense  8.4 µL per well

3. Add primers/barcodes  1.6 µL per well (from 384-well index plates, with 3.75µM/each forward/reverse primers; see **oligos**). Total reaction volume is now 20µl (10µl sample + 10µl PCR mix and primers).

4. Vortex plate, spin down and incubate in thermocycler with the following program:

| Step     | Temperature | Time  | Cycles |
|----------|-------------|-------|--------|
| Gap fill | 72°C        | 3 min | 1x     |

|                 |        |        |     |
|-----------------|--------|--------|-----|
| First denature  | 95°C   | 30 sec | 1x  |
| Denature        | 95°C   | 15 sec | 18x |
| Anneal          | 67°C   | 30 sec |     |
| Extend          | 72°C   | 45 sec |     |
| Final extension | 72°C   | 4 min  | 1x  |
|                 | 4-10°C | hold   |     |

9 Pool  1.5 µL from each well into a 1.5mL Eppendorf tube.

## 10 Library cleanup

We prepare SPRI-beads in 20% PEG-8000 solution as in:

[https://openwetware.org/wiki/SPRI\\_bead\\_mix#Ingredients\\_for\\_50\\_mL\\_2](https://openwetware.org/wiki/SPRI_bead_mix#Ingredients_for_50_mL_2)

1. Add 0.9X SPRI-beads in 20% PEG solution. Incubate for  00:05:00

 Room temperature

2. Place on magnetic rack  00:03:00

3. Remove supernatant


4. Add 1 volume 80% EtOH (fresh). Incubate for  00:00:30

5. Remove supernatant

6. Repeat EtOH wash

7. Air dry for  00:10:00 -  00:15:00


8. Re-suspend beads thoroughly in  100 µL EB or TE buffer

9. Repeat cleanup (from step 1-7) and elute in  30 µL EB or TE buffer

## 11 (optional) Quality control of DNA libraries

Using Agilent HS 5000 DNA chips (or equivalent)





Pooled (and dliuted) DNA-library from 384-well plate.

This library was sequenced on a NextSeq loading 2.5pM based on a peak of 334bp. Sequencing was paired-end 37bp, 8bp dual index.

## Reverse transcription and cDNA amplification

12 Following step 4, cytosolic/RNA fraction plate contains 2µl solution per well.

### Primer annealing

Thaw plate. Spin down. Incubate in thermocycler at 72 °C for 00:03:00 .


Remove to ice immediately.

13 **Prepare RT master-mix**

| Reagent                  | React ion conc. | µl per reacti on | 384 w plate |
|--------------------------|-----------------|------------------|-------------|
| SMARTScribe RT (100u/µl) | 15u/µl          | 0.475            | 199.5       |
| RNase Inhibitor (40u/µl) | 1.66u/ µl       | 0.125            | 52.5        |
| 5X First Strand buffer   | 1X              | 1                | 420         |
| DTT (100mM)              | 8.33 mM         | 0.25             | 105         |
| Betaine (5M)             | 1.66M           | 1                | 420         |
| MgCl2 (1M)               | 10mM            | 0.03             | 12.6        |
| TSO (100uM)              | 1.66µ M         | 0.05             | 21          |



|                                |   |      |      |
|--------------------------------|---|------|------|
| Nuclease free H <sub>2</sub> O | - | 0.07 | 29.4 |
|                                |   |      |      |
| Total                          |   | 3    | 1260 |

Dispense  2 µL per well

Cover plate with new film and spin down.

#### 14 Incubate in thermocycler

 42 °C  01:30:00

 70 °C  00:05:00

 4 °C hold

#### 15 cDNA preamplification

|                                  | Reaction conc. | µl per reaction | 384 w plate |
|----------------------------------|----------------|-----------------|-------------|
| Nuclease free H <sub>2</sub> O   | -              | 1.0688          | 470.25      |
| Kapa HiFi HotStart ReadyMix (2X) | 1X             | 6.25            | 2750        |
| IS_PCR primer (10µM)             | 0.1µM          | 0.125           | 55          |
| Lambda Exonuclease (10u/µl)      | 0,045u/µl      | 0.05625         | 24.75       |
|                                  |                |                 |             |
| Total                            |                | 7.5000          | 3300        |

Dispense  7.5 µL per well . Total reaction volume will be 12.5µl.

#### 16 Spin down. Cover with new lid. Incubate in thermocycler with the following program:

| Step                 | Temperature | Time   | Cycles |
|----------------------|-------------|--------|--------|
| Lambda exonuclease   | 37°C        | 30 min | 1x     |
| Initial denaturation | 95°C        | 3 min  | 1x     |
| Denaturation         | 98°C        | 20 sec | 18-24x |












|                  |      |        |  |
|------------------|------|--------|--|
| Annealing        | 67°C | 15 sec |  |
| Elongation       | 72°C | 4 min  |  |
| Final elongation | 72°C | 5 min  |  |
|                  | 4C   | Hold   |  |

#### Note

**NOTE:** The number of cycles of pre-amplification will be different for different cell types. We suggest running a pilot (ideally qPCR-monitored to determine inflection point, for example by using 1X dsGreen to the reaction above)

## 17 cDNA cleanup

Using 20% SPRI-bead solution (as in step10 for DNA library cleanup).

1. Add 0.7X volume of SPRI beads per well. Mix well by pipetting
2. Incubate  00:05:00  Room temperature
3. Place on magnetic stand for  00:03:00
4. Carefully remove supernatant
5. Add  40 µL 80% EtOH and incubate  00:00:30
6. Remove EtOH (without disturbing the beads)
7. Wash again with EtOH. Make sure to remove well.
8. Allow beads to air-dry for  00:10:00 -  00:15:00
9. Remove plate from magnetic stand
10. Elute beads in  15 µL EB or TE buffer Mix well by pipetting
11. Incubate  00:05:00  Room temperature
12. Place on magnetic plate for  00:03:00
13. *Optional: Carefully remove supernatant to the elution plate*


## 18 cDNA quantification

We measure concentration of random wells using Qubit HS dsDNA, adapted to a 96-well plate reader.

1. Add 98.5  $\mu\text{L}$  of 1X Qubit HS dsDNA solution (or mix dye and buffer separately) to a flat-bottom, black plate
2. Add 1.5  $\mu\text{L}$  of cDNA sample
3. Add Standards (NOTE: We make a 8-step ladder from 0ng/ $\mu\text{L}$   $\rightarrow$  10ng/ $\mu\text{L}$  Qubit Standard DNA in TE buffer)
3. Read in plate reader using 485nm excitation/528nm emission
4. Calculate cDNA concentration

## 19 (optional) cDNA quality control

Using Agilent HS 5000 DNA chips (or equivalent)



Example of a single immune (=small) cell cDNA profile (cytosolic fraction from DNTR protocol)

## 20 Make cDNA dilution plate

Dilute cDNA based on average concentration from Qubit measurements.

Target concentration 150 pg per  $\mu\text{L}$  in 15  $\mu\text{L}$  (optionally in same plate)

## cDNA fragmentation


### 21 Prepare Tn5 master mix

Let TAPS-PEG equilibrate at 37°C and mix well before use.

| Reagent | Reaction conc. | $\mu\text{L}$ per reaction | 384 w plate |
|---------|----------------|----------------------------|-------------|
|         |                |                            |             |





|   |   |       |     |
|---|---|-------|-----|
| Nuclease free H <sub>2</sub> O                              | -   | 1.050 | 525 |
| TAPS-PEG (50mM TAPS, 25mM MgCl <sub>2</sub> , 40% PEG-8000) | 10mM TAPS<br>5mM MgCl <sub>2</sub><br>8% PEG-8000 | 0.500 | 250 |
| psfTn5, loaded with 50μM MED S-A/B                          |   | 0.250 | 125 |
|   |   |       |     |
| Total   |   | 1.800 | 900 |

Dispense  1.8 μL per well in a new plate (**tagmentation plate**)

22 Add  0.7 μL cDNA (normalized to **150pg/μl**)

Mix well by vortexing plate. Cover with new lid and spin down.

23 Incubate in thermocycler at  55 °C  00:10:00

Remove immediately and stop reaction by adding  1 μL per well of 0.2% SDS.

Vortex, spin down and incubate  00:07:00 at  55 °C

## cDNA library PCR and barcoding


24 **Make PCR master-mix**

| Reagent          | Reaction conc. | μl per reaction | 384 w plate (420 x) |
|------------------|----------------|-----------------|---------------------|
| H <sub>2</sub> O | -              | 13.25           | 5565                |
| KAPA HiFi        | 1X             | 5               | 2100                |



|                           |                |      |     |
|---------------------------|----------------|------|-----|
| Buffer (5X)               |                |      |     |
| dNTP (10mM/each)          | 0.3mM/each     | 0.75 | 315 |
| KAPA enzyme (1u/ $\mu$ l) | 0.02u/ $\mu$ l | 0.5  | 210 |
|                           |                |      |     |
| Total                     |                | 19.5 |     |

Dispense  19.5  $\mu$ L per well to **tagmentation plate** (containing 3.5 $\mu$ L sample after step 23)


- 25 Add primers/barcodes  2  $\mu$ L per well (from 384-well index plates, with 3.75 $\mu$ M/each forward/reverse primers; see **oligos**).

Total reaction volume is 25 $\mu$ L (3.5 $\mu$ L sample + 21.5 $\mu$ L PCR mix and primers).

- 26 Vortex. Spin down and cover. Incubate in thermocycler as below:

| Step            | Temperature | Time   | Cycles |
|-----------------|-------------|--------|--------|
| Gap fill        | 72°C        | 3 min  | 1x     |
| First denature  | 95°C        | 30 sec | 1x     |
| Denature        | 95°C        | 15 sec | 12x    |
| Anneal          | 67°C        | 30 sec |        |
| Extend          | 72°C        | 45 sec |        |
| Final extension | 72°C        | 4 min  | 1x     |
|                 | 4-10°C      | hold   |        |

## cDNA library pooling and clean-up

- 27 Pool  2.5  $\mu$ L from each well to an 1.5ml Eppendorf tube

- 28 **Library cleanup (as for DNA libraries)**

We prepare SPRI-beads in 20% PEG-8000 solution as in:

[https://openwetware.org/wiki/SPRI\\_bead\\_mix#Ingredients\\_for\\_50\\_mL\\_2](https://openwetware.org/wiki/SPRI_bead_mix#Ingredients_for_50_mL_2)

1. Add 0.9X SPRI-beads in 20% PEG solution. Incubate for 00:05:00

Room temperature

2. Place on magnetic rack 00:03:00

3. Remove supernatant

4. Add 1 volume 80% EtOH (fresh). Incubate for 00:00:30

5. Remove supernatant


6. Repeat EtOH wash

7. Air dry for 00:10:00 - 00:15:00

8. Re-suspend beads thoroughly in 100  $\mu$ L EB or TE buffer

9. Repeat cleanup (from step 1-7) and elute in 30  $\mu$ L EB or TE buffer

## 29 Pooled library QC



Pooled cDNA library of 784 cells on HS D5000 Agilent tapestation

## Citations

### Step 7

Picelli S, Björklund AK, Reinius B, Sagasser S, Winberg G, Sandberg R. Tn5 transposase and tagmentation procedures for massively scaled sequencing projects.

<https://doi.org/10.1101/gr.177881.114>