



Mar 08, 2024

🌐 Dispensing (non-L1) C. elegans to 96 well imaging plates using Integra VIAFILL



Forked from [Dispensing C. elegans to 96 well tracking plate using Integra VIAFILL](#)

DOI

dx.doi.org/10.17504/protocols.io.rm7vzxwj2gx1/v1

e.warren Warren¹, Bonnie Evans^{1,2}

¹MRC Laboratory of Medical Sciences; ²Imperial College London

Behavioural Genomics



Bonnie Evans

MRC Laboratory of Medical Sciences, Imperial College London

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzxwj2gx1/v1>

Protocol Citation: e.warren Warren, Bonnie Evans 2024. Dispensing (non-L1) C. elegans to 96 well imaging plates using Integra VIAFILL. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.rm7vzxwj2gx1/v1>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 08, 2024

Last Modified: March 08, 2024

Protocol Integer ID: 96374

Keywords: interga viafill dispenser, imaging plate, using interga viafill protocol, interga viafill protocol, agar plates in advance, well plate, plate, agar plate, bleach

Abstract

Protocol for dispensing later larval stage *C. elegans* into 96 well plates using the Interga VIAFILL dispenser. A bleach synchronised *C. elegans* population should be "re-fed" onto agar plates in advance.

Materials

M9 Buffer [3 g KH_2PO_4 , 6 g Na_2HPO_4 , 5 g NaCl, 1 ml 1 M MgSO_4 , H_2O to 1 litre. Sterilize by autoclaving.]

Integra viafill

Integra Viafill Cassette: 8 Channel Small Bore. Part number: 5722

Troubleshooting

Configure the Integra VIAFILL

- 1 Insert a small cassette into the machine.
- 2 Create a new program (or use existing program) with the following settings:
 - 2.1 Mode: Repeat Dispense
Tubing: 8 channel small
Plate: 96
Volume: 10 μ l
Map: Full
Speed: Fast
 - 2.2 To adjust the height, put a plate into the stage and press the "Height" button on the display. The plate clearance, dispense and pre-dispense heights can all be easily adjusted using the arrows. Click the back arrow to return to the program screen.
 - 2.3 Save these settings.

Preparing *C. elegans*

- 3 Wash worms off plates with several ml of M9 solution and transfer to 15ml falcon tubes.
- 4 Centrifuge for 2 minutes at 1500 rpm (RCF:210, ascending 9; descending 7) – program 1
- 5 Remove supernatant with plastic Pasteur pipette.
- 6 Add 15 ml M9.
- 7 Repeat steps 3 to 6.
- 8 Transfer the worm solution to a sterile conical flask.



- 9 Gently mix the worm suspension by swirling the flask, and pipette ten 10 μ l droplets to a microscope slide or cover glass. Count the number of worms present in each droplet.
- 10 Dilute the worm suspension to achieve your desired worm density (30 to 60 worms per 100 μ l) by adding M9 buffer.

Dispense *C. elegans* with the VIAFILL

- 11 Add a sterile stir bar in the conical flask, and place on a magnetic stirrer.
- 12 Turn the magnetic stirrer onto the lowest setting to agitate the worm suspension.
- 13 Select the VIAFILL program.
- 14 Select the volume you wish to dispense (10 μ l).
- 15 Place the end of the tubing from the cassette into the worm solution.
- 16 Press "Prime" to prime the tubing. Check all tips of the cassette are dispensing correctly.
- 17 **Test dispense** 96 droplets onto a plate lid.
- 18 Count the number of worms in 10 random droplets. Dilute worm solution with M9 to the correct desired density again if necessary.
- 19 Press "Recover" and re-prime the tubing before dispensing onto imaging plates
- 20 Place a 96-well imaging plate in the stage and press "Start".



- 21 Once all imaging plates have been filled, dry with the lids off under laminar flow for 1 to 2 hours.

Cleaning the cassette

- 22 Press "Recover" to recover the worm solution into the bottle.

Note

Dispose of excess worms by adding Rely+On™ *Virkon* to 1%.

- 23 Place the tubing ends into water.
- 24 "Prime" the tubing several times so that the water runs through.
- 25 Remove the end of the tubing from the water and "Prime" with air to remove any remaining liquid.
- 26 Remove the cassette from the machine, place in an empty tip box, and autoclave.