



Jun 19, 2018

Dissection and Fixation of Drosophila CNS

DOI

dx.doi.org/10.17504/protocols.io.qjydupw

Geoffrey W. Meissner¹, Jonathan B. Grimm¹, Rebecca M. Johnston¹, Ben Sutcliffe², Julian Ng², Gregory S.X.E. Jefferis², Sebastian Cachero², Luke D. Lavis¹, Oz Malkesman¹

¹Janelia Research Campus, Howard Hughes Medical Institute, Ashburn, Virginia, USA;

²Division of Neurobiology, Medical Research Council Laboratory of Molecular Biology, Cambridge CB2 0QH, United Kingdom



Geoffrey Meissner

HHMI Janelia Research Campus

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.qjydupw>

External link: <https://www.janelia.org/project-team/flylight/protocols>

Protocol Citation: Geoffrey W. Meissner, Jonathan B. Grimm, Rebecca M. Johnston, Ben Sutcliffe, Julian Ng, Gregory S.X.E. Jefferis, Sebastian Cachero, Luke D. Lavis, Oz Malkesman 2018. Dissection and Fixation of Drosophila CNS. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.qjydupw>



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 30, 2018

Last Modified: June 19, 2018

Protocol Integer ID: 12632

Keywords: chemical tagging of neuron, fixation of drosophila cn, drosophila cn, drosophila brain, improved drosophila melanogaster transgene, drosophila melanogaster transgene, rapid labeling of specific cell, chemical fluorophore ligand, labeling technique with standard immunohistochemistry, new fluorophore ligand, neuron, chemical tagging, neural tissue, cell, rapid labeling, labeling tag, standard immunohistochemistry, janelia fluor, specific cell


Abstract

The use of genetically encoded 'self-labeling tags' with chemical fluorophore ligands enables rapid labeling of specific cells in neural tissue. To improve the chemical tagging of neurons, we synthesized and evaluated new fluorophore ligands based on Cy, Janelia Fluor, Alexa Fluor, and ATTO dyes and tested these with recently improved *Drosophila melanogaster* transgenes. We found that tissue clearing and mounting in DPX substantially improves signal quality when combined with specific non-cyanine fluorophores. We compared and combined this labeling technique with standard immunohistochemistry in the Drosophila brain.

Troubleshooting

Summary

1



De-wax Anesthetized Flies

- 2 Prepare three wells in a glass well plate: one with 70% ethanol and the other two with cold Schneider's Insect Medium (S2). Place a sticky dot next to the 70% well to clearly distinguish it from the other wells filled with S2.
- 3 Anesthetize flies with CO₂ or cold. Grasp a fly by the wings or legs with forceps and briefly submerge the anesthetized fly first in 70% ethanol (2 seconds) followed immediately by a brief dip in the first well of cold S2. Then submerge the fly in the second well of cold S2 where it will remain until dissected.

Do not use a transfer pipette to move flies between wells because this adds ethanol to the S2 wells. Extended exposure to ethanol will accelerate denaturation of the tissue and fluorescent proteins making the tissue unusable.

- 4 Keep the dish with flies submerged in S2 on ice until they are dissected.

Do not anesthetize and rinse more flies than can be dissected within 20 minutes. Flies kept submerged longer will die and their brains will be unusable due to postmortem changes

Dissection

- 5 Dissect in cold S2. Transfer the rinsed fly to a Sylgard-lined dish with cold S2 and dissect.

Replace your puddle of S2 with fresh cold S2 when it becomes littered with dissection debris.

Dissect only as many flies as you can comfortably complete in 20 minutes before transferring the dissected brains to fixative. All fixations are precisely timed (see next).

Fixation

- 6 **Timed Fixation.** Within 20 minutes of dissection, transfer tissue to a 2mL Protein LoBind tube with ~1.9 mL of 2% paraformaldehyde (PFA) in S2 at room temperature (RT) and incubate for 55 minutes at RT while on a nutator. Cover tissue to protect from light.

Optimal fixation time is 55 minutes for 2% PFA in S2. Use a timer to precisely time the duration of your fixation.

- 7 **Fix Removal - Rinse.** Place the tubes upright to allow the tissue to sink. Use a transfer pipetted to aspirate the fix and fill tube with phosphate buffered saline with 5% Triton X-100 (PBT) at RT. Invert the tube a few times. Let the tissue settle and aspirate the rinse solution with a transfer pipette.

To help the tissue sink, add a small amount (~15 μ L) of PBT at RT to the tube. The detergent (Triton) in the PBT displaces air trapped in the trachea, helping the tissue sink.

- 8 **Fix removal - Washes.** Add 1.75 mL of RT PBT and nutate for 10-20 minutes at RT. Repeat for a total of 4 washes. Protect samples from light during washes.

- 9 **Storage.** Store the tubes of tissue in 0.5% PBT at 4°C. Nutate or lay the tubes flat in a covered box on a rotator. Do not store upright. Protect from light.

Typically, these tubes will begin the IHC (immunohistochemistry) process the following day but can be stored for up to 3 days. If stored more than overnight, aspirate the old PBT and do a brief wash with 0.5% PBT before beginning IHC processing.