

Feb 11, 2019

Dissection and fixation of murine colonic tissue for myenteric plexus visualization

DOI

dx.doi.org/10.17504/protocols.io.xz6fp9e

Dante Heredia¹, Terence Smith¹

¹University of Nevada Reno, School of Medicine



Dante Heredia

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.xz6fp9e

Protocol Citation: Dante Heredia, Terence Smith 2019. Dissection and fixation of murine colonic tissue for myenteric plexus visualization. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.xz6fp9e>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 11, 2019

Last Modified: February 11, 2019

Protocol Integer ID: 20254

Keywords: Myenteric plexus, colon, large intestine, enteric

Abstract

Protocol for harvest of colonic intestinal tissue, with the intent of imaging the myenteric plexus.



Materials

Zamboni fixative can be purchased commercially or 1.5% picric acid and 2% paraformaldehyde in 0.1N PBS.
Aquamount mounting media is commercially available

Safety warnings

! Fixatives can be toxic. Exercise proper use of PPE when handling.

Before start

Prepare Zamboni Fixative

Prepare Krebs-Ringers Solution



- 1 A ventral midline incision is made and the whole colon is carefully excised into a Sylguard lined dissection dish.
- 2 Cut along the mesenteric border until the colonic tube is now a rectangular in shape.
- 3 Pin the colon at 110% of length and width mucosa side up. Gently remove the mucosal layer. Re-pin in a new Sylguard lined dish at 100% ensuring tissue is taught but not stretched.
- 4 Fix tissue using ice cold Zamboni fixative for 15 minutes.
- 5 Remove tissue from dish and wash 3 times in PBS for 15 minutes each wash.
- 6 Cut tissue into 3 sections; oral, middle and anal (in order to fit on slide) and mount on glass slide using Aquamount mounting media and a glass cover slip.
- 7 Image sample.