



Jan 21, 2026

Dissociated Cardioids Assembly Using 30% Matrigel Droplets

DOI

<https://dx.doi.org/10.17504/protocols.io.3byl489e2vo5/v1>

Harry Akligoh¹

¹Northeastern University - Konry lab



Harry Akligoh

Northeastern University - Konry lab

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl489e2vo5/v1>

Protocol Citation: Harry Akligoh 2026. Dissociated Cardioids Assembly Using 30% Matrigel Droplets. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.3byl489e2vo5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: January 21, 2026



Last Modified: January 21, 2026

Protocol Integer ID: 239041

Keywords: dissociated cardioids assembly, cardiac organoid, microfluidic way

Abstract


Matured cardiomyocytes are needed for disease modelling and drug screening. In this protocol we describe a microfluidic way of reassembling dissociated cardiac organoids back into cardiac organoids.

30% Matrigel Droplet Generation and Dissociated Cardioid Assembly

- 1 Thaw frozen stock of dissociated cardioids in a waterbath set to 37C for two minutes. Don't allow it to thaw completely.
- 2 Transfer into 15ml Falcon tube and resuspend in 5ml RPMI 1640+B27 supplement.
- 3 Centrifuge at 200xg for 5 min and decant the supernatant.
- 4 Add 1ml RPMI 1640+B27 to the cell pellet and resuspend.
- 5 Measure cardiomyocytes viability by taking 10ul of cell suspension and adding it to 10ul of Trypan blue (1:1). Count using an automated countess 3 haemocytometer.
- 6 Prepare 30% Matrigel by diluting down stock Matrigel with RPMI 1640+B27.
- 7 Resuspend 2-5 million dissociated cardioids in 30% Matrigel and use for droplet generation.
- 8 Generate about 300um droplets off-chip into 1.5ml eppendorf tube prefilled with 100-500ul of RPMI 1640+B27.
- 9 Aliquot 2ml of RPMI 1640+B27 into each well of 6-well plate.
- 10 Pipette 100ul of 30% Matrigel-cardiomyocytes droplets into each well of the 6-well plate.
- 11 Swell the plate gently making sure not to merge the droplets.
- 12 Incubate in a 37C/5% CO2 incubator and change media everyday.

13 Take images and videos everyday beginning Day 0 until Day 10. On each day check for viability by staining with Calcein-AM (Add 1.5ul of Calcein-AM to 1.5ml of RPMI 1640+B27 and use that to stain the droplets for 2 hours)

14



Schema of encapsulated dissociated cardioids in 30% Matrigel