



Mar 02, 2024

## DNA Concentration Protocol

DOI

<https://dx.doi.org/10.17504/protocols.io.bp2l6xprklqe/v1>

chyanne.rosenbaum<sup>1</sup>

<sup>1</sup>Eagle Fish Genetics Lab



**EagleFish GeneticsLab**

Eagle Fish Genetics Lab

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.bp2l6xprklqe/v1>

**Protocol Citation:** chyanne.rosenbaum 2024. DNA Concentration Protocol. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bp2l6xprklqe/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** November 21, 2023



**Last Modified:** March 02, 2024

**Protocol Integer ID:** 91296

**Keywords:** dna concentration protocol, dna concentration protocol this protocol, eagle fish genetics lab, extracted dna, dna, protocol, times the initial concentration, concentration, initial concentration, fish, lab

## Abstract

This protocol is how the **Eagle Fish Genetics Lab** takes extracted DNA in a 96 well plate and concentrates it to five times the initial concentration.

## Guidelines

This creates a 5x concentration.



## Materials

### Section 1 - Prepare Your Plate

- a. Corresponding Initial DNA Extraction plate(s)
- b. low profile non-skirted PCR tray
- c. p200 pipette (green) set to 25  $\mu$ L
- d. p200 pipette (green) tips
- e. Nexttec sealing tape (clear plastic seal)

Machine(s): Vortex, centrifuge, and PCR plate heat sealer

### Section 2 - Thermocycler

Machine(s): Thermocycler

### Section 3 - Rehydrate

Supplies needed:

- a. Nuclease free water
- b. Reservoir for the water
- c. p10 multichannel pipette (red) set to 5  $\mu$ L
- d. p10 pipette tips (red)
- e. Heat seal

Machine(s): Vortex, centrifuge, and PCR plate heat sealer



## Section 1 - Prepare Your Plate

- 1 Supplies needed per concentration plate:
  - a. Corresponding Initial DNA Extraction plate(s)
  - b. low profile non-skirted PCR tray
  - c. p200 pipette (green) set to 25  $\mu$ L
  - d. p200 pipette (green) tips
  - e. Nexttec sealing tape (clear plastic seal)
- 2 Sign up for a thermocycler (1-2 hours)
- 3 Label low profile non-skirted PCR tray: "Concentration", date, tray #, project initials and your initials.
- 4 Using the p200, take 25  $\mu$ L of initial DNA extract and dispense into the label low profile non-skirted PCR tray.
  - a. Use fresh tips for each sample to prevent contamination.
  - b. *Pro tip* – if you are not going to fill the whole tray, outline the working wells with sharpie.
- 5 Seal the plate with clear plastic seal.
- 6 Balance and spin down at maximum speed (3700) for 5 minutes on centrifuge.

## Section 2 – Thermocycler

- 7 Turn on thermocycler by using the switch on the back right of the machine.
- 8 Open the lid to the thermocycler.
- 9 Remove clear plastic seal from your plate and discard the seal.
- 10 Place your plate into the machine.

- a. Keep the lid of the machine OPEN.
- 11 Follow the string of commands below which corresponds to your thermocycler.
- a. PCT-100 (Machine #15)
    - i. Select "Run" → Press "Proceed" → Select "<MAIN>" → Press "Proceed" → Select "70EVAPO" → Press "Proceed" → Select "NO" → Press "Proceed" → the machine will start.
  - b. DNA Engine (Machine #16)
    - i. Select "Run" → Press "Proceed" → Select "<MAIN>" → Press "Proceed" → Select "EVAPO-85" → Press "Proceed" → Select "PLATE" → Press "Proceed" → next to the "VOLUME" enter "25" → Press "Proceed" → the machine will start.
- 12 Check after half an hour to see if all liquids have been evaporated, if not, continue to monitor.
- a. Avoid continuous heating of dry wells, do not over evaporate.
- 13 Turn off thermocycler with switch in the back and the close lid.

### Section 3 - Rehydrate

- 14 Supplies needed:
- a. Nuclease free water
  - b. Reservoir for the water
  - c. p10 multichannel pipette (red) set to 5  $\mu$ L
  - d. p10 pipette tips (red)
  - e. Heat seal
- 15 Fill a reservoir a quarter of the way with nuclease free water.
- 16 Dispense 5  $\mu$ L of nuclease free water into each of the evaporated working wells.
- a. Use a new tip for each well.
  - b. Do NOT use the hover techniques as some DNA may be aerosolized and can still contaminate the hovering tips.
- 17 Seal plate with heat seal
- 18 Vortex WELL, quickly spin down and let sit at room temperature for 5 min.

