



May 28, 2025

DNA extraction and DNA quality control from fish tissue using the QIAwave DNA Blood & Tissue Kit (QIAGEN) or NucleoSpin Tissue Kit (MACHEREY-NAGEL)



DOI

<https://dx.doi.org/10.17504/protocols.io.yxmvmmyoz9v3p/v1>

Hervé Rogissart^{1,2,3}, Cecile Chardon¹

¹UMR CARRETEL; ²INRAE; ³USMB



Hervé Rogissart

USMB-INRAE

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yxmvmmyoz9v3p/v1>

Protocol Citation: Hervé Rogissart, Cecile Chardon 2025. DNA extraction and DNA quality control from fish tissue using the QIAwave DNA Blood & Tissue Kit (QIAGEN) or NucleoSpin Tissue Kit (MACHEREY-NAGEL). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.yxmvmmyoz9v3p/v1>

Manuscript citation:

<https://www.qiagen.com/us/resources/resourcedetail?id=d155700e-edce-4fb1-9f97-8faa3265ec8e&lang=en>

<https://www.mn-net.com/media/pdf/5b/d0/d9/Instruction-NucleoSpin-Tissue.pdf>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 20, 2025

Last Modified: May 28, 2025

Protocol Integer ID: 124736

Keywords: DNA extraction, Fish, Salmonids, dna quality control from fish tissue, dna extraction from other fish species, genomic dna from fish tissue, dna extraction procedure, summary sheet of the dna extraction procedure, dna extraction, qiawave dna blood, extracting genomic dna, fish tissue, dna quality control, dna quality, nucleospin tissue kit, other fish species, dna, fish, extraction, tissue kit, salvelinus alpinus

Abstract

This protocol describes a column-based method for extracting genomic DNA from fish tissue (e.g., *Salvelinus alpinus*) using either the QIAwave DNA Blood & Tissue Kit (QIAGEN) or the NucleoSpin Tissue kit (MACHEREY-NAGEL), adapted from the manufacturers' recommendations. The protocol is also well-suited for DNA extraction from other fish species. A summary sheet of the DNA extraction procedure is included to support its practical implementation in the laboratory.

Attachments



Rogissart-Chardon-

DN...

198KB

Image Attribution

Hervé Rogissart



Materials

■ Samples

- Fish tissue
- Preserved in ethanol (final concentration > 90 % of ethanol) or at  -20 °C
- Weight needed = 2-10 mg

■ Reagents

To clean the workspace

- DNA/RNA-ExitusPlus
- Ethanol (70 %)

For prepare sample:

- DNA/RNA-ExitusPlus
- Ethanol 96 - 100 %
- Water

For DNA extraction

- DNA extraction kit: **Qiagen Kit** or **Macherey-Nagel Kit**
- Ethanol 96 - 100 %, molecular grade to prepare wash buffer
- Ultrapure water

■ Materials (excluding solutions preparation)

- Bunsen burner
- Chisel and lab pincer
- Balance
- Specific DNA-work station (sterile area equipped with air filtration)
- Microcentrifuge for 1.5 to 2 mL tubes (relative centrifugal force needed: 6000 to 19000 x g)
- Horizontal vortexer with microtube holder
- Vortexer
- Stove or dry bath (to have 56 °C)
- Pipettes: 100-1000 µL; 10-100 µL
- 2 trash cans: 1 for liquid and 1 for solid

■ Consumables

- Tips with filter:
 - > 1000 µL
 - > 100 µL
- 2 mL sterile microcentrifuge tube: 2 per sample
 - 1 to prepare and lyse tissue (steps 1-2-3) and 1 to elute DNA (step 6)
- Gloves



Protocol materials

✕ Qubit™ dsDNA BR Assay Kit **Thermo Fisher Scientific Catalog #Q32853**

✕ NucleoSpin Tissue **Macherey-Nagel Catalog #740952.250**

✕ DNA/RNA-ExitusPlus **Applychem Catalog #A7089**

✕ QIAwave DNA Blood & Tissue Kit **Qiagen Catalog #69556**

Safety warnings

- ⚠ Buffers used in these kits require safety precautions, mainly Buffer AL and Buffer AW1 or Buffer B3 and BW which contain guanidine salt.

Safety informations of all buffers are available at :

Qiawave DNA Blood & Tissue Kit (QIAGEN)

NucleoSpin Tissue (MACHEREY-NAGEL)

Before start

- Wear gloves throughout the extraction process

1. Preheat dry bath or stove to  56 °C

2. Clean the workspace with  DNA/RNA-ExitusPlus **Applychem Catalog #A7089** followed by

 70 % volume **ethanol**

3. Unpack  QIAwave DNA Blood & Tissue Kit **Qiagen Catalog #69556** (**Qiagen Kit**)

or

 NucleoSpin Tissue **Macherey-Nagel Catalog #740952.250** (**Macherey-Nagel Kit**)

4. Check/Prepare buffers

5.1. If using **Qiagen Kit**

- Buffer ATL** and **Buffer AL** may precipitate during storage

- Warm to  56 °C until fully dissolved

- Before first use, prepare **Buffer AW1**, **Buffer AW2** and **Buffer AE** as indicated on the manufacturer's protocol. The information below is for a 250 columns kit.

→ **Buffer AW1**: In a  250 mL sterile glass bottle, transfer the entire volume of **Buffer AW1/C** ( 98 mL) and add **ethanol**  96-100 % volume ( 130 mL). Cap the glass bottle tightly and mix thoroughly by inverting the bottle several times. To label the glass bottle, use the enclosed label and transfer it onto the glass bottle. The final volume of **Buffer AW1** is  228 mL . This solution is stable for 1 year when stored closed at room temperature.

→ **Buffer AW2**: In a  250 mL sterile glass bottle, transfer the entire volume of **Buffer AW2/C** ( 6 mL) and add **ultrapure water** ( 60 mL) + **ethanol**  96-100 % volume ( 160 mL). Cap the glass bottle tightly and mix thoroughly by inverting the bottle several times. To label the glass bottle, use the enclosed label and transfer it onto the glass bottle. The final volume of Buffer AW2 is  226 mL . This solution is stable for 1 year when stored closed at room temperature.

→ **Buffer AE**: In a 150mL sterile glass bottle, transfer the entire volume of buffer **AE/C**  10 mL (10 mL) and add **ultrapure water**  110 mL . Cap the glass bottle tightly and mix thoroughly by inverting the bottle several times. To label the glass bottle, use the enclosed label and transfer it onto the glass bottle. The final volume of Buffer AE is  120 mL .

- Prepare aliquots of each solution as needed



5.2. If using **Macherey-Nagel Kit**

- **Buffer T1** and **Buffer B3** may precipitate during storage

- Warm between  50 °C and  70 °C until fully dissolved

- Before first use, prepare **Buffer B5** and **Proteinase K solution** as indicated on the manufacturer's protocol.

→ **Buffer B5**: add appropriate volume of **ethanol**  96-100 % volume to **Buffer B5 Concentrate** as indicated on the bottle and write the date on the bottle. This solution is stable for 1 year when stored closed at

 Room temperature .

→ **Proteinase K solution**: add appropriate volume of Proteinase Buffer PB (provided) to lyophilized **Proteinase K**. Ensure all powder is dissolved and solution is homogeneous. Transfert **Proteinase K solution** in 1.5 or 2 mL sterile microtube. This solution is stable for 6 months at  -20 °C .

- *Prepare aliquots of each solution as needed*



Prepare sample:

- 1
 - Before start
 - For each sample, label a  2 mL sterile microcentrifuge tube (Eppendorf safe-lock) with a unique identifier
 - Decontaminate cutting tools: clean with DNA ExitusPlus, rinse with water and cover with ethanol
 - Sterilize cutting tools using bunsen burner
 - Cut up to approximately **2-10 mg de tissue**
 - Place tissue in annotated tube

Note

- Decontaminate or replace cutting tools for each sample to prevent cross-contamination
- Cutting the tissue into small pieces to enable more efficient lysis
- Store at  +4 °C for quick use (< 8h) or at  -20 °C for further use (> 8h)

Lyse sample 1/2:

- 2
 - Prepare the lysis mixture containing the following per sample

	A	B
	Qiagen kit	Macherey-Nagel kit
	180 µL buffer ATL	180 µL buffer T1
	20 µL Proteinase K	25 µL Proteinase K

Lysis mixture composition

Note

- Homogenize **Buffer** and **Proteinase K**
- When preparing the mix for multiple samples, allow a marge for 1-2 additional samples

- Add lysis mixture to each sample



A	B
Qiagen kit	Macherey-Nagel kit
200 µL	205 µL

Lysis mixture volume / sample

- Vortex briefly  00:00:05 by pulse-vortexing
- Centrifuge , 00:00:05 in a microcentrifuge

Note

Check that tubes are tightly sealed to prevent evaporation during incubation, and ensure that tissue samples are fully immersed in the lysis buffer. If this was not the case, a brief centrifugation was performed to re-establish contact between the tissue and the solution.

- Incubate at  56 °C for at least  02:00:00 in a dry bath or in a stove

Note

- If using a thermomixer set to  300 rpm
- Make sure the incubator is preheated to 56°C otherwise the lysis time must be increased
-  Overnight incubation is applied to ensure that digestion will be complete and/or to facilitate the organization of extraction.

Lyse sample 2/2 & Adjust DNA binding conditions:

- 3
 - To save time while waiting for the **lysis incubation**, per sample prepare and annotate:
 - spin columns for DNA purification placed in waste tubes (provided)
 -  2 mL microtubes for DNA elution
 - Few minutes before the end of incubation, prepare **buffer and ethanol mixture** containing the following per sample



	A	B
	Qiagen kit	Macherey-Nagel kit
	200 µL AL	200 µL B3

Buffer and ethanol mixture composition

Note

- Homogenize **Buffer** and **ethanol mixture**
- When preparing the mix for multiple samples, allow a marge for 1-2 additional samples

- Each sample must be vortexed briefly  00:00:15 and centrifuged briefly  00:00:05
- Add **buffer and ethanol mixture**

	A	B
	Qiagen kit	Macherey-Nagel kit
	400 µL	410 µL

Buffer and ethanol mixture volume / sample

- Immediately vortex for  00:00:15 to mix and centrifuge  , 00:00:05 in a microcentrifuge

Note

- It is important to immediately mix the samples by vortexing or pipetting after adding the **Buffer and ethanol mixture** to ensure a homogeneous solution. Failure to do so may affect DNA yield.

Bind DNA:

- 4
 - Transfer the mixture samples to **spin column** placed in a waste tube (provided and annotated)



- Centrifuge at Room temperature to **capture the DNA** on the column membranes. Follow the parameters indicated in the table below.

A	B	C
	Qiagen kit	Macherey-Nagel kit
Column name	DNeasy Mini Spin Column	NucleoSpin Tissue Column
Centrifugation	6000 x g - 1 min	11000 x g - 1 min

- Discard the flow-through and reuse the **Waste Tube**

Wash & Dry silica membrane:

- 5
- To **clean DNA**, do 2 washes as described in the table below
 - Between each centrifugation, discard the flow-through and reuse the **Waste Tube**
 - To **dry silica membrane**, do an additional centrifugation described in table below

A	B	C
	Qiagen	MN
1st wash	500 µL AW1 6000 x g - 1 min	500 µL BW 11000 x g - 1 min
2nd wash	500 µL AW2 19000 x g - 2 min	600 µL B5 11000 x g - 1 min
Dry silica membrane	19000 x g - 1 min	11000 x g - 1 min

Elute:

- 6
- Place the **Spin Column** in a 2 mL annotated clean tube
 - For each sample check the correspondence of Spin Column label with 2 mL tube label
 - Add 100 µL of elution buffer (provided in the kit) to the membrane of each column





A	B
Qiagen kit	Macherey-Nagel kit
Buffer AE	Buffer BE

- Incubate at Room temperature 00:01:00
- To elute DNA from the membrane, centrifuge by respecting these instructions

A	B
Qiagen kit	Macherey-Nagel kit
6000 x g - 1 min	11000 x g - 1 min

Safety information

- **Position the lid inward the center of the centrifuge**
- **Caution do not discard the eluate**
 - Eluting with 100 μ L increases the DNA concentration but reduces the overall DNA yield

DNA storage:

7 **DNA is ready** to use immediately or store at:

- 4 °C for use a few days
- -20 °C for use within a few weeks
- -80 °C for long-term storage

Note

- Before using, allow the sample to equilibrate to Room temperature
- Perform pulse-vortexing for 00:00:05 followed by microcentrifugation

DNA quantification and DNA quality control:

8 **NanoDrop measurement:**



- Use  2 μL of **used elution buffer** as blank as standard with the NanoDrop
- To validate blank, add  2 μL to **used elution buffer** into the surface of the NanoDrop and read

Equipment

NanoDrop™ One UV-Vis Spectrophotometer

NAME

spectrophotometer

TYPE

Thermo Scientific

BRAND

ND-ONE-W

SKU

<https://www.thermofisher.com/order/catalog/product/ND-ONE-W>

LINK

Sample Volume (Metric): Minimum 1 μL ; Spectral Bandwidth: ≤ 1.8 nm (FWHM at Hg 254 nm); System Requirements: Windows™ 8.1 and 10, 64 bit; Voltage: 12 V (DC); Wavelength Range: 190–850 nm

SPECIFICATIONS

- To measure DNA concentration of sample, add  2 μL of sample DNA into the surface of the NanoDrop and read
 - Note down the A260/280, A260/230 and concentration (ng/ μL)
 - or export the data to a USB storage device

Note

- Clean the surface between each read

9 **Qubit Measurement:**

- Ensure the Qubit is calibrated before starting

Equipment

Invitrogen™ Qubit™ 3 Fluorometer

NAME

Accurately measures DNA, RNA, and protein using the highly sensitive fluorescence-based Qubit quantitation assays

TYPE

Invitrogen™ Q33216

BRAND

Q33216

SKU

<https://www.fishersci.co.uk/shop/products/qubit-3-0-fluorometer/15387293>

LINK

■ Use the

 Qubit™ dsDNA BR Assay Kit Thermo Fisher Scientific Catalog #Q32853

- Analyze  1 µL of **extracted DNA**
- Note the concentration (ng/µl)