



Apr 19, 2025

DNA extraction and Nanopore library prep for ultralong-reads

 Forked from [DNA extraction and Nanopore library prep from 15-30 whole flies- V.3.2](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.14egn98wql5d/v1>


James Hemker¹, Hannah Gellert¹, Bernard Y Kim²

¹Stanford University; ²Princeton University



James Hemker

Stanford University



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn98wql5d/v1>

Protocol Citation: James Hemker, Hannah Gellert, Bernard Y Kim 2025. DNA extraction and Nanopore library prep for ultralong-reads. [protocols.io](#) <https://dx.doi.org/10.17504/protocols.io.14egn98wql5d/v1>

Manuscript citation:



License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 13, 2024

Last Modified: April 19, 2025

Protocol Integer ID: 112003

Keywords: Drosophila, nanopore, ligation, bead-free, HMW, ultra-long, typical drosophila nanopore, nanopore library prep for ultralong, genome assembly of drosophila species, drosophila species, nanopore library prep, drosophila, dna extraction, genome assembly, genome, dna, using ont promethion sequencer, ont promethion sequencer, read n50, reads this protocol, sequencing

Funders Acknowledgements:

Bernard Y Kim

Grant ID: NIGMS F32GM135998

Dmitri A Petrov

Grant ID: NIGMS R35GM118165

Abstract

This protocol is optimized for rapid and cost-effective (about \$150) genome assembly of *Drosophila* species from laboratory lines using ONT PromethION sequencers. Following this protocol, a typical Drosophila Nanopore sequencing run should have read N50 of 20-40kbp. Sequencing is halted at about 40-60X depth of coverage (10-14 Gbp on MinKNOW for most species, assuming ~20% of data is removed by a quality filter). Adapted from Kim et al. 2021



Guidelines

This protocol is used to prepare ~1-2 ug of Nanopore library from a single reaction. The amount loaded onto the flow cell depends on the quality of the library. Larger amounts of longer libraries should be loaded to keep the molar concentration of adapted ends consistent. However, longer libraries tend to clog the flow cell more quickly, necessitating frequent DNase flushing and reloading and reducing throughput. Two libraries with the same N50 but where one has a larger number of >100kb fragments will sequence differently.

Ballpark estimates of R10.4.1 library loads maintaining good pore occupancy are:

Read N50 1kb: 10-15 ng library

Read N50 5kb: 25 ng library

Read N50 10kb: 50 ng library

Read N50 20kb: 100 ng library

Read N50 30kb: 200 ng library

Read N50 40kb+: 300 ng library

To maximize read lengths, one should not wait until all active pores have been depleted to flush and reload. A DNase flush should take place as soon as sequencing throughput starts to decrease, or about every 8 hours. A flow cell with loaded library can be stored at 4C overnight with no ill effects.



Materials

MATERIALS

⊗ 10% SDS solution

⊗ NEBNext Companion Module for Oxford Nanopore Technologies Ligation Sequencing – 24 rxns **New**
England Biolabs Catalog #E7180S

⊗ Ligation sequencing kit 1D **Oxford Nanopore Technologies Catalog #SQK-LSK109**

⊗ Chloroform **Merck MilliporeSigma (Sigma-Aldrich) Catalog #CX1055-6**

⊗ Phenol Chloroform Isoamyl Alcohol (25:24:1) Tris-saturated (pH 8.0) **Fisher Scientific Catalog #BP1752I-400**

⊗ 3M sodium acetate

⊗ Proteinase K Solution (20 mg/mL) RNA grade **Thermo Fisher Scientific Catalog #25530049**

⊗ RNase A solution **Merck MilliporeSigma (Sigma-Aldrich) Catalog #R6148**

⊗ Tris-EDTA (TE) buffer pH 8.0 1X

⊗ Homogenization Buffer (HB) [0.1M NaCl 30mM Tris-HCl pH 8.0 10 mM EDTA 0.5% Triton X-100]

⊗ Lysis Buffer (LB) [0.1M Tris-HCl pH 8.0; 0.1M NaCl; 20mM EDTA]

⊗ Hydration Buffer (STE) [400mM NaCl 20mM Tris-HCl pH 8.0 30mM EDTA]

⊗ DNase wash buffer (DWB) [300mM KCl 2mM CaCl₂ 10mM MgCl₂ 15 mM HEPES pH 8.0]

⊗ Elution Buffer (EB) [10 mM Tris-HCl pH 8.0]

⊗ Short Read Eliminator (SRE) **Circulomics Catalog #SS-100-101-01**

DNA extractions are performed in Phase lock gel tubes to minimize handling and to maximize yield. A cheaper alternative to the official phase lock gel tubes is to put ~200uL of Dow Corning High Vacuum Grease into a 2.0 mL LoBind tube with a small syringe. Care should be taken with homebrew phase lock gel tubes as using too little grease will result in the phase lock layer collapsing during the chloroform extraction step.

Although less effective, a solution of [0.8M NaCl, 9% w/v PEG8000, 10mM Tris-Cl pH 8.0] can be substituted for the Short Read Eliminator. See John Tyson's "Rocky Mountain" protocol for more details (<https://www.protocols.io/view/rocky-mountain-adventures-in-genomic-dna-sample-pr-7euhjew>). The SRE XS or XL versions can be used if DNA is short or sufficiently long. This may require a bit of trial and error to figure out.

Equipment	
DNA LoBind tubes, 1.5 mL	NAME
Tubes	TYPE
Eppendorf	BRAND
022431021	SKU
https://online-shop.eppendorf.us/US-en/Laboratory-Consumables-44512/Tubes-44515/DNA-LoBind-Tubes-PF-56252.html	LIN K
1.5 mL	SPECIFICATIONS

Equipment	
DNA LoBind tubes, 1.5 mL	NAME
Tubes	TYPE
Eppendorf	BRAND
022431048	SKU
https://online-shop.eppendorf.us/US-en/Laboratory-Consumables-44512/Tubes-44515/DNA-LoBind-Tubes-PF-56252.html	LIN K
1.5 mL	SPECIFICATIONS

Equipment	
Large-orifice pipet tips, 200uL	NAME
Pipette tips	TYPE
Fisher	BRAND
02-707-134	SKU
https://www.fishersci.com/shop/products/fisherbrand-large-orifice-pipet-tips-1-200-l-packaging-hrs-10-x-96/02707134	LIN K
200 uL	SPECIFICATIONS

Equipment	
Dounce Homogenizer, 2mL	NAME
Tissue Grinder	TYPE
Kimble	BRAND
885300-0002	SKU
https://www.kimble-chase.com/advancedwebpage.aspx?cg=886&cd=4&SKUTYPE=202&SKUFLD=SKU&DM=1250&WEBID=6856	LIN K
2 mL with Pestles A and B	SPECIFICATIONS



Equipment

5PRIME Phase Lock Gel tube, light, 2mL	NAME
Quantabio	BRAND
2302820	SKU
https://www.quantabio.com/phase-lock-gel	LINK
Light	SPECIFICATIONS



Protocol materials



NEBNext Companion Module for Oxford Nanopore Technologies Ligation Sequencing – 24 rxns **New**
England Biolabs Catalog #E7180S



PCR Tubes, 0.2mL, flat cap, natural, PCR Tube; 0.2mL; Natural; w/flat cap; 1000/Pk. **Thermo**
Fisher Catalog #3412



Nuclease-free water or water filtered using a Milli-Q filtering system **Ambion Catalog #AM9932**



Oxford Nanopore Ligation Sequencing Kit **Oxford Nanopore Technologies Catalog #SQK-LSK110**



Ligation Sequencing Kit V14 **Oxford Nanopore Technologies Catalog #SQK-LSK114**



Large-Orifice Pipet Tips 200µL **Fisher Scientific Catalog #02-707-134**



Elution Buffer (EB) [10 mM Tris-HCl pH 8.0]



Short Read Eliminator (SRE) **Circulomics Catalog #SS-100-101-01**



Agencourt AmPure XP beads **Catalog #A63880**



DNA LoBind Tubes, 1.5 mL **Eppendorf Catalog #0030108051**



DNA Precipitation Buffer (PB) [0.8 M NaCl 9% w/v PEG 8000 10mM Tris-HCl pH 8.0]



ONT Flow Cell Wash Kit **Oxford Nanopore Technologies Catalog #EXP-WSH004**



SDS, 10% Solution **Life Technologies Catalog #AM9822**



1M Tris-HCl, pH 8.0 **Thermo Fisher Scientific Catalog #15568025**



Sodium Chloride **Fisher Scientific Catalog #S271**



Proteinase K **Thermo Fisher Scientific Catalog #100005393**



RNAse A **New England Biolabs Catalog #T30318-2**



10% Triton X-100 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #648463**



0.5M EDTA, pH 7.5 **bioworld Catalog #40120779**



Nuclease-Free Water **Thermo Fisher Scientific Catalog #AM9932**



Phenol Chloroform Isoamyl Alcohol (25:24:1) Tris-saturated (pH 8.0) **Fisher Scientific Catalog #BP1752I-400**



Chloroform **Merck MilliporeSigma (Sigma-Aldrich) Catalog #CX1055-6**



Sodium acetate **P212121**



10 mM Tris-HCL pH 8.0



Before start

This protocol is for DNA extraction from whole *Drosophila*. Before starting the protocol, up to ~400 whole female flies should be starved for 1 day. Fresh flies work best for generating ultra-long reads. The nuclear extraction steps require significant input material, hence the extremely high number of flies. Testing can be done for specific flies to determine the best number for your situation.



Fly collection

1d

- 1 The day before DNA extraction, sex and sort ~200 female flies and transfer to empty fly vials with no media for starvation. Cotton plug should be soaked with water to prevent dehydration.

Tissue homogenization

15m

- 2 Prepare Homogenization Buffer (HB), Lysis Buffer (LB), and Lysis Master Mix (LMM).

Molar recipe Homogenization Buffer:

- [M] 0.1 Molarity (M) NaCl
- [M] 30 millimolar (mM) Tris-Hcl pH 8.0
- [M] 10 millimolar (mM) EDTA
- [M] 0.5 % (v/v) Triton X-100

Recipe for  50 mL Homogenization Buffer :

-  45 mL  Nuclease-Free Water Thermo Fisher Scientific Catalog #AM9932
-  2.5 mL  10% Triton X-100 Merck MilliporeSigma (Sigma-Aldrich) Catalog #648463
-  1.5 mL  1M Tris-HCl, pH 8.0 Thermo Fisher Scientific Catalog #15568025
-  1 mL  0.5M EDTA, pH 7.5 bioworld Catalog #40120779
-  0.292 g  Sodium Chloride Fisher Scientific Catalog #S271

Molar recipe for Lysis Buffer:

- [M] 0.1 Molarity (M) Tris-Hcl pH 8.0
- [M] 0.1 Molarity (M) NaCl
- [M] 20 millimolar (mM) EDTA

Recipe for  50 mL Lysis Buffer :

-  43 mL  Nuclease-Free Water Thermo Fisher Scientific Catalog #AM9932
-  5 mL  1M Tris-HCl, pH 8.0 Thermo Fisher Scientific Catalog #15568025
-  2 mL  0.5M EDTA, pH 7.5 bioworld Catalog #40120779
-  0.292 g  Sodium Chloride Fisher Scientific Catalog #S271

Recipe for Lysis Master Mix per sample (not 50-fly batches) (multiply by n+1 for n samples)

- 500 µL Lysis Buffer
- 15 µL SDS, 10% Solution [Life Technologies Catalog #AM9822](#)
- 15 µL Proteinase K [Thermo Fisher Scientific Catalog #100005393](#)
- 5 µL RNase A [New England Biolabs Catalog #T30318-2](#)

2.1 Hold the LMM on a heat block at 50 °C while doing the extractions.

3 Before beginning DNA extraction, separate the flies into groups of 50 or less, as it will make grinding them in the Dounce homogenizer much easier. If too many flies are homogenized at the same time, it will be difficult to move the tighter pestle properly. Excessive force may result in spillage of homogenate or glass breakage. Split the flies into batches to avoid this issue.

Equipment

Dounce Tissue Grinder

Tissue Grinder

Kimble Kontes

885300-0002

<https://www.dwk.com/na/kimble-kontes-dounce-tissue-grinder-with-two-glass-pestles-and-tube-2-ml-885300-0002>

NAME

TYPE

BRAND

SKU

LINK

4 Note: Perform extractions On ice to try and limit activity of endogenous DNAses that will chew up DNA. You want to get the extracted material into LMM as soon as possible.

Add 50 flies to Dounce pestle along with 500 µL HB .

5 Dounce 8-10X with the loose pestle (A), and then follow with 8-10 strokes of the tighter pestle (B). Make sure that the fly bodies are crushed, however overgrinding will lead to shorter reads.



- 6 Use a wide-bore tip to move the homogenate into a new 1.5 mL LoBind tube.

Equipment

DNA LoBind tubes, 1.5 mL

NAME

Tubes

TYPE

Eppendorf

BRAND

022431021

SKU

<https://online-shop.eppendorf.us/US-en/Laboratory-Consumables-44512/Tubes-44515/DNA-LoBind-Tubes-PF-56252.html>

LINK

1.5 mL

SPECIFICATIONS

- 7 Repeat steps 5-6 for all 50-fly batches, each going into a new LoBind tube.

Clean or switch to new Dounce pestles between samples.

- 8 Pellet debris via centrifugation for 1 minute at  500 x g

- 9 Pipette the supernatant containing the nuclei into new LoBind tubes with a wide-bore tip. Make sure to get as much of the yellow-white layer on top of the dark debris as possible.

- 10 Pellet nuclei via centrifugation for 5 minutes at  2000 x g, 4°C . If 4C centrifuge is unavailable, room temperature is OK.

- 11 Decant off the supernatant, and wipe off the white lipid smear on the side of the tube. Can use a kimwipe wrapped around tweezers.

- 12 For each sample, resuspend the first batch tube with  100 µL LMM by pipetting vigorously using a P1000. Transfer the resuspended nuclei to the next batch tube and repeat, until all nuclei for one sample are gathered in a single tube.



- 13 Add another  400 μL LMM to the final tube, bringing the total to  500 μL LMM per sample.

Lysis

4h

- 14 Incubate lysis tube at  50 $^{\circ}\text{C}$ for ~  04:00:00 . Every 45 minutes, invert the tubes to make sure that the LMM has a chance to contact all pieces of material.

4h

Note

Sometimes a bit of vigorous shaking is needed, especially if there is a lot of material. This is OK - the end product needs to be sheared for library prep anyway. Qualitative observations suggest that thorough mixing improves DNA yield and purity substantially.

Phenol chloroform extraction

1h

- 15 Spin down 1 phase lock gel tube per sample at  15000 x g for  00:00:30 .

Note

Although not essential, phase lock gel tubes help minimize shearing and loss of yield caused by repeated pipetting. Dow Corning High Vacuum Grease is compositionally identical to the light phase lock gel material. You can buy the 5.3oz tube from Amazon and squeeze some into a 10mL BD syringe for dispensing. This size of tube/syringe fits well for minimal mess and hassle. Avoid overfilling and air bubbles. Autoclave but be warned this may cause a mess, so wrap the syringe in foil beforehand.

About  250 μL of grease is placed into a 2mL LoBind tube to make the homebrew phase lock gel tube.

IMPORTANT: If an insufficient amount of grease is applied, the phase lock layer will collapse during the chloroform extraction.

Reference: <https://bitesizebio.com/18944/diy-phase-separating-gel-clean-and-cheap/>



Safety information

WARNING: If you are using normal tubes in lieu of LoBinds, do not use polystyrene tubes for the phenol-chloroform extraction. They will melt and burst in the centrifuge. Polypropylene tubes do not melt.

Equipment

5PRIME Phase Lock Gel tube, light, 2mL	NAME
Quantabio	BRAND
2302820	SKU
https://www.quantabio.com/phase-lock-gel	LINK
Light	SPECIFICATIONS

16 Transfer all of the lysate solution from the heat block to the phase lock gel tube by pipetting with a wide-bore tip.

17 Add an equal volume (~  400 μ L) of

 Phenol Chloroform Isoamyl Alcohol (25:24:1) Tris-saturated (pH 8.0) **Fisher Scientific Catalog #BP1752I-400**

to the phase lock tube.

Safety information

This should be performed inside the fume hood.

Note

We recommend shaking the bottle of phenol first to remix, and then letting it sit for ~5 minutes. It will then separate into the phenol on the bottom and the tris buffer on top. Make sure to pipette from the phenol layer.



- 18 Mix by placing tubes on a rocker at medium speed for 00:08:00 .

8m

Note

We use a rocking platform, so the tubes are placed on their sides horizontally to maximize the surface area. When solution is well mixed, aqueous (top) layer will be a cloudy milky color.

- 19 Centrifuge the phase lock tube at 10000 x g for 00:08:00 . Phase lock layer should now separate aqueous and organic layers.

8m

- 20 Add an equal volume (400 μ L) of

Chloroform **Merck MilliporeSigma (Sigma-Aldrich) Catalog #CX1055-6** to the tube.

Safety information

This step should be performed inside the fume hood.

- 21 Mix by placing tubes on a rocker at medium speed for 00:08:00 .

8m

- 22 Centrifuge the phase lock tube at 15000 x g for 00:08:00 . Phase lock layer should now separate aqueous and organic layers.

8m

- 23 Quickly decant the aqueous (top) layer into a fresh 2.0 mL LoBind tube.



Note

Try to perform the decanting step in a few seconds, and don't tap/shake the phase lock tube to get the last drops out. Care must be taken as the chloroform significantly weakens the phase lock gel layer. If the phase lock tube is inverted for too long during decanting, the layer will collapse and everything will pour out. It's best to leave a couple of drops behind but avoid the hassle of cleaning this up.

IMPORTANT: It is highly recommended to use LoBind tubes in this and subsequent steps. The coating will prevent DNA sticking to the tube. This is helpful for maximizing yield and minimizing shearing.

Safety information

This step should be performed inside the fume hood.

Equipment

DNA LoBind Tube	NAME
2 mL	TYPE
Eppendorf	BRAND
022431048	SKU
https://www.eppendorf.com/us-en/Products/Laboratory-Consumables/Tubes/DNA-LoBind-Tubes-p-022431048	LINK

DNA precipitation, wash, and resuspension

1h 30m

- 24 Chill 100% ethanol on ice and additionally make fresh  1 mL 70% ethanol using nuclease-free water, per sample.

25 Add 10% volume (~  50 μL) of [M] 3 Molarity (M)  Sodium acetate **P212121** to each sample.

26 Add 2-2.5x volumes (typically  1000 μL) of cold 100% ethanol to the tube, and mix with careful swirling and gentle rocking.

Expected result

DNA should slowly precipitate into a single white stringy clump, and un-precipitated DNA should be visible as shimmering strands at the bottom of the tube that are attached to the white clump.

Note

If the extraction tube turns cloudy, it is likely salt precipitation because the solution is too nonpolar and not DNA. Add water dropwise with thorough mixing and the solution should clear up.

27 Using a P200 pipette and a wide-bore tip, transfer the stringy clump to a fresh 1.5 mL LoBind tube.

Note

This step can be somewhat tricky. The DNA clump can get stuck to your pipette tip and be very difficult to get off. Making sure that the DNA is fully precipitated (no un-precipitated strands present) and quick pipetting helps prevent this. Position the pipette tip right above the clump and aspirate quickly to bring the clump into the tip. Using a similar quick motion, dispense the DNA clump and liquid into the fresh 1.5 mL tube.

28 Add  150 μL (or enough to cover the DNA) of 70% ethanol.

29 Remove and discard the ethanol. Be careful not to discard any DNA.

30 Add  150 μL 70% ethanol to wash the DNA. Gently swirl to mix.



31 Centrifuge the tube at  2000 x g for  00:05:00 . 5m

32 While being careful not to disturb the pellet, pipette off the ethanol.

33 Add  175 μ L 70% ethanol

34 Spin at  10000 rcf for  00:01:00 . 1m

35 Being careful not to disturb the DNA pellet, remove the ethanol.

36 Add  200 μ L 70% ethanol

37 Spin at  10000 rcf for  00:01:00 . 1m

38 Remove the ethanol.

39 Spin at  10000 rcf for  00:01:00 .

40 Using a 10uL pipette, remove any excess ethanol.

41 Allow the DNA to air dry right until the moment it becomes translucent (~  00:05:00).

Do not over-dry the pellet.

Note

Oftentimes the whole DNA pellet will not become translucent but the edges of the pellet will.



- 42 Resuspend in  60 μL  10 mM Tris-HCL pH 8.0 and incubate at  50 $^{\circ}\text{C}$ for at least  01:00:00 .

1h

Note

Usually resuspend in  60 μL , but can range from  31 μL -  90 μL depending on the amount of material

- 43 Briefly spin down tube to gather any condensation and store at  4 $^{\circ}\text{C}$.

DNA resuspension

1w

- 44 Keep the DNA at 4C for at least a few days to obtain proper resuspension. If the DNA blob is not resuspending (which will be likely with very long reads), mix the DNA twice, very slowly, with a P1000 pipette tip. This will encourage DNA to resuspend and make sure it is adequately sheared for library prep.

Note

While it is possible to over-shear the DNA, under-shearing makes for much worse Nanopore runs because it makes the sample too viscous. More shearing may be needed depending on how much the sample was mixed during incubation, or how fresh the flies were.

If there is a lot of DNA that doesn't seem to be resuspending, shearing with a blunt-end needle may be necessary. When using a blunt-end needle, be sure to 'unstick' the plunger before using it on DNA. Aspirate DNA into the needle (not the plunger) 3x times maximum.



- 45 Check sample concentration and quality of  1 μL aliquots using Qubit and Nanodrop.

Note

Ideally, this should Qubit at $>75 \text{ ng}/\mu\text{L}$ and have Nanodrop ratios of $260/280 >1.8$ and $260/230 >2.0$. If sample is above $150 \text{ ng}/\mu\text{L}$, dilute with more 10mM Tris to $150 \text{ ng}/\mu\text{L}$ before SRE.

Short Read Elimination 1

- 46 Using a cut-off P200 tip (a wide bore will be too small to fit in the PCR tube), gently transfer  30 μL of sample to a 1.5 mL DNA LoBind tube. Add  30 μL SRE XL buffer. Using a wide-bore P200 tip, quickly but gently mix the tube. The precipitation buffer described here can be used in place of the SRE buffer but is not as effective at removing small DNA fragments as SRE.

Note

Dilute sample down to $150 \text{ ng}/\mu\text{L}$ as possible before mixing with buffer

 Short Read Eliminator (SRE) **Circulomics Catalog #SS-100-101-01**

 DNA LoBind Tubes, 1.5 mL **Eppendorf Catalog #0030108051**

 DNA Precipitation Buffer (PB) [0.8 M NaCl 9% w/v PEG 8000 10mM Tris-HCl pH 8.0]

 Large-Orifice Pipet Tips 200 μL **Fisher Scientific Catalog #02-707-134**

- 47 Centrifuge the sample at  10000 $\times g$ for  00:30:00 or until DNA has pelleted and solution is no longer viscous. Meanwhile, prepare  500 μL fresh 70% ethanol with nuclease-free water.



48 Pipette off the supernatant, taking care not to disturb the DNA pellet.

Note

We have increased our yield by leaving 10-15 μL of supernatant in the bottom of the tube going into the first wash. This is particularly important if the pellet is invisible.

49 Add  150 μL of 70% ethanol. Pipette slowly, with the tip touching the front wall of the tube, so that the pellet is not disturbed.

50 Centrifuge at  10000 x g for  00:02:00 .

51 Pipette off the supernatant, taking care not to disturb the DNA pellet. Make sure all the supernatant is removed and only the pellet remains.

52 Repeat wash: 

Note

The second centrifuge (step 43) can be shorter, ~1 minute.

53 Briefly spin sample and use a P10 to remove any remaining ethanol.

54 Resuspend pellet in  48.5 μL EB.

 Elution Buffer (EB) [10 mM Tris-HCl pH 8.0]

55 Incubate the tube on the heat block at  50 $^{\circ}\text{C}$ for at least  01:00:00 . Briefly spin down the tube to collect condensation. Incubate at least  48:00:00 and preferably for twice that time at  4 $^{\circ}\text{C}$.

Note

Before proceeding to the next stage, we recommend  1 μL Qubit to ensure there is still DNA.

DNA repair and end-prep

- 56 Thaw NEBNext repair and DNA-tailing mixes and buffers from the Nanopore Companion Module. Vortex buffers and flick mixes after thawing. Spin down tubes and keep chilled on ice.



NEBNext Companion Module for Oxford Nanopore Technologies Ligation Sequencing – 24 rxns **New England Biolabs Catalog #E7180S**

- 57 Add  3.5 μL of FFPE DNA Repair Buffer,  3.5 μL of End-Prep Reaction Buffer,  2 μL of FFPE DNA Repair Mix, and  3 μL of End-Prep Reaction Mix to a PCR tube. Add remaining  47.5 μL of DNA sample to the PCR tube using a cut-off P200 tip. Mix tube with gentle flicking (or very gentle pipetting with the cut-off P200 tip), and then briefly spin down.

Note

To increase efficiency and decrease amount of pipette tips needed, prep PCR tubes with buffers and mixes and add the HMW DNA sample last.

Note

We have found that half-reaction volumes for this step do not decrease yield. To follow half-reaction volumes, see the Single Fly Forked Protocol for more information.



PCR Tubes, 0.2mL, flat cap, natural, PCR Tube; 0.2mL; Natural; w/flat cap; 1000/Pk. **Thermo Fisher Catalog #3412**



Nuclease-free water or water filtered using a Milli-Q filtering system **Ambion Catalog #AM9932**



- 58 In a thermal cycler, incubate at 20 °C for 01:00:00 then 65 °C for 00:30:00 . After this, sample can be held at 4 °C temporarily until ready to proceed.

1h 30m

Bead Clean Up

- 59 Prepare 500 µL of 70% ethanol per sample.
- 60 Transfer sample from PCR tube to a LoBind tube using a cut-off P200 tip.
- 61 Add equal volume of AmPure XP beads (normally 59.5 µL) to the sample. Immediately use a P200 wide bore tip to mix 5x.

Note

This step must be performed quickly; otherwise, DNA will precipitate onto pipette tip and will result in sample loss.

Note

If needed, briefly spin down to ensure there are no bubbles or any sample on the wall of the LoBind tube.

Agencourt AmPure XP beads **Catalog #A63880**

- 62 Incubate at room temperature for 00:05:00 .
- 63 Place the LoBind tube on a magnetic rack and wait until solution is clear and the beads are pelleted.
- 64 Remove the supernatant by placing pipette tip on the wall of the LoBind tube opposite of the beads. Pipette very slowly to ensure no DNA is pulled off.

5m

**Note**

If DNA is pulled off, add supernatant back to tube and wait for solution to clear. Then try again.

- 65 Wash by adding  190 μL of 70% ethanol (enough to cover the beads on the wall of the LoBind tube).

Note

Work quickly to add the 70% ethanol at this step to avoid the beads drying out.

- 66 Remove and discard ethanol.

- 67 Wash again by adding  200 μL of 70% ethanol.

- 68 Remove and discard ethanol.

Note

Briefly spin and use a P10 pipette to remove any remaining excess of ethanol.

- 69 Resuspend sample in  32 μL of nuclease-free water.

- 70 Place tube on heat block at  50 $^{\circ}\text{C}$ until the pellet has dissolved.

Note

This step can take a long time. If there is concern about the DNA not resuspending off the beads, the tube can be stored at  4 $^{\circ}\text{C}$ overnight and then the sample removed from the beads the following morning.



71 Place LoBind tube on magnet rack until solution is clear.

72 Using a cut off P200 tip, remove the supernatant containing the aqueous DNA.

Note

Qubit  1 μL of sample to ensure DNA concentration before proceeding to next step.

Note

This is a safe stopping point. Sample can be stored at  4 $^{\circ}\text{C}$.

Adapter ligation

73 Thaw AMXF, Quick T4 Ligase, LNB, and LFB from the NEBNext Nanopore Companion Module and the Nanopore LSK110 kit or LA, Quick T4 Ligase and LNB from Nanopore LSK114 kit . Mix AMXF or LA, Quick T4 ligase, and LFB by flicking. Mix LNB by pipetting. Briefly spin the tubes down and keep chilled on ice.



NEBNext Companion Module for Oxford Nanopore Technologies Ligation Sequencing – 24 rxns **New England Biolabs Catalog #E7180S**



Oxford Nanopore Ligation Sequencing Kit **Oxford Nanopore Technologies Catalog #SQK-LSK110**



Ligation Sequencing Kit V14 **Oxford Nanopore Technologies Catalog #SQK-LSK114**

74 Add  30 μL prepared DNA sample (the extra  1 μL can be used to Qubit),  2.5 μL AMXF or LA, and  5 μL Quick T4 ligase to a fresh 1.5 mL DNA LoBind tube. Gently flick the tube to mix.

75 Add  12.5 μL LNB to the sample. Working quickly, mix by gentle pipetting with a wide-bore tip. DNA precipitation is normal, but if the DNA precipitates before you finish mixing it will stick to your pipette tip and you will lose a significant amount of library.



Large-Orifice Pipet Tips 200µL **Fisher Scientific Catalog #02-707-134**

76 Incubate the reaction mixture at room temperature for 00:20:00 . 20m

77 Add 20 µL of AmPure XP beads and mix quickly with wide-bore tip.

78 Incubate at room temperature for 00:05:00 . 5m

79 Place tubes on magnetic and wait for solution to clear.

80 On the magnet, use a cut-off P200 tip to pull the supernatant off the beads very slowly, then dispense the supernatant back onto the bead pellet slowly. Let the sample sit on the magnet for a few minutes.

81 Pipette off supernatant with a normal pipette tip, pipetting from the front of the tube away from the pellet.

82 Add 95 µL of LFB to the tube. SFB or a 1:1 dilution of PB can be used here.

Note

DO NOT USE ETHANOL TO WASH PREPARED LIBRARY. It will denature the motor protein.

Lightly tap the tube to encourage adapter on the beads to come off, but not necessarily for beads to resuspend.

83 Being careful not to disturb the pellet, pipette off all the supernatant.

84 Wash again using 105 µL of LFB. Pipette LFB on to the beads more quickly to get the pellet off the side of the tube. Lightly tap the tube to mix but not fully resuspend.



85 While on magnet remove LFB. Briefly spin and use a P10 pipette to remove any remaining excess of LFB.

86 Resuspend pellet in  21 μL EB for R9.4.1 sequencing or in  30 μL for R10.4.1 sequencing.

 Elution Buffer (EB) [10 mM Tris-HCl pH 8.0]

87 Incubate library on the heat block at  34 $^{\circ}\text{C}$ for  01:00:00 . Briefly spin down the tube to collect condensation then incubate for at least  48:00:00 before the next step.

2d 1h

88 Place sample on magnet wait until solution is clear. Use a cut off P200 tip to remove sample from beads and place in a new 1.5mL Lo Bind tube.

(Optional) Library size selection with SRE buffer

89 Quantify library concentration using  1 μL of the prepared library with Qubit. This step should not be performed unless library concentration is greater than 40 ng/uL. If the concentration is greater than 100ng/uL the library should be diluted to improve size selection performance.

90 Add an equal volume ( 20 μL) of SRE XL buffer to the library and gently pipette mix using a wide-bore tip.

 Short Read Eliminator (SRE) **Circulomics Catalog #SS-100-101-01**

 Large-Orifice Pipet Tips 200 μL **Fisher Scientific Catalog #02-707-134**

91 Centrifuge at  10000 x g, 00:30:00 .

92 Pipette off the supernatant, being careful not to disturb the DNA pellet at the bottom of the tube.

**Note**

Similar to previous SRE step, leave 10-15uL of supernatant in the bottom of the tube for the first wash.

- 93 Add 100 uL of LFB, SFB, or 1:1 diluted PB (similar to step 46) to wash the pellet. It does not really matter which one is used.

Note

DO NOT USE ETHANOL TO WASH PREPARED LIBRARY. It will denature the motor protein.

- 94 Centrifuge tube at  10000 x g, Room temperature for  00:02:00 .

- 95 Being careful not to disturb the pellet, pipette off all the supernatant.

- 96 Repeat wash: 

- 97 Resuspend pellet in  21 µL EB for R9.4.1 sequencing or in  30 µL for R10.4.1 sequencing.

 Elution Buffer (EB) [10 mM Tris-HCl pH 8.0]

- 98 Incubate the tube on the heat block at  37 °C for  01:00:00 . Briefly spin down the tube to collect condensation, and incubate at least  48:00:00 at  4 °C before sequencing.

Tips for sequencing the library-R10.4.1**30m**

- 99 Thaw 1 tube SB (LSK110), 1 tube LIS, 1 tube of FCF per sample, and 1 tube FCT. Mix SQB and FB by flicking.



100 Warm the FCF at 37 °C for 00:30:00

30m

101 Add 30 µL FCT to FCF and pipette 10x to ensure thorough mixing

Note

We recommend marking the top of the FCF tube after FCT has been added.

102 Follow the official instructions to prime the flow cell.

103 While the flow cell is priming, prepare the library by adding 70 µL of LIS, Sample library (LIS and library should total to 100uL), and 100 µL SB. Lightly tap to mix until swirls disappear but wait to pipette mix until just before loading.

104 Pipette mix prepared library 2x times and then following official instructions to load the flow cell.

105 Over the course of a sequencing run, pores will get clogged and become inactive. It is essential to flush the flow cell at 10-14 hour intervals to make these pores available again. We recommend Nanopore's Flow Cell Wash Kit (EXP-WSH00).

ONT Flow Cell Wash Kit **Oxford Nanopore Technologies Catalog #EXP-WSH004**



Protocol references

Kim, Bernard Y., Jeremy R. Wang, Danny

E. Miller, Olga Barmina, Emily Delaney, Ammon Thompson, Aaron A.

Comeault et al. "Highly contiguous assemblies of 101 drosophilid genomes." *Elife* 10 (2021): e66405.

Kim, B.Y., Gellert, H.R., Church, S.H., Suvorov, A., Anderson, S.S.,

Barmina, O., Beskid, S.G., Comeault, A.A., Crown, K.N., Diamond, S.E.

and Dorus, S., 2024. Single-fly genome assemblies fill major

phylogenomic gaps across the Drosophilidae Tree of Life. *PLoS Biology*, 22(7), p.e3002697