

Sep 12, 2022

DNA Extraction

DOI

<https://dx.doi.org/10.17504/protocols.io.5qpvojbw7l4o/v1>

klaus.stark¹, Kira Julia Stanzick¹

¹Department of Genetic Epidemiology, University of Regensburg



Kira Julia Stanzick

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5qpvojbw7l4o/v1>

Protocol Citation: klaus.stark , Kira Julia Stanzick 2022. DNA Extraction. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.5qpvojbw7l4o/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 19, 2022



Last Modified: September 12, 2022

Protocol Integer ID: 68916

Keywords: dna extraction dna extraction from frozen blood clot, gentra puregene blood kit for dna extraction, gentra puregene blood kit, dna extraction, dna extraction dna extraction, frozen blood clot, qiagen clotspin basket, ml whole blood, whole blood, dna, qiagen

Abstract

DNA extraction from frozen blood clots is challenging. Here we applied QIAGEN Clotspin Baskets and the Gentra Puregene Blood Kit for DNA extraction from 5.5 ml whole blood without anticoagulating additives.

Protocol materials

- ⊗ RBC Lysis Solution **Qiagen**
- ⊗ Cell Lysis Solution **Qiagen**
- ⊗ Proteinase K **Qiagen**
- ⊗ Protein Precipitation Solution **Qiagen**
- ⊗ Isopropanol
- ⊗ Proteinase K (20 mg/ml)
- ⊗ Isopropanol
- ⊗ ethanol
- ⊗ DNA Hydration Solution **Qiagen**



Blood clot preparation and red blood cell lysis

40m 19s

- 1 Frozen blood clots in 15 ml tubes were transferred from -20°C to a warming cabinet 55 °C for 00:10:00 and thereafter immediately placed On ice 10m
- 2 The tube was inverted to loosen the clot. The blood clot was completely poured with 5 mL RBC Lysis Solution **Qiagen** into the Clotspin Basket placed on a 50 ml tube
- 3 To disperse the clot, centrifugate 2000 x g, 00:05:00 5m
- 4 The remaining clot material from the Clotspin basket was transferred through the basket to the filtrate with 00:10:00 RBC Lysis Solution **Qiagen** and the basket was discarded. 10m
- 5 To completely disperse the clotted material, the filtrate was **vortexed** vigorously for 00:00:03 and placed for 00:05:00 at Room temperature on a circulating shaker (**250 1/min**). 5m 3s
- 6 The tubes were again **vortexed** vigorously for 00:00:03 and centrifugated 2000 x g, 00:05:00 at 2000 x g for 5 min 5m 3s
- 7 The supernatant was carefully discarded, taking care that the pellet remains in the tube. If no pellet was visible, about 0.5 mL of the supernatant was kept in the tube
- 8 The tube was **vortexed** rigorously for 00:00:10 and additional 5 mL RBC Lysis Solution **Qiagen** was added to the pellet, followed by **vortexing** for 00:00:03 and incubation on a circulating shaker (**250 1/min**) for 00:05:00 at Room temperature 5m 13s

White blood cell lysis

1d 0h 5m 30s

- 9 Centrifugation 2000 x g, 00:05:00 to pellet the DNA-containing white blood cells 5m



10 Discard supernatant carefully, leaving about 0.2 ml of residual liquid. Vortex rigorously for 00:00:10 10s

11 Adding of 5 mL Cell Lysis Solution **Qiagen** and 25 μ L of Proteinase K **Qiagen** [M] 20 mg/mL as followed by rigorously **vortexing** for 00:00:10 10s

12 For complete lysis of DNA containing cells, incubate Overnight at 55 °C 10s

Protein precipitation

17m 20s

13 samples were cooled down On ice for 00:05:00 and 1.7 mL Protein Precipitation Solution **Qiagen** was added, followed by rigorous **vortexing** for 00:00:20 5m 20s

14 After centrifugation 2000 x g, 00:10:00 , incubate the samples for 00:02:00 On ice . 12m

DNA precipitation

3m

15 For precipitation of DNA, the supernatant was **carefully** transferred in a 50 ml tube, containing 5 mL Isopropanol

16 The samples were mixed by gently **inverting the tube 50 times** and centrifugated 2000 x g, 00:03:00 3m

17 The supernatant was **carefully** discarded, and the tube was drained on a clean piece of absorbent paper without losing the pellet. The DNA pellet was inspected visually and with two following options

STEP CASE

big red to brown pellet 9 steps



Additional washing step to remove remaining erythrocytes

DNA washing

4h 23m 20s

- 18 2 mL of Cell Lysis Solution **Qiagen** and 10 μ L of Proteinase K (20 mg/ml) were added, followed by incubation at 55 $^{\circ}$ C for 02:00:00 in a warming cabinet **or** at Room temperature Overnight on a circulating shaker (**250 1/min**) 4h
- 19 samples were cooled down On ice for 00:05:00 and 1 mL Protein Precipitation Solution **Qiagen** was added, followed by rigorous **vortexing** for 00:00:20 5m 20s
- 20 After centrifugation 2000 x g, 00:10:00 , the samples were incubated for 00:02:00 On ice 12m
- 21 DNA precipitation was done with 3 mL Isopropanol . The samples were mixed by gently **inverting the tube 50 times** and centrifuged 2000 x g, 00:03:00 3m
- 22 The supernatant was carefully discarded. The tube was drained on a clean piece of absorbent paper, taking care that the pellet remained in the tube
- 23 5 mL of [M] 70 % (v/v) ethanol were added immediately and the tubes were inverted until the pellet was detached
- 24 After centrifugation 2000 x g, 00:03:00 at 2000 x g for 3 min, the supernatant was carefully discarded, and the DNA pellet was air dried at room temperature for 10 min or until the pellet got glassy 3m

DNA hydration

1h 2m

- 25 Addition of 0.5 mL DNA Hydration Solution **Qiagen** for a **large** pellet **or** 0.3 mL for a **smaller** pellet was followed by incubation at 65 $^{\circ}$ C in a warming 1h 2m



cabinet for  01:00:00 . To fully dissolve the DNA, the sample were put on a **shaker**

 Overnight at  Room temperature

- 26 Samples were centrifugated briefly and the solved DNA was transferred to a 2 ml cup and **stored** at  -20 °C .