

Oct 13, 2022

DNA extraction for fermented plant based foods

DOI

<https://dx.doi.org/10.17504/protocols.io.36wgqj9wkvk5/v1>

Anique Ahmad¹, Arya Gautam¹, Tsedenia Deneke¹, Ahmed Shibl¹, Aashish Jha¹

¹New York University, Abu Dhabi

Anique Ahmad: aa5506@nyu.edu

Aashish Jha: jhaar@nyu.edu



Genetic Heritage Group NYU Abu Dhabi

NYU Abu Dhabi

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.36wgqj9wkvk5/v1>

Protocol Citation: Anique Ahmad, Arya Gautam, Tsedenia Deneke, Ahmed Shibl, Aashish Jha 2022. DNA extraction for fermented plant based foods. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.36wgqj9wkvk5/v1>

Manuscript citation:

Anique Ahmad, Arya Gautam, Tsedenia Deneke, Ahmed, Shibl, Aashish Jha (2022). DNA extraction from fermented plant based foods. protocols.io

dx.doi.org/10.17504/protocols.io.36wgqj9wkvk5/v1



License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: September 19, 2022

Last Modified: October 13, 2022

Protocol Integer ID: 70218

Keywords: based foods dna extraction, foods dna extraction, dna extraction, fermented plant, dna, based food

Funders Acknowledgements:

The Bill and Melinda Gates Foundation

Abstract

DNA Extraction for fermented plant based foods

Guidelines

Disinfect mortar and pestle with 70% ethanol in between samples

If the goal is metagenomic sequencing, we recommend bypassing the mortar and pestle step. Simply suspend the samples in 1 ml of ZymoBIOMICS™ Lysis Solution, shake vigorously or vortex, and transfer 750 uL of the supernatant to a ZR BashingBead™ Lysis Tube (0.1 & 0.5 mm) and proceed to Step 2.

Materials

ZymoBIOMICS DNA Extraction Kit

Fecal samples

Mortar and pestle

Liquid nitrogen

Ice box with dry ice

1.5 ml microcentrifuge tubes

Vortex

Omni Bead Ruptor Elite bead beater

Centrifuge

Safety warnings

! Handle liquid nitrogen with proper training and PPE



Before start

- Label the following:
 - ZR BashingBead™ Lysis Tubes (0.1 & 0.5 mm)
 - ZymoSpin™ III-F Filter in Collection Tube
 - Collection Tube
 - ZymoSpin™ II-CR Column in Collection Tube
 - ZymoSpin™ III-HRC Filter in Collection Tube
 - 2 Sets of 1.5 ml microcentrifuge tubes (not provided with kit)
- Spin the labelled ZR BashingBead™ Lysis Tubes (0.1 & 0.5 mm) for 10 seconds in a mini-centrifuge to ensure that the beads have settled at the bottom
- Include 1 control per batch and assign its position randomly.



- 1 Take  250 mg of food sample while on dry ice and transfer to mortar. Add liquid nitrogen until food sample is completely immersed. Using pestle vigorously homogenize sample until a powder forms. Transfer all of the sample to a *ZR BashingBead™ Lysis Tube (0.1 & 0.5 mm)*. Add  750 μL ZymoBIOMICS™ Lysis Solution to the tube. Cap tightly.
- 2 Secure in a Omni Bead Ruptor Elite bead beater fitted with a 2 ml tube holder assembly and process at max speed (30 m/s) for  00:05:00 . Rest for  00:05:00 . 10m
- 3 Centrifuge the *ZR BashingBead™ Lysis Tubes (0.1 & 0.5 mm)* in a microcentrifuge at $\geq$  10000 x g, 00:01:00 . 1m
- 4 Transfer up to  600 μL supernatant to the *Zymo-Spin™ III-F Filter* in the *labelled Collection Tube* and centrifuge at  8.000 x g, 00:01:00 . Discard the Zymo-Spin™ III-F Filter. 1m
- 5 Add  600 μL of ZymoBIOMICS™ DNA Binding Buffer to the filtrate in the *labelled Collection Tube* from Step 4. Mix well.
- 5.1 Repeat Step 5 so that the final volume of ZymoBIOMICS™ DNA Binding Buffer added is  1200 μL .
- 6 Transfer  800 μL of the mixture from Step 5 to a *Zymo-Spin™ IICR Column* in a *Collection Tube* and centrifuge at  10.000 x g, 00:01:00 . 1m
- 6.1 Discard the flow through from the *Collection Tube* and repeat Step 6.
- 7 Add  400 μL ZymoBIOMICS™ DNA Wash Buffer 1 to the *Zymo-Spin™ IICR Column* in a *new Collection Tube* and centrifuge at  10.000 x g, 00:01:00 . Discard the flow-through. 1m
- 8 Add  700 μL ZymoBIOMICS™ DNA Wash Buffer 2 to the *Zymo-Spin™ IICR Column* in a *Collection Tube* and centrifuge at  10.000 x g, 00:01:00 . Discard the flow-through. 1m



- 9 Repeat with  200 μL ZymoBIOMICS™ DNA Wash Buffer 2 to the *Zymo-Spin™ IICR Column* in a Collection Tube and centrifuge at  10.000 x g, 00:01:00 . 1m
- 10 Transfer the *Zymo-Spin™ IICR Column* to a clean 1.5 ml microcentrifuge tube and add  100 μL (50 μL minimum) ZymoBIOMICS™ DNase/RNase Free Water directly to the column matrix and incubate for  00:05:00 . Centrifuge at  10.000 x g, 00:01:00 to elute the DNA 6m
- 11 Place a *Zymo-Spin™ III-HRC Filter* in a new *Collection Tube* and add  600 μL ZymoBIOMICS™ HRC Prep Solution. Centrifuge at  8.000 x g, 00:03:00 . 3m
- 12 Transfer the eluted DNA (Step 10) to a prepared *Zymo-Spin™ III-HRC Filter* in a clean 1.5 ml microcentrifuge tube and centrifuge at exactly  16.000 x g, 00:03:00 . 3m
- 13 The filtered DNA is now suitable for PCR and other downstream applications. Eluted DNA should be frozen (–30 to –15°C or –90 to –65°C)