

Jul 21, 2023

## DNA Extraction for Formalin Specimen

DOI

[dx.doi.org/10.17504/protocols.io.e6nvwdk97lmk/v1](https://dx.doi.org/10.17504/protocols.io.e6nvwdk97lmk/v1)

heeseung<sup>1</sup>

<sup>1</sup>Rouse lab



heeseung

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvwdk97lmk/v1>

**Protocol Citation:** heeseung 2023. DNA Extraction for Formalin Specimen. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.e6nvwdk97lmk/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** In development

We are still developing and optimizing this protocol

**Created:** July 21, 2023



**Last Modified:** July 21, 2023

**Protocol Integer ID:** 85314

**Keywords:** dna extraction for formalin specimen dna extraction, formalin specimen dna extraction, qiaamp dna ffpe advanced kit, dna extraction, formalin specimen, qiagen kit, formalin, extraction, dna

## Abstract

DNA Extraction for Formalin Specimen (Qiagen Kit: QIAamp DNA FFPE Advanced Kit)

## Troubleshooting



## Deparaffinization process

- 1 Place the tissue in a 1.5 ml or 2 ml microcentrifuge tube. Add **300 µl Deparaffinization Solution**, **vortex** vigorously for **10 s**, and centrifuge briefly to bring the sample to the bottom of the tube.

**Note:** For specimens stored only in alcohol (=not fixed in paraffin), "deparaffinization step" is not required. In this case, instead of 300 µl Deparaffinization Solution, add **300 µl ATL**, **300 µl AL**, and **300 µl ethanol (96–100%)**, and proceed to the next **step 3**, without step 2 (incubation).

- 2 Incubate at **56°C** for **3 min**, then allow to **cool to room temperature**.

## Preparation for Extraction 1

- 3 Add **25 µl Buffer FTB**, **55 µl RNase-free Water**, and **20 µl Proteinase K**. Mix by vortexing. Briefly centrifuge the tube to spin down any tissue that sticks to the tube wall or under the cap of the tube after vortexing.

## Lysis

- 4 Incubate for **2 h** at **56°C** and **1000 rpm**.

**Note:** In general, more than one hour is required, and the results were most successful when the process was conducted for more than 2 hours. For relatively older samples, it can be proceed "overnight".

- 5 Incubate for **1 h** at **90°C** without shaking.

## Preparation for Extraction 2

- 6 Carefully remove and discard the upper **blue phase**. Keep the lower aqueous lysate, add **150 µl RNase-free Water**, then vortex.

**Note:** If the deparaffinization process (step 1 and 2) is omitted, remove the **900 µl of the upper lysate**.

- 7 Add **2 µl RNase A**, vortex, and incubate for **2 min** at **room temperature** on the bench.



- 8 Add **20 µl Proteinase K**, vortex, and incubate for **15 min** at **65°C** and **450 rpm**.
- 9 Add **250 µl Buffer AL** and **250 µl ethanol (96–100%)** to each sample and mix thoroughly by vortexing

## Extraction

- 10 Transfer **450 µl lysate** to the QIAamp UCP MinElute column (in a 2 ml collection tube), and centrifuge at **15,000 rpm** for **30 s**.  
**Note:** Maximum speed is recommended.
- 11 Transfer the residual lysate to the same QIAamp UCP MinElute column, and centrifuge at **15,000 rpm** for **1 min**. Discard the flow-through from step 10 and 11 and reuse the collection tube.  
**Note:** Maximum speed is recommended.
- 12 Add **500 µl Buffer AW1** to each spin column, and centrifuge at **15,000 rpm** for **30 s**. Discard the flow-through and reuse the collection tube.  
**Note:** Maximum speed is recommended.
- 13 Add **500 µl Buffer AW2** to each spin column, and centrifuge at **15,000 rpm** for **30 s**. Discard the flow-through and reuse the collection tube.  
**Note:** Maximum speed is recommended.
- 14 Add **250 µl ethanol (96–100%)** to the spin column, and centrifuge at **15,000 rpm** for **30 s**. Discard the flow-through and collection tube. Place the spin column into a new 2 ml collection tube and centrifuge for **3 min at full speed** to remove any residual liquid to dry the membrane.  
**Note:** Maximum speed is recommended.
- 15 Place the QIAamp UCP MinElute column into a clean 1.5 ml microcentrifuge tube, and discard the collection tube containing the flow-through. Open the lid of the QIAamp MinElute column and apply **20 µl Buffer ATE** to the center of the membrane.
- 16 Close the lid and incubate at **room temperature** for **1 min**, then centrifuge at **full speed** for **1 min** to elute the DNA. **Step 15** is repeated **twice** more so that final volume is **60 µl (20 µl+20 µl+20 µl)**.



**Note:** The final volume can be up to 100  $\mu$ l, but the sequence result in the case of 60  $\mu$ l was the most successful.