


May 16, 2018

DNA Extraction for plant samples by CTAB

 [GigaScience](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.pzqdp5w>



Ting Yang¹, Chenyu Wu¹

¹BGI-Shenzhen

GigaScience Press

BGI



Ting Yang

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.pzqdp5w>

External link: <https://doi.org/10.1093/gigascience/giz007>

Protocol Citation: Ting Yang, Chenyu Wu 2018. DNA Extraction for plant samples by CTAB. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.pzqdp5w>

**Manuscript citation:**

Huan Liu, Jinpu Wei, Ting Yang, Weixue Mu, Bo Song, Tuo Yang, Yuan Fu, Xuebing Wang, Guohai Hu, Wangsheng Li, Hongcheng Zhou, Yue Chang, Xiaoli Chen, Hongyun Chen, Le Cheng, Xuefei He, Hechen Cai, Xianchu Cai, Mei Wang, Yang Li, Sunil Kumar Sahu, Jinlong Yang, Yu Wang, Ranchang Mu, Jie Liu, Jianming Zhao, Ziheng Huang, Xun Xu, Xin Liu, Molecular digitization of a botanical garden: high-depth whole-genome sequencing of 689 vascular plant species from the Ruili Botanical Garden, *GigaScience*, Volume 8, Issue 4, April 2019, giz007, <https://doi.org/10.1093/gigascience/giz007>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 08, 2018

Last Modified: May 21, 2018

Protocol Integer ID: 12048

Keywords: dna extraction for plant sample, ctab dna extraction, dna extraction, leaves by ctab, plant sample, stored leaf, extraction, leaf, ctab, dna

Abstract

DNA extraction of -80°C stored leaves by CTAB.

Materials

MATERIALS

-  PVP **Merck MilliporeSigma (Sigma-Aldrich)**
-  20 mM EDTA **Merck MilliporeSigma (Sigma-Aldrich)**
-  100 mM Tris-HCl **Merck MilliporeSigma (Sigma-Aldrich)**
-  1.5 M NaCl **Merck MilliporeSigma (Sigma-Aldrich)**
-  2% CTAB **Merck MilliporeSigma (Sigma-Aldrich)**
-  1% β -mercaptoethanol **BBI Biotech**
-  Chloroform: Isoamyl alcohol **Sangon Biotech**
-  isopropanol **Sangon Biotech**
-  sodium acetate **Qiagen**
-  ethanol **BBI Biotech**
-  TE buffer **Thermo Fisher Scientific**

STEP MATERIALS

-  PVP **Merck MilliporeSigma (Sigma-Aldrich)**
-  20 mM EDTA **Merck MilliporeSigma (Sigma-Aldrich)**
-  100 mM Tris-HCl **Merck MilliporeSigma (Sigma-Aldrich)**
-  1.5 M NaCl **Merck MilliporeSigma (Sigma-Aldrich)**
-  2% CTAB **Merck MilliporeSigma (Sigma-Aldrich)**
-  1% β -mercaptoethanol **BBI Biotech**
-  Chloroform: Isoamyl alcohol **Sangon Biotech**
-  isopropanol **Sangon Biotech**
-  sodium acetate **Qiagen**
-  ethanol **BBI Biotech**
-  TE buffer **Thermo Fisher Scientific**
-  PVP **Merck MilliporeSigma (Sigma-Aldrich)**
-  20 mM EDTA **Merck MilliporeSigma (Sigma-Aldrich)**
-  100 mM Tris-HCl **Merck MilliporeSigma (Sigma-Aldrich)**
-  1.5 M NaCl **Merck MilliporeSigma (Sigma-Aldrich)**
-  2% CTAB **Merck MilliporeSigma (Sigma-Aldrich)**
-  1% β -mercaptoethanol **BBI Biotech**



✕ Chloroform: Isoamyl alcohol **Sangon Biotech**

✕ isopropanol **Sangon Biotech**

✕ sodium acetate **Qiagen**

✕ ethanol **BBI Biotech**

✕ TE buffer **Thermo Fisher Scientific**

Protocol materials

✕ 100 mM Tris-HCl **Merck MilliporeSigma (Sigma-Aldrich)**

✕ 1% β -mercaptoethanol **BBI Biotech**

✕ isopropanol **Sangon Biotech**

✕ ethanol **BBI Biotech**

✕ 2% CTAB **Merck MilliporeSigma (Sigma-Aldrich)**

✕ 20 mM EDTA **Merck MilliporeSigma (Sigma-Aldrich)**

✕ TE buffer **Thermo Fisher Scientific**

✕ PVP **Merck MilliporeSigma (Sigma-Aldrich)**

✕ 1.5 M NaCl **Merck MilliporeSigma (Sigma-Aldrich)**

✕ sodium acetate **Qiagen**

✕ ethanol **BBI Biotech**


✕ TE buffer **Thermo Fisher Scientific**

✕ isopropanol **Sangon Biotech**

✕ Chloroform: Isoamyl alcohol **Sangon Biotech**

✕ Chloroform: Isoamyl alcohol **Sangon Biotech**

- 1 Grind -80°C stored leaves (100-200 mg) to fine powder using liquid nitrogen in the presence of 50 mg PVP (PolyVinyl Pyrrolidone, Mr 10,000). Note: It's very important that the tissue doesn't thaw.



 **PVP Merck MilliporeSigma (Sigma-Aldrich)**

- 2 Add 1.5 mL of freshly prepared 2x CTAB buffer and incubate at 60°C for 45 min with intermittent shaking every 10 minutes.

20 mM EDTA (PH 8.0, 0.5M)	40 ml
100 mM Tris-HCl (PH 8.0, 1M)	100 ml
1.5 M NaCl	87.6 g
2% CTAB (w/v)	20 g
1% β -mercaptoethanol	10 ml (add before use)

Then add distilled water to make it up to 1000 ml.

00:45:00

 20 mM EDTA **Merck MilliporeSigma (Sigma-Aldrich)** 100 mM Tris-HCl **Merck MilliporeSigma (Sigma-Aldrich)** 1.5 M NaCl **Merck MilliporeSigma (Sigma-Aldrich)** 2% CTAB **Merck MilliporeSigma (Sigma-Aldrich)** 1% β -mercaptoethanol **BBI Biotech**

- 3 Centrifuge at 12,000 rpm for 15 min at room temperature.

00:15:00

- 4 Carefully transfer the aqueous phase into a new tube. Note: Use wide-bore tips for transferring the aqueous phase to avoid mechanical damage to DNA.

- 5 Add double volume of Chloroform: Isoamyl alcohol (24 : 1) under fume hood, and invert gently 15 to 20 times and centrifuge at 12,500 rpm for 15 min. Note: If the aqueous layer appears translucent, repeat the step until the solution is transparent. Be careful to avoid transferring any chloroform.

00:15:00

 Chloroform: Isoamyl alcohol **Sangon Biotech**

- 6 Add double volume of chilled isopropanol followed by 1/3rd volume of 3 M sodium acetate (pH 5.2) and keep at -20°C for one hour to precipitate the DNA. Note: The longer the chilled incubation, the more the precipitation.

01:00:00

 isopropanol **Sangon Biotech** sodium acetate **Qiagen**

- 7 Centrifuge at 12,000 rpm for 15 min and discard the supernatant.



00:15:00

- 8 To the pellet, add 70% chilled ethanol and spool out the pellet carefully and centrifuge again at 12,000 rpm for 15 min, and repeat the step again.

00:15:00

 ethanol **BBI Biotech**

- 9 Discard the supernatant and vacuum dry or air dry the pellet at room temperature. Note: Make sure that there is no residual ethanol, this is very critical especially if the DNA is to be used directly for PCR. Overdrying should also be avoided as it makes the pellet difficult to resuspend.
- 10 Add 40-50 μ L TE buffer (10 mM Tris-HCl, 1 mM EDTA, pH 8), and quantify the DNA using Qubit flurometer/Nanodrop, and run an 0.8% Agarose gel electrophoresis.

 TE buffer **Thermo Fisher Scientific**

- 11 Store at -20°C/-40°C till further use.