

Jan 28, 2022

DNA extraction from colonial tunicates

DOI

dx.doi.org/10.17504/protocols.io.b33sqqne



Marta Wawrzyniak¹, Simon Blanchoud¹

¹University of Fribourg

Blanchoud lab, UNIFR



Marta Wawrzyniak

University of Fribourg

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.b33sqqne>

Protocol Citation: Marta Wawrzyniak, Simon Blanchoud 2022. DNA extraction from colonial tunicates. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.b33sqqne>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 20, 2022

Last Modified: January 28, 2022

Protocol Integer ID: 57170

Keywords: DNA extraction, colonial tunicates, ascidians, dna extraction from colonial tunicate, hotphenol dna extraction protocol, botrylloides diegensis, dna extraction, dna, colonial tunicate

Abstract

This protocol has been successfully used with *Botrylloides diegensis* and was adapted to our needs based on the HotPhenol DNA extraction protocol.

Guidelines

Change gloves frequently, particularly as the protocol progresses from crude extracts to more purified materials. Use sterile tubes. Perform all steps on ice and use RNase-free and DNase-free water unless otherwise stated.

Materials

Heat bath setup at 70 C (If working with fresh samples: glass beads 0.1mm and eppendorf thermal shaker)
phenol pH 8 (4 C)

Ph-Ch-IA: phenol:chloroform:isoamyl alcohol (25:24:1, best prepared fresh)

Ch-IA: chloroform:isoamyl alcohol (24:1)

SDS-Lysis buffer : 10mL Lysis buffer, 4mL SDS 10%

Lysis buffer 50mL: 12.3g 3M Sodium acetate (pH 5.2), 7.3g 0.5M EDTA , Nuclease-free water

12.3g 3M sodium acetate in 50mL Nuclease-free water

80% Ethanol

ultra pure water

Troubleshooting



- 1 This protocol was developed to extract both RNA and DNA in parallel (See [RNA extraction from colonial tunicates](#) , steps 1-7) using the same samples. However it could be run directly on fresh samples (steps 1.1-1.4).
- 1.1 Clean the slide from which you will take the colony of your interest. See [Cleaning colonial ascidians](#).
- 1.2 Isolate a cleaned colony composed of approx. 20 zooids.
- 1.3 Transfer to a tube and spin at maximum speed for 00:02:00 . 2m
- 1.4 Remove the excess water.


Cell Lysis

4m 30s

- 2 Add 500 μ L of Ph-Ch-IA solution and 350 μ L of SDS-Lysis buffer.
- 3 Heat the tube at 70 $^{\circ}$ C for 00:02:00 . 2m
- 3.1 If working with fresh samples, add glass beads 0.1mm to the tube and shake on eppendorf thermal shaker at 900rpm.
- 4 Mix by vortexing at maximum speed for 00:00:45 . 45s
- 5 Cool for 00:01:00 on ice. 1m
- 6 Mix again by vortexing at maximum speed for 00:00:45 . 45s
- 7 Heat the tubes at 70 $^{\circ}$ C for 00:10:00 , mix regularly by inversion. 10m




7.1 If working with fresh samples, shake the tubes on eppendorf thermal shaker at 900rpm.

8 Mix again by vortexing at maximum speed for  00:00:45 .

45s

9 Cool for  00:01:00 on ice.



1m

10 Mix again by vortexing at maximum speed for  00:00:45 .


45s


DNA extraction

3m

11 Centrifuge at  Room temperature at maximum speed for  00:03:00 .


3m

12 Transfer  400 μL of the upper aqueous phase to a new tube.


13 Add  400 μL of Ph-Cl-IA solution.


14 Shake the tube by inversion for  00:00:30 .

30s

15 Centrifuge at maximum speed for  00:03:00 .


3m

16 Transfer  300 μL of the upper aqueous phase to a new tube.

17 Add  300 μL of Ph-Cl-IA solution.


18 Shake the tube by inversion for  00:00:30 .


30s

19 Centrifuge at maximum speed for  00:03:00 .

3m



20 Transfer  200 μL of the upper aqueous phase to a new tube.

21 Add  200 μL of CI-IA solution.

22 Shake the tube by inversion for  00:00:30 .

30s


23 Centrifuge at maximum speed for  00:03:00 .


3m

24 Transfer aqueous phase to a new tube.

DNA precipitation

3h 30m



25 Add 2 volumes of [M] 100 % volume Ethanol (typically  300-400 μL).

26 Add 0.1 volume of [M] 3 Molarity (M) sodium acetate (typically  15-20 μL).

27 Mix by inversion.


28 Incubate at  -20 $^{\circ}\text{C}$ for  03:00:00 .

3h

29 Centrifuge at maximum speed for  00:20:00 at  4 $^{\circ}\text{C}$.

20m

30 Discard the supernatant.

31 Add  450 μL of cold [M] 80 % volume ethanol.



- 32 Centrifuge at maximum speed for 00:05:00 at Room temperature .
- 33 Discard the supernatant.
- 34 Add 200 μL of cold 80 % volume ethanol.
- 35 Centrifuge at maximum speed for 00:05:00 at Room temperature .
- 36 Discard the supernatant.
- 37 Resuspend the pellet in ultra pure water (typically 20-100 μL).
- 38 Measure the DNA concentration using the NanoDrop.
- 39 Store at -20 °C for short storage or at -80 °C for long storage.

5m

5m