

Nov 09, 2022

DNA extraction from fecal samples

DOI

dx.doi.org/10.17504/protocols.io.5jyl8jnorg2w/v1

Anique Ahmad¹, Arya Gautam¹, Tsedenia Denekew¹, Aashish Jha¹

¹New York University, Abu Dhabi

Aashish Jha: jhaar@nyu.edu



Genetic Heritage Group NYU Abu Dhabi

NYU Abu Dhabi

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5jyl8jnorg2w/v1>

Protocol Citation: Anique Ahmad, Arya Gautam, Tsedenia Denekew, Aashish Jha 2022. DNA extraction from fecal samples. [protocols.io](https://dx.doi.org/10.17504/protocols.io.5jyl8jnorg2w/v1) <https://dx.doi.org/10.17504/protocols.io.5jyl8jnorg2w/v1>

Manuscript citation:

Anique Ahmad, Arya Gautam, Tsedenia Denekew, Aashish Jha (2022). DNA extraction from fecal samples. [protocols.io](https://dx.doi.org/10.17504/protocols.io.5jyl8jnorg2w/v1) dx.doi.org/10.17504/protocols.io.5jyl8jnorg2w/v1

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 16, 2022

Last Modified: November 09, 2022

Protocol Integer ID: 70128

Keywords: extraction from fecal samples dna extraction, fecal samples dna extraction, dna extraction, fecal sample, extraction, dna, sample, dna extraction from fecal samples dna extraction

Abstract

DNA Extraction from fecal samples

Guidelines

Sample Type	Maximum Input
Feces	200 mg
Soil	250 mg
Liquid Samples ¹ and Swab Collections ²	250 µl
Cells (isotonic buffer, e.g. PBS)	50-100 mg (wet weight) (10 ⁹ bacterial and 10 ⁸ yeast cells)
Samples in DNA/RNA Shield™, ³	≤ 1 ml

Table 1. Sample type and maximum input for DNA extraction using this protocol

Materials

ZymoBIOMICS DNA Extraction Kit

Fecal samples

Homogenizing pestle

Liquid nitrogen or dry ice

Ice box with ice

1.5 ml microcentrifuge tubes

Vortex

Omni Bead Ruptor Elite bead beater

Centrifuge

Troubleshooting



Before start

- Label the following:
 - ZR BashingBead™ Lysis Tubes (0.1 & 0.5 mm)
 - ZymoSpin™ III-F Filter in Collection Tube
 - Collection Tube
 - ZymoSpin™ II-CR Column in Collection Tube
 - ZymoSpin™ III-HRC Filter in Collection Tube
 - 2 Sets of 1.5 ml microcentrifuge tubes (not provided with kit)
- Spin the labelled ZR BashingBead™ Lysis Tubes (0.1 & 0.5 mm) for 10 seconds in a mini-centrifuge to ensure that the beads have settled at the bottom
- Include 1 control per batch and assign its position randomly.



- 1 Transfer fecal sample to *ZR BashingBead™ Lysis Tube (0.1 & 0.5 mm)* and add ZymoBIOMICS™ Lysis Solution:
 - 1.1
 - **Flash frozen fecal samples** will be solid and *an aliquot* ( 15 mg) of feces is expected in a 2mL container (Table 1). Add  250 µL of ZymoBIOMICS™ Lysis Solution to the fecal sample and homogenize evenly using the homogenizing pestle while submerged in liquid nitrogen (or on dry ice). Once the sample is homogenized, add another  500 µL of ZymoBIOMICS™ Lysis Solution such that the final volume is  750 µL . Mix well by vortexing.
 - Spin the sample container for  00:00:10 in a minicentrifuge.
 - Transfer the entire  750 µL of the sample mix into the *ZR BashingBead™ Lysis Tube (0.1 & 0.5 mm)* and proceed to Step 2.
 - 1.2 **Fecal samples collected using Zymo DNA/RNA Shield™** will exist as a solution (Table 1). Mix the sample by shaking well or vortexing. Add up to  250 µL of such sample to a *ZR BashingBead™ Lysis Tube (0.1 & 0.5 mm)*. Adjust final volume to  1 mL with ZymoBIOMICS™ Lysis Solution. Cap the tube tightly and proceed to Step 2.
 - 1.3 **Fecal samples collected using other stabilization kits** may exist as a solution (Table 1). Mix the sample by shaking well or vortexing. Add  250 µL of such sample to a *ZR BashingBead™ Lysis Tube (0.1 & 0.5 mm)*. Add  750 µL ZymoBIOMICS™ Lysis Solution to the tube. Cap tightly and proceed to Step 2.
- 2 Secure in a Omni Bead Ruptor Elite bead beater fitted with a 2 ml tube holder assembly and process at max speed (30 m/s) for  00:05:00 . Rest for  00:05:00 .
- 3 Centrifuge the *ZR BashingBead™ Lysis Tubes (0.1 & 0.5 mm)* in a microcentrifuge at ≥  10000 x g, 00:01:00 .
- 4 Transfer up to  600 µL supernatant to the *Zymo-Spin™ III-F Filter* in the *labelled Collection Tube* and centrifuge at  8.000 x g, 00:01:00 . Discard the Zymo-Spin™ III-F Filter.
- 5 Add  600 µL of ZymoBIOMICS™ DNA Binding Buffer to the filtrate in the *labelled Collection Tube* from Step 4. Mix well.
- 5.1 Repeat so that the final volume of ZymoBIOMICS™ DNA Binding Buffer added is  1200 µL .

10s

10m

1m

1m



- 6 Transfer  800 μL of the mixture from Step 5 to a *Zymo-Spin™ IICR Column* in a Collection Tube and centrifuge at  10.000 x g, 00:01:00 . 1m
- 6.1 Discard the flow through from the Collection Tube and repeat step 6.
- 7 Add  400 μL ZymoBIOMICS™ DNA Wash Buffer 1 to the *Zymo-Spin™ IICR Column* in a new *Collection Tube* and centrifuge at  10.000 x g, 00:01:00 . Discard the flow-through. 1m
- 8 Add  700 μL ZymoBIOMICS™ DNA Wash Buffer 2 to the *Zymo-Spin™ IICR Column* in a Collection Tube and centrifuge at  10.000 x g, 00:01:00 . Discard the flow-through. 1m
- 9 Repeat with  200 μL ZymoBIOMICS™ DNA Wash Buffer 2 to the *Zymo-Spin™ IICR Column* in a Collection Tube and centrifuge at  10.000 x g, 00:01:00 . 1m
- 10 Transfer the *Zymo-Spin™ IICR Column* to a clean 1.5 ml microcentrifuge tube and add  100 μL (50 μL minimum) ZymoBIOMICS™ DNase/RNase Free Water directly to the column matrix and incubate for  00:05:00 . Centrifuge at  10.000 x g, 00:01:00 to elute the DNA 6m
- 11 Place a *Zymo-Spin™ III-HRC Filter* in a new *Collection Tube* and add  600 μL ZymoBIOMICS™ HRC Prep Solution. Centrifuge at  8.000 x g, 00:03:00 . 3m
- 12 Transfer the eluted DNA (Step 10) to a prepared *Zymo-Spin™ III-HRC Filter* in a clean 1.5 ml microcentrifuge tube and centrifuge at exactly  16.000 x g, 00:03:00 . 3m
- 13 The filtered DNA is now suitable for PCR and other downstream applications. Eluted DNA should be frozen (–30 to –15°C or –90 to –65°C)