


Feb 11, 2020

DNA extraction from whole blood using simple salting out procedure

DOI

dx.doi.org/10.17504/protocols.io.bbzqip5w



Schuele Lab¹, Faria Zafar¹

¹Department Pathology, Stanford University School of Medicine

Neurodegeneration Method Development Community
Tech. support email: ndcn-help@chanzuckerberg.com



Schuele Lab

Stanford University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bbzqip5w

Protocol Citation: Schuele Lab, Faria Zafar 2020. DNA extraction from whole blood using simple salting out procedure. protocols.io <https://dx.doi.org/10.17504/protocols.io.bbzqip5w>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol in our group and it is working.

Created: January 31, 2020

Last Modified: February 11, 2020

Protocol Integer ID: 32528

Keywords: DNA, DNA extraction, blood, whole blood, human blood, salting out, nuclei lysis, molecular weight dna extraction from whole blood, dna extraction, molecular weight dna extraction, extraction, dna, genetic analysis, range sequencing, procedure the purpose, procedure,

Abstract

The purpose of this protocol is high-molecular weight DNA extraction from whole blood for genetic analyses, including no amplification long-range sequencing.

Citation

S. A. Miller, D. D. Dykes, H. F. Polesky (1988)
. A simple salting out procedure for extracting DNA from human nucleated cells. Nucleic Acids Research.

[10.1093/nar/16.3.1215](https://doi.org/10.1093/nar/16.3.1215)

[LINK](#)

Attachments



DNA extraction from ...

212KB



Guidelines

Storage Vial labeling:

3 lines, 6pt bold font; 1.5" label

a) Subject ID#: four digit SUSL ID (e.g. SUSL-2345), Sample type: *WBC*

b) Date of blood draw: MM/DD/YY

c) Study: [Name Study], Initials



Materials

MATERIALS

- ☒ Sodium chloride **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S3014**
- ☒ Ammonium chloride (≥ 99.5 %) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A9434**
- ☒ Potassium bicarbonate (≥ 99.5 %) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #90339**
- ☒ EDTA 500 mM Solution pH 8.0 ULTROL® Grade **Merck Millipore (EMD Millipore) Catalog #324504**
- ☒ Trizma® hydrochloride / Tris-HCl **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T5941**
- ☒ Ethyl Alcohol 200 Proof (GR ACS) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #EX0276**
- ☒ 2-Propanol (99.5 %) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #278475**
- ☒ Proteinase K from Tritirachium album **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P2308**

Note

If using different vendors, choose molecular biology grade reagents.

Buffers and Solution Recipes

10X Red blood cell (RBC) lysis buffer (1l)

NH ₄ Cl	82.91 g
KHC O ₃	10.01 g
EDTA 0.5 M	20 ml
complete to 1000 ml with autoclaved dH ₂ O	

* dilute to 1x before use with autoclaved dH₂O

1 M Tris-HCl, pH 8 (used to make nucleolysis buffer and TE)

Tris - HCl	15.7 g
Bring volume to ~ 95ml with autoclaved ddH ₂ O	
measure & adjust pH to	

8
Bring final volume up to 100 ml with autoclaved ddH ₂ O

10X Nucleus lysis buffer (100ml)

NaCl	23.38 g
Tris-HCl (1 M, pH = 8)	10 ml
EDTA 0.5 M	4 ml
complete to 100 ml with autoclaved ddH ₂ O	

* dilute to 1x before use with autoclaved ddH₂O

20 % SDS (used end of Day 1 for DNA extraction to poke holes in cell membrane)

SDS	20 g
autoclaved ddH ₂ O	80 ml

5 M NaCl

NaCl	146.1 g
Bring volume to 500 ml with autoclaved dH ₂ O	

70 % EtOH (used Day 2 of DNA extraction to "clean" DNA)

100 % EtOH	35 ml
ddH ₂ O	15 ml

TE (used at end of DNA extraction to re-suspend DNA and prepare for storage)



1 M Tris- HCl	1 ml
0.5 M EDTA	0.2 ml
Bring volume to 100 ml with autoclaved ddH ₂ O	

Proteinase K reconstitution:

Solutions can be prepared in 25 millimolar (mM) Tris-HCl buffer , pH 8.0 , containing 1 millimolar (mM) Calcium chloride . Amounts for typical usage is 50 µg/ml – 200 µg/ml . Solutions are stable at pH 8 at 4 °C . It is stable in a broad range of environments: pH, buffer salts, detergents (SDS), and temperature. In the presence of 0.1 – 0.5 % SDS, proteinase K retains activity and will digest a variety of proteins and nucleases in DNA preparations without compromising the integrity of the isolated DNA.

Materials

- 10 ml EDTA collection tube (K₂EDTA, 5.4mg; 16 × 100mm; 10mL; Pink)
- 50 ml conical tube (Fisher catalog # 12565268)
- 1.5 ml clear Eppendorf tube (Fisher catalog # 5408129)
- 1000 µl filter pipette tips
- 5.0 ml serological pipette (Fisher catalog # 13-678-11D)
- 50 ml serological pipette (Fisher catalog # 13-678-11F)
- Gloves
- Bleach (e.g. Clorox, as 20 % fresh solution)

Equipment

- Centrifuge: Sorvall Legend XTR, or SorvallIT 6000B, H10000B
- Centrifuge: Eppendorf 5417R
- Rotator

Safety warnings

⚠ Please see SDS (Safety Data Sheet) for hazards and safety warnings.

Before start

Prepare buffers and solutions (for recipes see 'Materials').



- 1 If you use FROZEN blood, start with step-case 'FROZEN blood'. If FRESH blood is used, select step-case 'FRESH blood'.

STEP CASE

FROZEN blood

37 steps

- 2 In the morning take out tubes of blood from the -80 °C freezer. Record ID numbers and put in fridge (4 °C) to thaw until later in the afternoon.
- 3 Once blood has thawed, invert several times or place on rotator for a couple of minutes.
- 4 Decant the blood sample into 50 ml conical tube.
- 5 Rinse blood tube with 1 mL of 1x RBC lysis buffer and add to 50ml conical tube.
- 6 Add 1x RBC lysis buffer up to 50 mL .
- 7 Place on shaker at Room temperature for 00:30:00 .

Note

Note: Blood solution becomes transparent.

- 8 Centrifuge 00:15:00 at 2000 rpm at Room temperature .
- 9 Discard supernatant in container with 20 % bleach .

Note

Note: The pellet contains the nuclei.

**Note**

If direct processing with Nuclei lysis, day 1, then continue **after this step** with Nuclei lysis (**step 12**).

- 10 Add  700 μ L 1x RBC lysis buffer and close the tube tightly. Pipette up and down to resuspend the pellet (cleans up the nuclei pellet from additional RBCs and allows for storage in small tubes). 
- 11 Transfer to 1.5 ml **pre-labeled** clear Eppendorf tube.
- 12 Centrifuge  00:05:00 at  2000 rpm at  Room temperature . 
- 13 Using a 1000 μ L pipette tip, discard the supernatant in container with  20 % bleach .
- 14 Pellets can be frozen at this point at  -80 $^{\circ}$ C .

Nuclei lysis (Day 1, afternoon)

- 15 Add  3 mL Nucleus Lysis Buffer and close tube tightly. 
- 16 Pipette up and down to resuspend pellet. 
- 17 Add  300 μ L of  10 % SDS and  70 μ L Proteinase K ( 10 mg/ml). 
- 18 Mix by gently swirling tube. 
- 19 Incubate  Overnight in water bath at  55 $^{\circ}$ C . 



DNA extraction (Day 2, morning)

20 Add  1 mL of  5 Molarity (M) NaCl to the tubes. 

21 Close caps tightly and **shake vigorously**.

Note

Note: Very important to **shake thoroughly** at this step to precipitate proteins.

22 Centrifuge for  00:20:00 at  3000 rpm at  Room temperature . 

23 While centrifuging, label 2× 15ml conical tubes per sample.

24 Transfer supernatant to a 15 ml conical tube and then repeat centrifugation for

 00:20:00 at  3000 rpm

Note

Note: The supernatant contains the DNA.

25 Transfer supernatant to next 15 ml conical tube.

26 Add equal volume of **ice-cold isopropanol** (usually about  4 mL). 

27 Invert tubes several times and DNA will precipitate.

Note

Note: You should be able to see **white strings** (this is the DNA!).



28 Centrifuge for  00:05:00 at  3000 rpm .



29 Discard supernatant.

30 Add  1 mL  70 % ethanol to each tube.



31 Centrifuge for  00:05:00 at  3000 rpm .



32 Carefully pipette off *ethanol* to not dislodge the DNA pellet.



33 Leave DNA to dry in **uncapped tubes** overnight.

Dissolve DNA

34 Add  250 μ L TE and allow DNA to dissolve in 15 ml conical tube stored in fridge ( 4 °C).



35 Once **completely dissolved**, transfer to labeled sterile Eppendorf tube.

36 Label top with **subject ID**. Label side with **date of extraction** and **concentration**.

37 Spec all samples using **Nanodrop**. Dilute samples to stock at  250 ng/ μ l –  350 ng/ μ l and to working dilution at  10 ng/ μ l .



38 Transfer samples into barcode tubes and log into DNA bank.



Citations

S. A. Miller, D. D. Dykes, H. F. Polesky. A simple salting out procedure for extracting DNA from human nucleated cells

[10.1093/nar/16.3.1215](https://doi.org/10.1093/nar/16.3.1215)