



Feb 06, 2023

Version 1

DNA extraction - Zooplankton - 96 wells V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlk4x7wg5r/v1>

coline.royaux^{1,2,3}, Nicolas Rabet^{1,2,3}, Céline Bonillo^{1,2,3}

¹Sorbonne Université; ²Muséum National d'Histoire Naturelle; ³UMR BOREA



Coline Royaux

Université Pierre et Marie Curie (Paris VI), Muséum National...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlk4x7wg5r/v1>

Protocol Citation: coline.royaux , Nicolas Rabet, Céline Bonillo 2023. DNA extraction - Zooplankton - 96 wells. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.j8nlk4x7wg5r/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 15, 2021



Last Modified: February 06, 2023

Protocol Integer ID: 54125

Keywords: parts of zooplanktonic freshwater crustacean, zooplanktonic freshwater crustacean, dna extraction, branchiopoda, zooplankton, copepoda, dna, extraction, new caledonia

Abstract

This protocol was used to extract DNA from whole or parts of zooplanktonic freshwater crustaceans (Copepoda, Branchiopoda, ...) from New Caledonia.



- 1 Prepare your 96-well extraction plate with one individual per well. Alternate genus in the wells to detect eventual contamination between wells. 1d
 - 1.1 Collect one individual from a sample 1m
 - 1.2 Note its genus and determine its sex with a binocular microscope 1m
 - 1.3 For big individuals (more than 5 mm), dissect a few legs and put it in the well. Be careful not to damage the rest of the body and put it in a tagged Eppendorf tube. For little individuals (less than 5 mm), put the whole body. 3m
- Note

If necessary, use alcohol to get the biological material to fall at the bottom of the well
- 1.4 When all 96 wells are filled, the biological material has to dry to go to lysis 12h
- Note

If necessary, use a micropipette to empty an excess of alcohol in the well
- Safety information

Make sure the plate is closed when you want to transport it elsewhere
- 2 Prepare the lysis 15m
 - 2.1 Mix  18 mL T1 buffer and  2.5 mL K proteinase in a Multi-Channel Reservoir and distribute  200 µL of the mix with a multimicropipette in each well 10m



Note

🧪 180 μ L T1 buffer and 🧪 25 μ L K proteinase in each well

- 2.2 Close your extraction plate with a heated aluminium foil and an adhesive plastic film 3m
- 3 Put your plate in a proofer at 🌡 56 °C ⌚ Overnight (6h or more) to lyse the tissues 6h
- 4 Perform the DNA extraction with a DNA extraction robot
- 4.1 Remove the adhesive film and aluminium foil from the plate and put it in the robot
- 4.2 It will deposit 🧪 200 μ L BQ1 buffer and 🧪 200 μ L ethanol in each wells and mix it
- 4.3 Then, 🧪 600 μ L of the wells content (lysate, BQ1, ethanol) are transfered on the tissue binding plate. Reagents excess are emptied in a waste container.
- 4.4 The tissue binding plate is then dried by a ⌚ 00:05:00 aspiration to bind DNA to the silica membrane of the binding plate 5m
- 4.5 The silica membrane is then washed with 🧪 600 μ L BW buffer and twice with 🧪 900 μ L B5 buffer per well. Each wash is intercalated by a ⌚ 00:05:00 aspiration dry. 5m
- 4.6 The waste container is then removed from under the binding plate which is dried again by a ⌚ 00:10:00 aspiration 10m
- 4.7 An empty extraction plate is placed under the binding plate to retrieve the genomic DNA from it



- 4.8 DNA is eluted from the tissue binding plate with  100 μ L BE buffer in each well and is collected in the new extraction plate underneath
- 4.9 After a  00:03:00 rest, the binding plate is dried for  00:02:00 and the elution is repeated with  100 μ L BE buffer
- 4.10 Retrieve the new extraction plate containing the genomic DNA and discard the rest

5m