

May 28, 2024

DNA Extraction: Zymo Research Quick-DNA Fecal/Soil Microbe Midiprep Kit (Cat #: D6110)

DOI

dx.doi.org/10.17504/protocols.io.3byl49b2ogo5/v1

Jenna Brown¹, Mirayana Marcelino Barros², Tiago Pereira³, Alejandro De Santiago Perez², Hunter Powell⁴

¹University of Georgia- Holly Bik; ²University of Georgia; ³Department of Marine Sciences, University of Georgia; ⁴UGA

Bik Lab UGA



Jenna Brown

University of Georgia- Holly Bik

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl49b2ogo5/v1>

Protocol Citation: Jenna Brown, Mirayana Marcelino Barros, Tiago Pereira, Alejandro De Santiago Perez, Hunter Powell 2024.

DNA Extraction: Zymo Research Quick-DNA Fecal/Soil Microbe Midiprep Kit (Cat #: D6110). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.3byl49b2ogo5/v1>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 02, 2024

Last Modified: May 28, 2024

Protocol Integer ID: 99117

Keywords: dna extraction of soil sample, soil microbe midiprep kit, dna extraction, dna fecal, zymo research quick, d6110, soil sample, dna, bik lab, kit, zymo, extraction

Funders Acknowledgements:

Holly Bik

Abstract

DNA Extraction: Zymo Research Quick-DNA Fecal/Soil Microbe Midiprep Kit (Cat #: D6110); kit and protocol used for DNA extraction of soil samples in the Bik Lab.

Troubleshooting



DNA Extraction: Zymo Research Quick-DNA Fecal/Soil Microbe Midiprep Kit (Cat #: D6110)

15m

- 1 Add **beta-mercaptoethanol** (not included in kit) to the **Genomic Lysis Buffer (Cat #: D3004-1-100)** to a final dilution of **0.5 % (v/v)** (v/v) *i.e.*, 2.5 ml per 500 ml for optimal performance.
- 2 Add **5 g** of soil sample to the bead/filter of a **ZR Bashing Bead Lysis/Filtration Tube (50 mL w/ 0.5 mm Beads; Cat #: S6010)**.
- 3 Add **6 mL** of **BashingBead Buffer (Cat #: D6001-3-150)** to the sample. Cap tube tightly and secure with parafilm.
 - 3.1 Place tubes in a bead beater (Bullet Blender) and process the samples for 10 minutes at speed 8.
Tip: Store tubes sideways until they go into the bead beater to prevent liquid from prematurely going through filter.
Tip: Select "Zymo MidiPrep" pre-loaded protocol on Bullet Blender in Bik Lab.
 - 3.2 Centrifuge the **ZR BashingBead Lysis/Filtration Tube** in a centrifuge at **5000 rpm, 00:05:00**.
 - 3.3 Remove bead/filter chamber from the top of the **ZR BashingBead Lysis/Filtration Tube** and transfer supernatant from the bottom of the tube (using a 10 mL pipette) to a new 50 mL tube.
- 4 Add **18 mL** of **Genomic Lysis Buffer** to the supernatant. Mix well.
Tip: Centrifuge for a few seconds if there is too much foam.
 - 4.1 Filter the entire mixture using a **Zymo-Spin V-E Column/Zymo-Midi Filter (Cat #: C1021-25)** mounted on a vacuum manifold, at ≥ 600 mm Hg.
 - 4.2 Disconnect the **Zymo-Spin V-E Column/Zymo-Midi Filter** and transfer the **Zymo-Spin V-E Column** to a **Collection Tube (Cat #: C1001-50)**. Spin the column at **10000 rpm, 00:01:00** in a microcentrifuge.
















5m



1m





- 5 Add  300 μL of **DNA Pre-Wash Buffer (Cat #: D3004-5-15)** to the column and spin at  10000 rpm, 00:01:00 . Discard the flow through.
Tip: Add reagent by inserting the pipette tip at an angle while pressing it against the wall of the column and without touching the filter to avoid overflow. 1m 
- 6 Add  400 μL of **g-DNA Wash Buffer (Cat #: D3004-2-50)** to the column and centrifuge at  10000 rpm, 00:01:00 . Discard the flow through. 1m 
- 7 **REPEAT WASH STEP:** Add  400 μL of **g-DNA Wash Buffer** to the column and centrifuge at  10000 rpm, 00:01:00 . Discard the flow through. 1m
- 8 Transfer the **Zymo-Spin V-E Column** to a 1.5 mL microcentrifuge tube. Add  150 μL of **DNA Elution Buffer (Cat #: D3004-4-16)** directly to the column matrix. Wait for 1 minute and then centrifuge at  10000 rpm, 00:01:00 to elute the DNA. 1m 
- 9 Place the **Zymo-Spin III-HCR Filter** in a clean **Collection Tube**. Add  600 μL of **Prep Solution (Cat #: D6035-1-30)**. Centrifuge at  8000 rpm, 00:03:00 . 3m 
- 9.1 Place the **Zymo-Spin III-HCR Filter** in a clean 1.5 mL microcentrifuge tube. Transfer the eluted DNA to the prepared **Zymo-Spin III-HCR Filter** and centrifuge at exactly  16000 rpm, 00:03:00 . The filtered DNA is now suitable for PCR and other downstream applications. 3m 