


Feb 12, 2024

DNA isolation from cattle tissues: blood, semen or any kind of tissues, including ear biopsies

DOI

<https://dx.doi.org/10.17504/protocols.io.kqdg3xojeg25/v1>



Cecile CG Grohs¹

¹Université Paris-Saclay, INRAE, AgroParisTech, GABI, 78350, Jouy-en-Josas, France



Cecile Grohs

Université Paris-Saclay, INRAE, AgroParisTech, GABI, 78350, ...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kqdg3xojeg25/v1>

Protocol Citation: Cecile CG Grohs 2024. DNA isolation from cattle tissues: blood, semen or any kind of tissues, including ear biopsies. [protocols.io https://dx.doi.org/10.17504/protocols.io.kqdg3xojeg25/v1](https://dx.doi.org/10.17504/protocols.io.kqdg3xojeg25/v1)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: January 30, 2024

Last Modified: February 22, 2024

Protocol Integer ID: 94392

Keywords: dna isolation from cattle tissue, isolate dna, routine method for isolate dna, dna isolation, cattle tissue, available frozen semen straw, dna, ear biopsy, genomic dna, including ear biopsy, removal of protein, precipitation of genomic dna, biopsy, washing of sample, semen, other tissue, tissue, different kinds of tissue

Abstract

Here we describe a routine method for isolate DNA from different kinds of tissues: commercially available frozen semen straws, blood, ear biopsies, or other tissues.

This protocol is based on a salting-out method and uses several commercially available solutions.

It consists of several steps: washing of samples, lysis, removal of proteins and precipitation of genomic DNA. This protocol was used to isolate hundreds of samples for years.

Guidelines

For recovering DNA from blood, use K3-EDTA tubes.

Salting out is a good method to obtain DNA, as it avoids using phenol/chloroform and allows to recover DNA from different qualities of blood.



Materials

☒ Isopropanol

☒ 70% ethanol

☒ Puregene Tissue Kit **Qiagen Catalog #158063**

☒ Puregene blood kit **Qiagen Catalog #158106**

☒ Proteinase K, 2mL **Qiagen Catalog #19131**

☒ Tris(2-carboxyethyl)phosphine hydrochloride solution **Merck MilliporeSigma (Sigma-Aldrich) Catalog #646547-10X1ML**

☒ Buffer RLT **Qiagen Catalog #79216**

☒ 1X PBS (Phosphate-buffered saline)

☒ DNA LoBind Tubes 2.0 mL **Eppendorf Catalog #30108078**

2mL tubes

2 X50mL tubes for each blood sample

Centrifuge for 2mL and 50mL tubes, blood tubes

3D heating rocker

Safety warnings

! See Safety Data Sheets for warnings and safety hazards.

Before start

As we use commercial sperm straws to perform our extractions, we do not always know the composition of these straws, the quantity of material contained, the nature of the diluents and preservatives used. This is why it is sometimes necessary to use several straws to obtain enough material for sequencing. It is also sometimes wise to perform several washes (see step 3) to eliminate contaminants from diluents and preservatives.



Prepare reagents

10m

1

For semen straws, Immediately before use, prepare a mix containing **RLT** buffer (Qiagen) and **TCEP** [Tris(2-carboxyethyl)phosphine hydrochloride] to a final volume of 500µL per sample as follow:

- 450 µL RLT
- 50 µL TCEP

Prepare samples

10m

2 Use DNA LoBind tubes at this stage (strongly recommended for sperm, not mandatory for other tissues).

2.1 For semen:

- Empty the 200 µL Straw in a 2 mL tube by cutting the two ends of the straw
- Rinse the straw with 200 µL 1X **PBS** at Room temperature

For ear biopsy:

- Open the device
- Transfer the biopsy in a 2 mL tube
- Remove extra preservative liquid and the plastic ball

For blood:

- Pour the 5 mL blood in a 50 mL tube

For other tissues:

- Cut a 150 mg piece and put it in a 2 mL tube

2.2 Wash step

30m

For semen:

- Add 800 µL more **PBS** (up to 1 mL 1X **PBS**)
- Pellet 1000 x g at Room temperature during 00:05:00
- Discard the supernatant





Second wash is optional (no significant impact observed)

- Re-suspend in  1 mL 1X **PBS**
- Pellet  1000 x g at  Room temperature during  00:05:00
- Discard the supernatant

For blood:

- Add  15 mL (i.e. three times the initial volume) **RBC** reagent from Puregene blood kit
- Shake vigorously and incubate  00:05:00 at  Room temperature
- Pellet white blood cells  3000 x g, 10°C, 00:15:00

Lysis

3h 30m

3 *Continue with Qiagen Puregene kit adapted as follow*

5m

Step one

For semen:

- Add  500 µL of **RLT/TCEP** mix
- Vortex by pulsing at max speed
- Incubate  00:05:00 on ice
- Add  500 µL **CLS**

For tissues, including ear biopsies:

- Add  800 µL **CLS**
- Add  60 µL **Proteinase K**

For blood:

- Add  5 mL **CLS** (i.e. the initial blood volume)

4 **Step 2**

- Mix by inversion (about 25 inversions)
- Incubate from  01:00:00 to  03:00:00 , depending on the lysis process,
 55 °C
- on a rotating shaker at  200 rpm

4h



Check that no undigested parts remain after lysis



Protein precipitation

30m 30s

5 For semen and tissues:

30m 30s

- Add  200 μL of **Protein Precipitation Solution**
- Vortex  00:00:15 and incubate  00:05:00 on ice
- Pellet  16000 x g, 4°C, 00:05:00

For blood:

- Add  1.65 mL of **Protein Precipitation Solution** (for 5mL blood)
- Vortex  00:00:15 and incubate  00:05:00 on ice
- Pellet  3000 x g, 10°C, 00:15:00

DNA precipitation

6 ▪ Transfert the supernatant to a new tube containing

5m

-  600 μL of Isopropanol (semen) -  2 mL tube
-  800 μL of Isopropanol (tissues) -  2 mL tube
-  5 mL of Isopropanol (blood) -  50 mL tube

- Carrefully invert the tube 25-50X times to form the pellet
- Incubate  00:05:00  Room temperature

For semen and tissues:

- Centrifuge as previously
- Discard supernatant

For blood:

- Pick up the pellet with a lab pipet and transfer it into a  2 mL tube
- Discard supernatant left

Wash DNA

20m

7 ▪ Add  600 μL 70% Ethanol

20m

- Centrifuge  16000 rpm, 4°C, 00:05:00
- Discard supernatant





- Allow the ethanol to evaporate, without drying the pellet  00:15:00

DNA resuspension

20m

- 8
- Add **TE** or **EB** buffer, depending on the subsequent analysis and usual procedure
- We recommend  50 μL for semen or tissue pellets, and  100 μL to  150 μL for pellets from 5mL blood
- Do not hesitate to heat for  01:00:00 at  50 $^{\circ}\text{C}$ on gentle agitation  100 rpm

1h

- 9
- DNA quantity control
- Measure the O.D. with a Nanodrop® device to obtain the DNA concentration
 - Dilute with extra **TE/EB** if you choose to standardize concentrations
 - Measure until the concentration is as expected

1h

DNA quality control

- Load the DNA on a 1% agarose gel in 0.5X TBE with Ethidium Bromide
- Perform an electrophoresis  01:00:00 at 70 V

Protocol references

See protocole for DNA isolation from cattle semen for long read sequencing at DOI dx.doi.org/10.17504/protocols.io.j8nlkw1qwl5r/v1