

Dec 19, 2022

## DNA Isolation from Snake Skin Shed

DOI

[dx.doi.org/10.17504/protocols.io.bp2l699ezlqe/v1](https://dx.doi.org/10.17504/protocols.io.bp2l699ezlqe/v1)



Alexis Lindsey<sup>1</sup>, Tonia S Schwartz<sup>1</sup>

<sup>1</sup>Auburn University

Schwartz Lab of Ecologi...



Alexis Lindsey

Auburn University

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l699ezlqe/v1>

**Protocol Citation:** Alexis Lindsey, Tonia S Schwartz 2022. DNA Isolation from Snake Skin Shed. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bp2l699ezlqe/v1>



**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** November 07, 2022

**Last Modified:** December 19, 2022

**Protocol Integer ID:** 72424

**Keywords:** snake, skin shed, DNA isolation, snake shed, shed, dna isolation from snake skin shed purpose, quality dna from shed reptile skin, shed dna isolation, dna isolation, shed reptile skin, snake skin shed purpose, protocol for skin, louisiana pine snake breeding project, dna, quality dna, skin, isolation

## Abstract

Purpose:

This protocol was developed for the Memphis Zoo's Louisiana Pine Snake Breeding Project. The protocol for skin shed DNA isolation was adapted from Fetzner (1999).

The time estimates assumes you are processing 24 samples and you are well practiced.

References: James W Fetzner (1999) Extracting High-Quality DNA from Shed Reptile Skins: A Simplified Method. BioTechniques 26:6



## Materials

### Equipment

- Sterile razor blades or scissors
- Dissection boards for cutting up the skin sheds

### Consumables

- Filtered micropipette tips (p1000, p200)
- 1.5-mL microcentrifuge tubes (VWR Catalog Number 76332-068)
- 2-mL microcentrifuge tubes ( VWR Catalog Number 20170-170)
- 15 ml conical tube for making 70% Ethanol.
- latex or nitrile gloves

### Reagents

- RNase AWAY, which also degrades DNA (Molecular BioProducts Catalog Number 7002)
- Proteinase K (20 mg/mL) (IBI Science Product Number IB05406)
- Lysis Buffer (10mM Tris-base, 10mM EDTA, 2% sodium dodecyl sulfate (SDS), pH 8.0)\*
- TE Buffer (10mM Tris-base, 0.1 mM EDTA, pH 8.0; Growcells.com Catalog No: MRGF-4240)
- 70% ethanol (200 proof ETOH, Deacon Labs Product Number 3916EA in Molecular Grade Water (see below)
- 5M aqueous solution ammonium acetate (ThermoFisher Scientific (Alfa Aesar) Product Number J60688)
- Isopropanol (ThermoFisher Scientific (Acros Organics) Catalog Number AC327272500)
- Molecular Grade Water (QualityBiological Catalog Number 351-029-131)

\*Lysis Buffer Recipe: For a final volume of 100 mL:

95 mL molecular-grade water

1 mL 1M Tris-base (VWR Product Number E199)

2 grams SDS (VWR Product Number 0227)

2 mL 0.5M EDTA (VWR Code E177)

## Troubleshooting



## Set Up

1h 10m

- 1 Turn on shaking incubator and set to 55 °C  
Obtain Ice  
Print list of samples 10m
- 2 Clean the razors, scissors, and dissection boards with RNAase away; rinse with molecular-grade water. 30m
- 3 Make fresh 70% ETOH. For example, to make 15 ml in a 15 ml tube use a 10 ml sterile pipette to take 10.5 mL of [M] 100 % volume (200 proof) ETOH, and another 10 ml sterile pipette to add 4.5 mL of molecular grade water. 10m
- 4 Set out and label 1.5 mL microcentrifuge tube for each sample to use in Step 6 (Digestion).  
Add 900 µL **Cell Lysis Buffer**  
Add 10 µL of **proteinase K**( 20 mg/mL). 10m
- 5 Set out and label a 2 mL tube for each sample to use in Step 13 (DNA precipitation).  
Add containing 900 µL **isopropanol**. 10m

## Isolation: Lysis

4h 5m

- 6 Cut up 1 in<sup>2</sup> piece of shed into smaller pieces with a sterile razor blade or scissors (change gloves and utensils between sheds).  
Put shed pieces in the labeled 1.5 mL microcentrifuge tube containing the Cell Lysis Buffer and Proteinase K (prepared in Step 4). 1h
- 7 Place in shaking incubator 300 rpm, 55°C 3 hours to overnight Vortex occasionally for the first few hours. 3h
- 8 Cool samples to room temperature and vortex. 5m

## Isolation: Precipitate Proteins

25m

- 9 Add 550 µL **5M ammonium acetate** to each tube and vortex for 10 seconds. 5m



- 10 Place samples on ice for Room temperature for 10 minutes. 10m
- 11 Centrifuge samples 17,000 x g, 00:05:00 to pellet protein and debris. 5m
- 12 Draw off as much supernatant as possible with a filtered tip into put in to a new 1.5 mL labeled tube. 2m
- 13 Centrifuge the supernatant at second time at 20,000 x g, 00:03:00 to pellet any residual protein and debris. 3m


## Isolation: Precipitate DNA

1h 13m

- 14 With a filtered tip, transfer supernatant from the second spin into the prepared 2mL tubes containing the isopropanol (prepared in Step 4). 2m
  - 15 Mix the supernatant with the isopropanol by inverting 50 times. If there is a lot of DNA you can see the strands condensing at this step (looks like thin white threads). 2m
  - 16 Place each tube into the centrifuge with the hinge facing out so the DNA pellet forms on that side of the tube. Centrifuge samples at 16,000 x g, 00:02:00 to pellet the DNA. 2m
- Expected result


There should be a small, white/clearish pellet of DNA present towards the bottom of the tube on the side of the hinge.
- 17 Pour off isopropanol into waste container. 1m
  - 18 Wash the DNA pellet by adding 500  $\mu$ L of **70% ethanol**. 2m
  - 19 Centrifuge at 16,000 x g, 00:02:00 . The DNA pellet should still be visible. 2m



- 20 Without dislodging the DNA pellet, carefully pour the supernatant out into waste container with as little movement as possible. 1m
- 21 Centrifuge the tubes again for  16,000 x g, 00:01:00 , and use a 10 µl tip to remove the residual ethanol. This will make the next step go faster. 1m \*
- 22 Invert tubes on a paper towel, with the tops open, until ethanol has completely evaporated. 1h


## Resuspension

1h 1m

- 23 Once the ethanol has evaporated (but the DNA pellet is not over dry), resuspend samples in  50 µL of **TE buffer**. 1m

### Expected result

This produces on average this produced 500 ng of DNA per µl. This volume of TE buffer can be increased or lowered if you want a higher concentration or if you started with a smaller amount of shed.

- 24 Sit at  300 rpm, 37°C, 01:00:00 or leave in 4°C overnight to fully resuspend the DNA. 1h

## Check Quality

- 25 Run 5 µL of resuspended DNA on 1% agarose gel to visual the quality and estimate quality. DNA can be quantified with the Nanodrop. For sensitive procedures (DNA sequencing library preparation) we recommend using the Agilent TapeStation or BioAnalyzer.