

Sep 28, 2023

Version 2

DNA Isolation (Gel Clean-up) V.2

DOI

dx.doi.org/10.17504/protocols.io.81wgbxnk3lpk/v2

NUS iGEM¹

¹National University of Singapore



NUS iGEM

National University of Singapore

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.81wgbxnk3lpk/v2>

Protocol Citation: NUS iGEM 2023. DNA Isolation (Gel Clean-up). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.81wgbxnk3lpk/v2> Version created by **NUS iGEM**

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 27, 2023



Last Modified: September 28, 2023

Protocol Integer ID: 88520

Keywords: DNA Isolation, DNA, Gel, Gel Electrophoresis, Buffer QG, Buffer PE, dna fragments from the agarose gel, dna isolation, dna fragment, gel electrophoresis, agarose gel, dna, gel clean, isolation, singapore igem team, nus

Abstract

2023 NUS-Singapore iGEM team followed this protocol to isolate the DNA fragments from the agarose gel after the gel electrophoresis.

Protocol materials

⊗ Buffer QG **Qiagen Catalog #19063**

⊗ Buffer PE **Qiagen Catalog #19065**

⊗ Buffer PE **Qiagen Catalog #19065**

⊗ Buffer PE **Qiagen Catalog #19065**

Troubleshooting

Safety warnings










- ❗
 - Proper lab PPE must be worn at all times.
 - When using the LED transilluminator, wear eyewear that is designed to block blue light to protect the eyes.











- 1 Prepare and label an Eppendorf tube.
- 2 Place the agarose gel onto the LED Transilluminator (blue light) to observe the DNA band(s).

Safety information

- Wear protective eyewear to protect the eyes from blue light.
- Turn off the LED transilluminator immediately when it is not in use.

- 3 Cut out the target DNA band from the agarose gel.
- 4 Put the gel piece into the Eppendorf tube.
- 5 Add  450 μL of  Buffer QG **Qiagen Catalog #19063** into the Eppendorf tube.
- 6 Heat the Eppendorf tube at  55 $^{\circ}\text{C}$ for  00:20:00 in the Thermo-Shaker. 20m
- 7 Add  150 μL of 100% isopropanol (IPA) into the Eppendorf tube and shake the tube to mix the solution well.
- 8 Transfer the whole solution into a QIAquick Spin Column (purple tube with a maximum volume of  750 μL).
- 9 Centrifuge it for  13 rpm, 00:01:00 . 1m
- 10 Discard the flow-through and place the QIAquick column back into the same tube.
- 11 Add  700 μL of  Buffer PE **Qiagen Catalog #19065** into the QIAquick column.



- 12 Centrifuge it for  13 rpm, 00:01:00 .
- 13 Discard the flow-through and place the QIAquick column back into the same tube.
- 14 Add  700 μL of  Buffer PE **Qiagen Catalog #19065** again into the QIAquick column.
- 15 Centrifuge it for  13 rpm, 00:01:00 .
- 16 Discard the flow-through and place the QIAquick column back into the same tube.
- 17 Centrifuge the emptied QIAquick column at  13 rpm, 00:01:00 to remove residual  Buffer PE **Qiagen Catalog #19065** . 1m
- 18 Transfer the QIAquick column into the newly labelled Eppendorf tube.
- 19 Add  30 μL of DI water into the QIAquick column.
- 20 Centrifuge the tube at  13 rpm, 00:01:00 , ensuring that the direction of the Eppendorf tube's cap is the same as the direction of spinning to avoid breaking. 1m
- 21 Discard the QIAquick column, the solution left in the Eppendorf tube contains the DNA fragment of interest.
- 22 Use the Nanodrop to measure and record the purity and concentration of the DNA fragment.

Equipment

NanoDrop™ One/OneC Microvolume UV-Vis Spectrophotometer		NAME
UV-Vis Spectrophotometer		TYPE
Thermo Scientific		BRAND
ND-ONE-W		SKU

23 Keep the isolated DNA fragment in the 🌡️ Room temperature .