



Nov 05, 2024

DNA Library Prep with MGIEasy FS DNA Library Prep Set

DOI

dx.doi.org/10.17504/protocols.io.261gerpkol47/v1

Qianyue Ji¹, Guoqiang Mai², Xiaohan Wang¹, Shanshan Liu³, Mo Han¹

¹BGI Research; ²China National GeneBank; ³MGI Tech



Qianyue Ji

BGI

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261gerpkol47/v1>

Protocol Citation: Qianyue Ji, Guoqiang Mai, Xiaohan Wang, Shanshan Liu, Mo Han 2024. DNA Library Prep with MGIEasy FS DNA Library Prep Set. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.261gerpkol47/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 27, 2024



Last Modified: November 05, 2024

Protocol Integer ID: 110978

Keywords: fs dna library prep set, dna library prep, wgs libraries construction for mgi, metagenomic, ng genomic dna, mgi, genomic dna, library prep set, library prep, mgi tech, sequencing platform series, fidelity pcr enzyme, dna, wgs libraries construction

Abstract

The MGIEasy FS DNA Library Prep Set(1000006987,MGI Tech) is specifically designed for WGS libraries construction for MGI high-throughput sequencing platform series. This library prep set is optimized to convert 5-400 ng genomic DNA into a customized library. This set uses advanced Adapter Ligation technology and High-fidelity PCR Enzymes to significantly increase library yield and conversion rate. This library prep set is applicable for samples from all common animals, plants, fungus, bacteria, etc., including human, mice, rice, Arabidopsis, yeast, E. coli, Metagenomics.

Materials

Materials

MGIEasy Dual Barcode Circularization Kit (Cat. No.: 1000020570)

Qubit ssDNA Assay Kit

MGIEasy FS DNA Library Prep Set

Supplies:

Pipette tips (assorted volumes)

5ml round-bottom polystyrene tubes, PCR tube, centrifuge tube

Equipment:

Qubit Fluorometer

Vortex

Micropipettes, various volumes

Microcentrifuge

Thermocycler

Magnetic rack

Troubleshooting

Fragmentation

1 Preparation

1.1 Mix the reagents before using and store the remaining reagents immediately after use.

Reagent	Requirements
Dilution Buffer	User-supplied:place at room temperature (RT).
Frag Buffer II	Thaw at RT, mix well, centrifuge briefly, and place on ice.
Frag Enzyme II	Keep on ice.

preparing the reagents

2 Fragmentation

2.1 Add 45 µL of genomic DNA to a new 0.2mL PCR tube. Add dilution buffer to make a total volume of 45 µL if the fragmentation gDNA volume not enough. Place the tube(s) On ice .

Reagent	Volume
gDNA	X µL
Dilution buffer	45 - X µL
Total	45 µL

Normalization of gDNA

2.2 Set and run program as table below. The thermocycle will perform the first step reaction described in table below and be kept at 4 °C .

Temperature	Time
70 °C Heated lid	On



Temperature	Time
4°C	Hold
30°C	8 min
65°C	15 min
4°C	Hold

Fragmentation reaction conditions

2.3 Mix the Frag Enzyme II by inverting 10 times and flicking the bottom gently, ensure that no residual reagent is left at the bottom. Centrifuge briefly and place On ice .

2.4 According to the desired reaction number, prepare the fragmentation mixture in a 0.2mL PCR tube On ice . Mix it well by vortexing 3 times, centrifuge briefly, and place On ice .

Reagent	Volume per reaction
Frag Buffer II	10µL
Frag Enzyme II	5 µL
Total	15 µL

Fragmentation mixture

2.5 Add 15 µL of fragmentation mixture to each sample tube. Mix by vortexing 3 times, centrifuge briefly, and place On ice .

2.6 Make sure the thermocycler has cooled to 4 °C . Place the PCR tube(s) into the thermocycler, and skip the 4°C Hold step to start the reaction at 30 °C .


















2.7 When the program is completed, centrifuge the PCR tube(s) briefly and place On ice .

3 Double size selection.

The following sample use a 36 µL 1st bead selection and a amount 12 µL 2nd bead selection to target a 330bp size fragment from fragmented DNA(60 µL).



- 3.1 Transfer all fragmented DNA to a new 1.5mL centrifuge tube. Add TE buffer for a final volume of  60 μL .
- 3.2 Add  36 μL of DNA Clean Beads to each sample tube. Mix with a vortexer until all beads are suspended.
- 3.3 Incubate at  Room temperature for  00:05:00 .
- 3.4 Centrifuge the tube(s) briefly and place on the magnetic rack for  00:05:00 until the liquid is clear. Then, carefully transfer the supernatant to a new 1.5mL centrifuge tube. Retain the supernatant and discard the beads.
- 3.5 Add  12 μL of DNA Clean Beads to the centrifuge tube with  94 μL supernatant. Mix with a vortexer until all beads are suspended.
- 3.6 Incubate at  Room temperature for  00:05:00 .
- 3.7 Centrifuge the tube(s) briefly and place on the magnetic rack for  00:05:00 until the liquid is clear. Carefully remove and discard all the supernatant.
- 3.8 While keeping the centrifuge tube on the magnetic rack, add  200 μL of 80% ethanol to each tube to wash the beads and tube wall. Wait for 30 sec. Carefully remove and discard the supernatant.
- 3.9 Repeat step 3.8 and try to remove all of the liquid from the tube.
- 3.10 Keep the tube(s) on the magnetic rack. Open the tube cap and air-dry the beads at  Room temperature until no wetness or glossiness is visible on the beads' surface. There should be no visible cracking on the surface of the beads.
- 3.11 Remove the tube(s) from the magnetic rack and add  43 μL of TE to elute the DNA. Mix with a vortexer until all beads are suspended.
- 3.12 Incubate the tube(s) at  Room temperature for  00:05:00 .

5m



5m



5m





5m



5m





- 3.13 Centrifuge the tube(s) briefly and place on the magnetic rack for  00:05:00 until the liquid is clear. Carefully transfer  41 μL of supernatant to a new 1.5mL centrifuge tube.
- 3.14 Quantify the purified size selection products with Qubit dsDNA Assay Kit.

5m



End repair


4 End repair

- 4.1 Transfer  100 ng DNA sample to a new 0.2mL PCR tube and add TE buffer for a total volume of  40 μL .

- 4.2 Prepare the End repair Mixture  On ice.

Reagent	Volume per reaction
ERAT Buffer	7.1 μL
ERAT Enzyme Mix	2.9 μL
Total	10 μL

End repair Mixture

- 4.3 Transfer  10 μL of the mixture to the 0.2mL PCR tube from step 3. Vortex 3 times and centrifuge briefly.

- 4.4 Place the 0.2mL PCR tube into the thermocycler and run the program.






Temperature	Time
70 °C Heated lid	On
37 °C	30 min
65°C	15 min
4°C	Hold

End repair program

- 4.5 When the program completed, briefly centrifuge to collect the solution at the bottom of the tube.



Adapter ligation

5 Adapter ligation

- 5.1 Add  5 µL of adapter to the corresponding sample tube. Vortex it 3 times (3 sec each), centrifuge briefly, and place  On ice .
- 5.2 Prepare the adapter ligation mixture. Vortex it 6 times (3 sec each), centrifuge briefly, and place  On ice .

	Reagent	Volume per reaction
	Ligation Buffer	23.4 µL
	DNA Ligase	1.6 µL
	Total	25 µL

adapter ligation mixture

- 5.3 Slowly pipette  25 µL of adapter ligation mixture to each sample tube. Vortex it 6 times (3 sec each), centrifuge briefly, and place  On ice .

The adapter ligation mixture is highly viscous. Pipette slowly and carefully.

- 5.4 Place the PCR tube(s) into the thermocycler. Run the program with the following conditions.

	Temperature	Time
	30 °C Heated lid	On
	23 °C	30 min
	4°C	Hold

adapter ligation program





- 5.5 When the program is completed, centrifuge the PCR tube(s) briefly and place On ice .
- 5.6 Add 20 μL TE buffer, for a total volume of 100 μL and transfer to a new 1.5mL centrifuge tube.

Cleanup of Adapter-ligated DNA

- 6 Cleanup of Adapter-ligated DNA
- 6.1 Add 50 μL of DNA Clean Beads to each sample tube. Mix with a vortexer until all beads are suspended.
- 6.2 Incubate the sample tube at Room temperature for 00:05:00 .
- 6.3 Centrifuge the sample tube briefly and place on the magnetic rack for 00:05:00 until the liquid is clear. Carefully remove and discard all the supernatant.
- 6.4 Keep the tube on the magnetic rack, add 200 μL of 80% ethanol to each tube to wash the beads and tube wall. Wait for 30 sec. Carefully remove and discard the supernatant.
- 6.5 Repeat step 6.4 once. Try to remove all liquid from the tube. If some liquid remains on the tube wall, centrifuge the tube briefly and place it on the magnetic rack for separation. Remove all liquid by using a low-volume pipette.
- 6.6 Keep the tube on the magnetic rack. Open the tube cap and air-dry the beads at Room temperature until no wetness or glossiness is visible on the beads' surface. There should be no visible cracking on the surface of the beads.
- 6.7 Remove the tube from the magnetic rack and add 21 μL of TE buffer to elute the DNA. Mix with a vortexer until all beads are suspended.
- 6.8 Incubate the tube(s) at Room temperature for 00:05:00 .

5m





5m



5m





- 6.9 Centrifuge the tube briefly and place on the magnetic rack for  00:05:00 until the liquid is clear. Carefully transfer  19 μL of supernatant to a new 0.2 mL PCR tube.

5m


PCR Amplification

7 PCR Amplification

- 7.1 Prepare the PCR mixture in a 0.2 mL PCR tube  On ice . Vortex it 3 times (3 sec each), centrifuge briefly, and place  On ice .

Reagent	Volume per reaction
PCR Enzyme Mix	25 μL
PCR Primer Mix	6 μL
Total	31 μL

PCR mixture

- 7.2 Add  31 μL of PCR mixture to each sample tube. Vortex 3 times (3 sec each) and centrifuge briefly to collect the solution at the bottom of the tube.

- 7.3 Place the PCR tube into the thermocycler. Run the program with the following conditions.



Temperature	Time	Cycles
Temperature	Time	Cycles
105 °C Heated lid	on	
95 °C	3min	1
98 °C	20s	8
60 °C	15s	
72 °C	30s	
72 °C	10min	1
4 °C	Hold	1













PCR program

- 7.4 When the program is completed, centrifuge the tube(s) briefly and transfer to a new 1.5mL centrifuge tube.

Cleanup of PCR Product

8 Cleanup of PCR Product

- 8.1 Add  50 μL of DNA Clean Beads to each sample tube from 7.4. Mix with a vortexer until all beads are suspended.
- 8.2 Incubate the sample tube at  Room temperature for  00:05:00 .
- 8.3 Centrifuge the sample tube briefly and place on the magnetic rack for  00:05:00 until the liquid is clear. Carefully remove and discard all the supernatant.
- 8.4 Keep the tube on the magnetic rack, add  200 μL of 80% ethanol to each tube to wash the beads and tube wall. Wait for  00:00:30 . Carefully remove and discard the supernatant.
- 8.5 Repeat step 8.4 once. Try to remove all liquid from the tube. If some liquid remains on the tube wall, centrifuge the tube briefly and place it on the magnetic rack for separation. Remove all liquid by using a low-volume pipette.
- 8.6 Keep the tube on the magnetic rack. Open the tube cap and air-dry the beads at  Room temperature until no wetness or glossiness is visible on the beads' surface. There should be no visible cracking on the surface of the beads.
- 8.7 Remove the tube from the magnetic rack and add  32 μL of TE buffer to elute the DNA. Mix with a vortexer until all beads are suspended.
- 8.8 Incubate the tube(s) at  Room temperature for  00:05:00 .

5m





5m



30s

5m



8.9 Centrifuge the tube briefly and place on the magnetic rack for  00:05:00 until the liquid is clear. Carefully transfer  30 μL of supernatant to a new 1.5mL centrifuge tube.

5m






8.10 Quality control of PCR product
Quantify the dsDNA with Qubit dsDNA-HS Assay Kit. Yield for PCR products: >1 pmol.
Assess the size range of purified PCR products with Agilent 2100 Bioanalyzer.

Denaturation and single-stranded circularization

20m 30s

9 Denaturation


9.1 Based on the PCR products concentration, add  300 ng of PCR products into a new 0.2 mL PCR tube. If the volume is less than  48 μL , add TE Buffer to make a total volume  48 μL .

9.2 Place the PCR tube(s) into the thermocycler. Run the program with the following conditions.




	Temperature	Time
	100 °C Heated lid	On
	95 °C	3 min
	4 °C	10 min

Denaturation reaction conditions (Volume: 48 μL)

9.3 After the reaction, centrifuge the tube briefly and place  On ice.





10 Single-stranded circularization

10.1 According to the desired reaction number, prepare the circularization reaction mixture in a new 0.2 mL PCR tube on ice. Vortex it 3 times (3 sec each), centrifuge briefly, and place  On ice.

Reagent	Volume per reaction
Dual Barcode Splint Buffer	11.6 µL
DNA Rapid Ligase	0.5 µL
Total	12.1 µL

Circularization reaction mixture


10.2 Add  12.1 µL of circularization reaction mixture to each sample tube . Vortex it 3 times (3 sec each), centrifuge briefly, and place  On ice .

10.3 Place the PCR tube(s) into the thermocycler. Run the program with the following conditions.







Temperature	Time
42 °C Heated lid	On
37 °C	10 min
4 °C	Hold

Single-stranded DNA circularization reaction conditions (Volume: 60 µL)

10.4 When the program is completed, place the PCR tube(s)  On ice , centrifuge briefly, and immediately proceed to the next step.

11 Digestion

11.1 According to the desired reaction number, prepare the digestion mixture in a 0.2 mL PCR tube  On ice . Vortex it 3 times (3 sec each), centrifuge briefly, and place  On ice .

11.2 Add  4 µL of digestion mixture to each sample tube. Vortex it 3 times (3 sec each), centrifuge briefly, and then place  On ice .



11.3 Place the PCR tube(s) into the thermocycler. Run the program with the following conditions.














	Temperature	Time
	42 °C Heated lid	On
	37 °C	10 min
	4 °C	Hold



Digestion reaction conditions (Volume: 64 µL)

- 11.4 When the program is completed, centrifuge the tube briefly and immediately add  7.5 µL of Digestion Stop Buffer to each sample tube. Vortex it 3 times (3 sec each), centrifuge briefly, and place  On ice .

12 Cleanup of digestion product



- 12.1 Mix the DNA Clean Beads thoroughly. Add  130 µL of DNA Clean Beads to each sample tube. Mix with a vortexer until all beads are suspended.
- 12.2 Incubate at  Room temperature for  00:05:00 . 5m 
- 12.3 Centrifuge the tube(s) briefly and place on the magnetic rack for  00:05:00 until the liquid is clear. Carefully remove and discard the supernatant. 5m
- 12.4 While keeping the tube(s) on the magnetic rack, add  160 µL of 80% ethanol to each tube to wash the beads and tube wall. Wait for  00:00:30 . Carefully remove and discard the supernatant. 30s
- 12.5 Repeat step 12.4. Try to remove all liquid from the tube. If some liquid remains on the tube wall, centrifuge the tube briefly and place it on the magnetic rack for separation. Remove all liquid by using a low-volume pipette.
- 12.6 Keep the tube(s) on the magnetic rack. Open the tube cap and air-dry the beads at  Room temperature until no wetness or glossiness is visible on the beads' surface. There should be no visible cracking on the surface of the beads.
- 12.7 Remove the tube(s) from the magnetic rack and add  25 µL of TE Buffer to elute the DNA. Mix with a vortexer until all beads are suspended.



12.8 Incubate at  Room temperature for  00:05:00 .

5m



12.9 Centrifuge the tube briefly and place on the magnetic rack for  00:05:00 until the liquid is clear. Carefully transfer  24 μL of supernatant to a new 1.5 mL centrifuge tube.

5m



13 QC of digestion product

13.1 Quantify the ssCir with Qubit ssDNA Assay Kit. The final Enzymatic Digestion products (ssDNA, ng) / input products of PCR (dsDNA, 300 ng) should be $\geq 7\%$.