

Sep 09, 2020

DNA Quantification -- CHEM 584



Forked from [DNA Concentration Measurement \(Protocol for Thermo Scientific NanoDrop™ 1000 Spectrophotometer\)](#).

DOI

dx.doi.org/10.17504/protocols.io.bk4xkyxn

Alba Balletbó¹, Ken Christensen²

¹Wageningen University; ²Brigham Young University



Ken Christensen

Brigham Young University

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bk4xkyxn

Protocol Citation: Alba Balletbó, Ken Christensen 2020. DNA Quantification -- CHEM 584. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bk4xkyxn>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: September 09, 2020

Last Modified: September 09, 2020

Protocol Integer ID: 41847

Abstract

Protocol adapted from the NanoDrop Spectrophotometer User's Manual.





Materials

TE Buffer or MQ Water

DNA Sample



- 1 Start a new experiment using the Nanodrop touchscreen.
- 2 Click on the corresponding application module (e.g., dsDNA).
- 3 Open the sampling arm and load a blank sample (e.g., TE Buffer, MQ Water, etc.).
 2 μL
- 4 Close the sampling arm on the machine to cover the blank sample.
- 5 Click "OK" to read the blank if the Nanodrop is not already in automatic mode.
- 6 Open the sampling arm and clean the blank off the upper and lower pedestals using a Kim Wipe.
- 7 Load your DNA sample and close the sampling arm.
 2 μL
- 8 Click "Measure sample" if the Nanodrop is not in automatic mode.
- 9 Take a picture of the screen that shows your results and the corresponding spectra for your records.
- 10 Clean the nucleic acid sample off of the upper and lower pedestals using a Kim Wipe.
- 11 Repeat for all samples.
- 12 Close program.