



Nov 02, 2018

DNA recovery from filter paper

DOI

dx.doi.org/10.17504/protocols.io.u8iezue

Diego Bonatto¹

¹Centro de Biotecnologia da UFRGS



Diego Bonatto

Centro de Biotecnologia da UFRGS

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.u8iezue

External link: <https://help.addgene.org/hc/en-us/articles/206127485-DNA-on-filter-paper-or-tube>

Protocol Citation: Diego Bonatto 2018. DNA recovery from filter paper. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.u8iezue>

Manuscript citation:

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 01, 2018

Last Modified: November 02, 2018

Protocol Integer ID: 17386

Abstract

Protocol to recover plasmids from filter papers - copy from Addgene webpage (<https://help.addgene.org/hc/en-us/articles/206127485-DNA-on-filter-paper-or-tube>)

Materials

MATERIALS

 TE buffer

STEP MATERIALS

 TE buffer Invitrogen - Thermo Fisher





Protocol materials

 TE buffer

 TE buffer Invitrogen - Thermo Fisher

 TE buffer Invitrogen - Thermo Fisher



- 1 To recover the plasmid, use a clean razor blade to cut out one of the circles containing your DNA.
- 2 Immerse the circle in  30 μL of  TE buffer **Invitrogen - Thermo Fisher** and pipette to mix.
- 3 After waiting for at least  00:10:00 , use  2 μL to transform competent bacteria.