


Dec 03, 2019

Domestication of L0 parts for Loop type IIS (Bsal and SapI)



Forked from [Domestication of L0 parts for Loop type IIS \(Bsal and SapI\)](#)



DOI

<https://dx.doi.org/10.17504/protocols.io.92jh8cn>

Eftychis Frangedakis¹, marta tomaselli², Marius Rebmann², Susana Sauret-Gueto²

¹University of Cambridge; ²Plant Sciences, University of Cambridge, OpenPlant

OpenPlant Project



Susana Sauret-Gueto

Plant Sciences, University of Cambridge, OpenPlant

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.92jh8cn>

Protocol Citation: Eftychis Frangedakis, marta tomaselli, Marius Rebmann, Susana Sauret-Gueto 2019. Domestication of L0 parts for Loop type IIS (Bsal and SapI). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.92jh8cn>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 03, 2019

Last Modified: December 03, 2019

Protocol Integer ID: 30507

Keywords: cloning of I0 part, I0 part, loop type ii, puap4

Abstract

Domestication and cloning of L0 parts into pUAP4 for **Loop type IIS assembly**.

Materials


MATERIALS

 Phusion high-fidelity PCR kit **Thermo Scientific Catalog #F553S**

 MinElute Gel Extraction Kit **Catalog #28606**

Domestication Guide.

1



Domestication of internal *BsaI* and *SapI* sites is necessary when a DNA sequence of interest contains a *BsaI* or *SapI* recognition site. For domestication, overlapping PCR can be used. Briefly, two PCR fragments are amplified upstream (fragment 1, with primers 1 and 2) and downstream (fragment 2, with primers 3 and 4) of the sequence to domesticate, which will be part of the reverse primer (primer 2) of fragment 1 and of the forward primer (primer 3) of fragment 2. Both primer 2 and 3 will be specially designed with a single nucleotide mismatch to alter the sequence to domesticate (taking care to not alter amino acid composition if the region is in the coding sequence). Primers will also contain *SapI* recognition sites to allow the two fragments to be ligated together into pUAP4 using *SapI* type IIS assembly. Primer 1 and primer 4 are designed according to the guidelines for cloning of LO parts into pUAP4. Blue arrows: *SapI* recognition site. Blue dashed lines: *SapI* cleavage site. Red arrows: *BsaI* recognition site. *CmR*: chloramphenicol bacterial resistance cassette. *LacZ*: *lacZ* α cassette for blue-white screening of colonies.



PCR and amplified fragments purification

- 2 Amplify the parts using your preferred polymerase (e.g. Phusion, Thermo Scientific™) following the manufacturer instructions.
- 3 Run your PCR product on agarose gel.
- 4 Extract the amplified fragments using a PCR purification kit (e.g. Qiagen MinElute™).

Fragments Assembly

- 5 Set up the assembly of the fragments as in [cloning of L0 parts into pUAP4](#)

[pUAP4]= plasmid length/200

[fragment]= fragment length/100

The formulas above will give you the required concentration in ng/μl.