

Dec 29, 2023

Version 1

DoTA-seq V3.1 V.1

Forked from a private protocol

DOI

dx.doi.org/10.17504/protocols.io.81wgbxx3ylpk/v1

freeman.lan¹

¹UW Madison



freeman.lan

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.81wgbxx3ylpk/v1>

Protocol Citation: freeman.lan 2023. DoTA-seq V3.1. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.81wgbxx3ylpk/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: August 24, 2023

Last Modified: December 29, 2023

Protocol Integer ID: 86909

Keywords: seq target primer, process of dota, sequencing library, multiplex dota, sequencing, seq generating, pcr step, seq, dota, cell suspension, generating large molecular weight primer dimer, single cell, large molecular weight primer dimer, pcr, cell

Abstract

This protocol describes the process of DoTA-seq generating a single cell sequencing library from a cell suspension. This workflow can be performed in two days, with the PCR step happening overnight. Before beginning this workflow make sure to have:

1. The necessary microfluidics devices prepared and ready to go
2. The multiplex DoTA-seq target primers validated to work together without generating large molecular weight primer dimers.

Please read the publication for further details.

Guidelines

Strongly recommend all pre-PCR steps (setting up reagents, washing gels) to be done in a PCR Clean hood. This has two purposes:

1. Reduce PCR contamination of templates which can strongly effect single-cell PCR reactions.
2. Reduce dust contamination of reagents which can clog devices and cause failures.



Materials

- ⊗ ddPCR Supermix for probes (no dUTP) **Bio-Rad Laboratories Catalog #1863024**
- ⊗ MetaPolyzyme **Merck MilliporeSigma (Sigma-Aldrich) Catalog #MAC4L-5MG**
- ⊗ Lysozyme from chicken egg white **Merck MilliporeSigma (Sigma-Aldrich) Catalog #L6876**
- ⊗ HFE 7500 Perfluorinated Oil
- ⊗ Perfluorooctanol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #370533**
- ⊗ TCEP-HCl **Gold Biotechnology Catalog #TCEP**
- ⊗ NN'-Bis(acryloyl)cystamine **Santa Cruz Biotechnology Catalog #sc-215506**
- ⊗ Ammonium persulfate **Catalog #A3678** ⊗ Acrylamide **P212121**
- ⊗ TEMED (Tetramethyl-ethylenediamine) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T9281**
- ⊗ Biorad Evagreen Droplet Oil **Bio-Rad Laboratories Catalog ##1864005**
- ⊗ DNA Clean & Concentrator™-5 **Zymo Research Catalog #D4003**
- ⊗ DNA Clean & Concentrator™-5 **Zymo Research Catalog #D4003**
- ⊗ Axygen® 0.2 mL Maxymum Recovery® Thin Wall PCR Tubes **Corning Catalog #PCR-02-L-C**
- ⊗ NEBNext Library Quant Kit for Illumina - 100 rxns **New England Biolabs Catalog #E7630S**
- ⊗ SYBR Green **Thermo Fisher Scientific**
- ⊗ Proteinase K solution, 20 mg ml – 1 **Ambion Catalog #AM2546**
- ⊗ Lysozyme from chicken egg white **Merck MilliporeSigma (Sigma-Aldrich) Catalog #L6876**
- ⊗ Pre-injection buffer 10mM HEPES pH 7.5 Pluronic 0.1%

















Safety information

Unpolymerized Acrylamide is toxic, handle with care and dispose according to regulations



Protocol materials

- SYBR Green **Thermo Fisher Scientific**
- NN'-Bis(acryloyl)cystamine **Santa Cruz Biotechnology Catalog #sc-215506**
- Ammonium persulfate **Catalog #A3678**
- TEMED (Tetramethyl-ethylenediamine) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T9281**
- MetaPolyzyme **Merck MilliporeSigma (Sigma-Aldrich) Catalog #MAC4L-5MG**
- Biorad Evagreen Droplet Oil **Bio-Rad Laboratories Catalog ##1864005**
- DNA Clean & Concentrator™-5 **Zymo Research Catalog #D4003**
- Proteinase K solution, 20 mg ml – 1 **Ambion Catalog #AM2546**
- Pre-injection buffer 10mM HEPES pH 7.5 Pluronic 0.1%
- TCEP-HCl **Gold Biotechnology Catalog #TCEP**
- Acrylamide **P212121**
- ddPCR Supermix for probes (no dUTP) **Bio-Rad Laboratories Catalog #1863024**
- Axygen® 0.2 mL Maxymum Recovery® Thin Wall PCR Tubes **Corning Catalog #PCR-02-L-C**
- NEBNext Library Quant Kit for Illumina - 100 rxns **New England Biolabs Catalog #E7630S**
- Lysozyme from chicken egg white **Merck MilliporeSigma (Sigma-Aldrich) Catalog #L6876**
- Lysozyme from chicken egg white **Merck MilliporeSigma (Sigma-Aldrich) Catalog #L6876**
- Perfluorooctanol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #370533**
- DNA Clean & Concentrator™-5 **Zymo Research Catalog #D4003**
- HFE 7500 Perfluorinated Oil
- PBS 0.1% Tween20
- PBS 0.1% Tween 20
- SYBR Green **Thermo Fisher Scientific**
- Cellbrite Fix 555 **Biotium Catalog ##30088**
- Acrylamide **P212121**
- NN'-Bis(acryloyl)cystamine **Santa Cruz Biotechnology Catalog #sc-215506**
- Ammonium persulfate **Catalog #A3678**
- Biorad Evagreen Droplet Oil **Bio-Rad Laboratories Catalog ##1864005**
- TEMED (Tetramethyl-ethylenediamine) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T9281**
- Biorad Evagreen Droplet Oil **Bio-Rad Laboratories Catalog ##1864005**
- Perfluorooctanol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #370533**

 Acetone Isopropanol PBS EDTA 1mM Tween20 2% w/v 1X PBS (Phosphate-buffered saline) Lysozyme from chicken egg white **Merck MilliporeSigma (Sigma-Aldrich) Catalog #L6876** MetaPolyzyme **Merck MilliporeSigma (Sigma-Aldrich) Catalog #MAC4L-5MG** 1X PBS (Phosphate-buffered saline) PBS 0.1% Tween20 Proteinase K solution, 20 mg ml – 1 **Ambion Catalog #AM2546** 10% SDS **Bio-Rad Laboratories Catalog #161-0146** TE Buffer (Tris 10mM EDTA 1mM pH 8) TE 2% Tween-20 TE 2% Tween-20 Pre-injection buffer 10mM HEPES pH 7.5 Pluronic 0.1% Pre-injection buffer 10mM HEPES pH 7.5 Pluronic 0.1% ddPCR Supermix for probes (no dUTP) **Bio-Rad Laboratories Catalog #1863024** DTT **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D0632** Biorad Evagreen Droplet Oil **Bio-Rad Laboratories Catalog ##1864005** Axygen® 0.2 mL Maxymum Recovery® Thin Wall PCR Tubes **Corning Catalog #PCR-02-L-C** EDTA (0.5 M), pH 8.0 **Life Technologies Catalog #AM9260G** HFE 7500 Perfluorinated Oil Perfluorooctanol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #370533** TCEP-HCl **Gold Biotechnology Catalog #TCEP** DNA Clean & Concentrator™-5 **Zymo Research Catalog #D4003** SPRIselect reagent kit **Beckman Coulter Catalog #B23317** High Sensitivity D1000 ScreenTape **Agilent Technologies Catalog #5067-5584** NEBNext Library Quant Kit for Illumina - 100 rxns **New England Biolabs Catalog #E7630S** HFE 7500 Perfluorinated Oil

Troubleshooting



Preparing Cells

10m

- 1 Prepare a cell suspension by washing twice in 1mL of PBS 0.1% Tween20 by spinning down at 5000 x g, 00:01:00
- 2 Resuspend cells in 100 μ L PBS 0.1% Tween 20
- 3 Add 1 μ L SYBR Green Thermo Fisher Scientific 10,000X dye to the cells to stain them
- 4 Count cells using a hemacytometer using the SYBR signal, calculate concentration of the cell suspension.
- 5 (Optional) Stain with Cellbrite Fix 555 Biotium Catalog ##30088 to get a cell membrane/wall stain

Preparing Gel

30m

- 6 Make 200 μ L Hydrogel Precursor Solution - Mix together in a tube:
 100 μ L Acrylamide P212121 monomer in water 25 Mass / % volume
 15 μ L
 NN'-Bis(acryloyl)cystamine Santa Cruz Biotechnology Catalog #sc-215506 in Methanol 5 Mass / % volume
 10 μ L Ammonium persulfate Catalog #A3678 10 Mass / % volume
 75 μ L Cell suspension diluted in PBS (**a total of 7e6 cells** to achieve a final concentration of **3.5e7 cells/mL** in the total solution)

Vortex Vigorously to Mix

Generate Gel Droplets

10m

- 7 Prepare and Load the Syringes with the gel sample and 600 μ L
 Biorad Evagreen Droplet Oil Bio-Rad Laboratories Catalog ##1864005 and connect to the microfluidic devices by following this protocol



Protocol

NAME

Loading Syringes to Inject into Microfluidics Device

CREATED BY

freeman.ian

[Preview](#)

- 8 Connect the syringes to the inlets of the **DoTA-seq Step 1** Microfluidics Device and Run the syringe pumps at **500uL/hr** for the gel syringe, and **900uL/hr** for the oil syringe.


- 9 Collect gel droplets for  00:20:00 in a 1mL tube.



20m

Note

Sometimes the initial droplet formation produces polydisperse droplets. In this case, wait 2 min for the bad emulsion to leave the outlet tubing into a waste tube, then begin collecting in the collection tube.

- 10 Make  200 μ L Gel Polymerization Oil - mix together in a tube:

 195 μ L Biorad Evagreen Droplet Oil **Bio-Rad Laboratories Catalog ##1864005** 5 μ L TEMED (Tetramethyl-ethulenediamine) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T9281**

- 11 Add the Gel polymerization oil to the collected droplets, invert slowly 3 times to mix, and Incubate the tube containing droplets at  37 °C for  00:10:00 to complete polymerization of the gel matrix.

10m



Note

You can now look at the emulsion under the microscope using

Countess slides **Thermo Fisher Scientific Catalog #C10228** to determine the encapsulation ratio of your cells. SYBRGreen and CF555 signal should be concordant and correspond to cells.

Breaking out gels from emulsion

12 Pulse spin the emulsion in a centrifuge to close pack the emulsion and drain the oil to the bottom of the tube.

13 Use a pipette to remove the oil from the bottom of the tube, leaving just the emulsion

14 Add 200 μL

1m

Perfluorooctanol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #370533** to break the emulsion

Vortex, then Wait 00:01:00 for the emulsion to break.

15 Pulse spin again and remove the oil in the bottom of the tube with a pipette.

16 Add 1000 μL of HFE 7500 Perfluorinated Oil to the tube.

10s

Mix by inverting 5 times, Wait 00:00:10 then remove with a pipette.

17 Add 1000 μL of Acetone to the tube.

10s

Mix by inverting 5 times, Wait 00:00:10 then remove with a pipette.

The gels should begin to flocculate and dehydrate.

18 Add 1000 μL of Isopropanol




10s

Mix by inverting 5 times, Wait 00:00:10 then remove with a pipette.

The gels should dehydrate and become hard.


Note: Do not wait too long as it could cause the gels to irreversibly aggregate into clumps.



- 19 Resuspend in  1000 μ L  PBS EDTA 1mM Tween20 2% w/v
- The gels can be stored at  4 $^{\circ}$ C for several days without changing DoTA-seq results.

Note







You can now look at the gels under the microscope using

 Countess slides **Thermo Fisher Scientific Catalog #C10228** to determine the encapsulation ratio of your cells.

For **SYBR Staining**, you should first wash a small aliquot in PBS 2% Tween to remove background SYBR+ oil droplets before visualization on the microscope.

You should see some loss in CF555 signal as the acetone and alcohol wash removes some bacterial membranes.

Lysing Bacteria

- 20 Wash gels 3 times in  1000 μ L  1X PBS (Phosphate-buffered saline) **(No Tween)** by centrifugation at  500 x g, 00:00:30 each time 30s
- 21 Make a Enzymatic Lysis Solution by adding:
-  20 mg
-  Lysozyme from chicken egg white **Merck MilliporeSigma (Sigma-Aldrich) Catalog #L6876**
-  100 μ L **[M] 1 mg/mL**
-  MetaPolyzyme **Merck MilliporeSigma (Sigma-Aldrich) Catalog #MAC4L-5MG**
-  900 μ L  1X PBS (Phosphate-buffered saline) **(No Tween)**
- 22 Resuspend the gels in this lysis solution. Incubate at  37 $^{\circ}$ C for  02:00:00 2h
- 23 Wash the gels 3 times in  1000 μ L  PBS 0.1% Tween20
- 24 Make a SDS Lysis solution by adding:
-  20 μ L **[M] 20 mg/mL**
-  Proteinase K solution, 20 mg ml – 1 **Ambion Catalog #AM2546**

 100 μL 10% SDS **Bio-Rad Laboratories Catalog #161-0146** 880 μL

TE Buffer (Tris 10mM EDTA 1mM pH 8)

25 Resuspend the gels in this SDS Lysis solution, incubate at 55 °C for 01:00:00

1h

26 Wash the gels three times in 1000 μL TE 2% Tween-20

Note: Use **2% Tween 20**, not 0.1% Tween

These gels can be stored at 4 °C in TE 2% Tween-20 for several days without impacting DoTA-seq results.

Note

You can now look at the gels under the microscope using

Countess slides **Thermo Fisher Scientific Catalog #C10228** to determine the encapsulation ratio and lysis efficiency of your cells.

It is advised to restain with SYBR and CF555 to get best signal. Lysed cells should exhibit SYBR signal but no CF555 Signal.

Barcoding the Cells

7m

27 Wash the gels three times in 1000 μL

Pre-injection buffer 10mM HEPES pH 7.5 Pluronic 0.1%

28 Resuspend gels in 100 μL

Pre-injection buffer 10mM HEPES pH 7.5 Pluronic 0.1%

29 Load the gels into a syringe following the protocol described in this excellent visual protocol.



Citation

Demaree B, Weisgerber D, Lan F, Abate AR (2018)
. An Ultrahigh-throughput Microfluidic Platform for Single-cell Genome Sequencing..
Journal of visualized experiments : JoVE.

<https://doi.org/10.3791/57598>

LINK

Alternatively, for a simpler version, you can also use a P200 pipette to directly pipette the gel into a syringe backfilled with HFE7500.

- 30 Generate a PCR Master Mix (This mix gives about ~10,000 cells per library - Scale up as required)

🧪 25 μ L

🧬 ddPCR Supermix for probes (no dUTP) **Bio-Rad Laboratories Catalog #1863024**

🧪 0.4 μ L [M] 50 micromolar (μ M) P7 Primer mixed with Barrev primer at

[M] 5 micromolar (μ M)

🧪 0.4 μ L [M] 50 micromolar (μ M) P5 Primer with appropriate I5 index

🧪 0.2 μ L Variable [M] 10 micromolar (μ M) DoTA-seq multiplex primer mix (**10-20nM final concentration per primer**)

🧪 0.2 μ L Variable [M] 10 micromolar (μ M) 16S DoTA-seq primers (**10-20nM final concentration**)

🧪 0.5 μ L Variable [M] 1 picomolar (pM) Freshly diluted from 500pM stock Barcode

Oligo

🧪 0.25 μ L [M] 500 millimolar (mM) Single use aliquots

🧬 DTT **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D0632**

Note

The ratio of 16S to DoTA-seq target primers mix can be varied depending on the relative amplification efficiencies. The best way to determine is to start from equal concentrations, then adjust based on the sequencing results (do most cells contain more 16S reads than target reads?)

20m



Note

IMPORTANT - BARCODE CONCENTRATIONS MAY NEED TO BE MEASURED


Typically, 0.5uL of 1pM barcode will give approximately 1 barcode for every 10 droplets. However, it is best to measure the barcode encapsulation rate by making PCR droplets containing the barcodes at the expected dilution and P7 and Barrev primers targeting the barcode for amplification. Visualize the resulting PCR emulsion using SYBRgreen staining under the microscope to obtain the real encapsulation ratio. Typically, the real barcode concentration can be ~5 fold off from the expected concentration based on manufacturer's labelling.

Note

Barcode oligos should always be freshly diluted from 500pM to 1pM before use, as we have found gradual loss of barcodes over time in a 1pM solution.

- 31 Load the PCR mastermix into the syringe following this protocol

Load  500 μ L of

 Biorad Evagreen Droplet Oil **Bio-Rad Laboratories Catalog ##1864005** into a syringe following this protocol

Protocol

NAME

Loading Syringes to Inject into Microfluidics Device

CREATED BY

freeman.lan

[Preview](#)

- 32 Connect the syringes to the **DoTA-seq Step 2** microfluidics device.

7m

Run the syringe pumps at **200uL/hr** for the gel and PCR mastermix, and **800uL/hr** for the oil syringe.



Collect droplets in an



Axygen® 0.2 mL Maxymum Recovery® Thin Wall PCR
Tubes **Corning Catalog #PCR-02-L-C**

for 00:07:00 for every 25 µL of PCR mastermix or until the PCR mastermix runs out.

33 Use a pipette to remove the oil in the PCR tube, leaving just the emulsion layer (it's okay to have a little bit of oil remaining).

34 Thermocycle the PCR emulsion as follows:

4h

95 °C 5 min

40 cycles of:

95 °C 30s

72 °C 10s

60 °C 5 min

72 °C 30s

Final incubation of:

72 °C 10min

12 °C Hold

All ramp times are at 1 °C per second

PCR Cleanup

1m


35 Keep the emulsion on ice to prevent polymerase activity

Add 25 µL EDTA (0.5 M), pH 8.0 **Life Technologies Catalog #AM9260G** to
the emulsion
Vortex the emulsion to mix

36 Add 25 µL HFE 7500 Perfluorinated Oil to the emulsion



Add  25 μL

 Perfluorooctanol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #370533** to the emulsion

Vortex the emulsion to mix


37 Wait  00:01:00 , then pulse centrifuge to separate the PCR mix from the oil.

1m

Note: If you do not see a clear separation of two clear phases and no emulsion remaining, repeat step 35.


Transfer the top aqueous phase to a new 1mL tube.

38 Add  20 μL [M] 1 Molarity (M)

 TCEP-HCl **Gold Biotechnology Catalog #TCEP** to the tube and vortex to completely decrosslink any remaining gels.


Note

You should be unable to obtain any "jellyish" substance by centrifugation! If there is any jellyish substance left it is not fully de-crosslinked. **ADD MORE TCEP.**

39 Clean up the PCR reaction using the  DNA Clean & Concentrator™-5 **Zymo Research Catalog #D4003** kit.

10m

Elute in  50 μL Elution Buffer.

40 Remove primer dimers and free barcodes using the  SPRIselect reagent kit **Beckman Coulter Catalog #B23317** with 0.7X volume of beads.

10m

41 Check the resulting library for primer dimers using

10m



Equipment

TapeStation

NAME

Agilent

BRAND

G2991AA

SKU

<https://www.agilent.com/en/product/tapestation-automated-electrophoresis/tapestation-instruments/4200-tapestation-system-228263>


LINK

with a

 High Sensitivity D1000 ScreenTape **Agilent Technologies Catalog #5067-5584**


Other high sensitivity capillary electrophoresis methods will also work.

There should be minimal primer dimers on the trace. Below is an example of an acceptable trace.



Example of an acceptable TapeStation trace.

42 Quantify the library using a qPCR library quantification kit such as


 NEBNext Library Quant Kit for Illumina - 100 rxns **New England Biolabs Catalog #E7630S**

1h



Note

Note that you must use a PCR based library quantification kit as not all amplicons contain all the adaptors for sequencing and therefore will throw off sequence non-specific forms of quantification!

- 43 Sequence the library on an Illumina sequencer using Custom Sequencing Primers listed here.  DoTA-seq-Oligo-Sequences.xlsx

Citations

Step 29

Demaree B, Weisgerber D, Lan F, Abate AR. An Ultrahigh-throughput Microfluidic Platform for Single-cell Genome Sequencing.

<https://doi.org/10.3791/57598>