

Oct 18, 2019

Double Digestion and Dephosphorylation of Plasmid

DOI

dx.doi.org/10.17504/protocols.io.79phr5n

Igem Dusseldorf¹

¹Heinrich-Heine Universität Düsseldorf



Igem Dusseldorf

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.79phr5n

Protocol Citation: Igem Dusseldorf 2019. Double Digestion and Dephosphorylation of Plasmid. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.79phr5n>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 16, 2019

Last Modified: October 18, 2019

Protocol Integer ID: 28687

1

Mix the following components gently:

Component	ng/ μ L
Plasmid DNA Tube X	
Component	Volume (μ L)
Sterile MilliQ Water	Fill up to 30.0 μ L
10x FastDigest Buffer	3.0
Plasmid DNA	Total ng of DNA/Plasmid DNA concentration ng/ μ L
FastDigest Enzyme 1	1.0
FastDigest Enzyme 2	1.0
FastAP Alkaline Phosphatase	1.0
Total Volume	30.0



Total ng of DNA	1000. 0
-----------------------	------------

- 2 Incubate at enzyme suitable conditions  00:30:00 -  01:00:00 and temperature

Note

Find out suitable conditions on website of particular science company

- 3 Inactivate the enzymes at suitable conditions and temperature

Note

Find out suitable conditions on website of particular science company

- 4 Run gel at 80-150 V until the dye line is approximately 75-80% of the way down the gel. A typical run time is about 1-1.5 hours, depending on the gel concentration and voltage.
- 5 Cut out digested band (~150 mg) and transfer to pre-weighed tube. Check weight.
- 6 Purify the band