

Oct 16, 2019

Double Digestion of Insert DNA

DOI

dx.doi.org/10.17504/protocols.io.79ehr3e

Igem Dusseldorf¹

¹Heinrich-Heine Universität Düsseldorf



Igem Dusseldorf

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.79ehr3e

Protocol Citation: Igem Dusseldorf 2019. Double Digestion of Insert DNA. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.79ehr3e>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 16, 2019

Last Modified: October 16, 2019

Protocol Integer ID: 28678

Materials

Sterile MiliQ water

10x Fast Digest Buffer

PCR-amplified DNA

Fast Digest Enzyme 1

Fast Digest Enzyme 2



Double digest of insert DNA

- 1 *Combined enzyme volume should not exceed 1/10 of the total reaction volume.*
For **double** digestion with **two** FastDigest restriction enzymes, mix:

1m

Component	Amount
10x FastDigest Buffer	2 µl
DNA	400 ng
FastDigest Enzyme 1	1 µl
FastDigest Enzyme 2	1 µl
Sterile MilliQ Water	Fill up with Sterile MilliQ Water to 20 µl

Note

Total Volume can vary. Other components need to be adjusted.

- 2 Mix gently and spin down



3 Incubate at enzyme suitable conditions

4m

Note

Find out suitable conditions on website of particular science company

4 Inactivate the enzyme at suitable conditions

1m

Note

Find out suitable conditions on website of particular science company

5

Note

Optional: Load on 1% agarose gel and run (duration and voltage vary)