

Jul 04, 2025

Double immunofluorescence – Paraffin-embedded sections

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlky5b6g5r/v1>

Raheleh Roudi¹

¹University of Stanford

Raheleh



Raheleh Roudi

University of Stanford

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlky5b6g5r/v1>

Protocol Citation: Raheleh Roudi 2025. Double immunofluorescence – Paraffin-embedded sections. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.j8nlky5b6g5r/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: July 04, 2025

Last Modified: July 04, 2025

Protocol Integer ID: 221773

Keywords: double immunofluorescence procedure, double immunofluorescence, different antigens in the same sample, different antigen, primary antibody, different species, paraffin, same sample, antibody

Funders Acknowledgements:

Raheleh Roudi

Grant ID: U24CA264298

Heike E. Daldrup-Link

Grant ID: U24CA264298

Abstract

In order to be able to examine the co-distribution of two (or more) different antigens in the same sample, a double immunofluorescence procedure can be carried out. Primary antibodies raised in different species can be used either in parallel (in a mixture) or in a sequential way.

Preparation of slides and samples

- 1
 - 1.1. Deparaffinize sections in xylene 2×5 min.
 - 1.2. Hydrate with 100% ethanol 2×3 min.
 - 1.3. Hydrate with 95% ethanol 1 min.
 - 1.4. Rinse in distilled water

2. Fixation

- 2
 - 2.1 Fix the samples either in ice-cold methanol, acetone (1-10 min) or in 3-4% paraformaldehyde in PBS pH 7.4 for 15 min at room temperature.
 - 2.2 Wash the samples twice with ice cold PBS.

3. Permeabilization

- 3
 - 3.1 Incubate the samples for 10 min with PBS containing 0.25% Triton X-100
 - 3.2 Wash cells in PBS three times for 5 mins.

4. Blocking and simultaneous incubation

- 4
 - 4.1. Incubate tissues with 1% BSA in PBST for 30 min to block unspecific binding of the antibodies
 - 4.2. Incubate tissues in **the mixture of two primary antibodies** in 1% BSA in PBST in a humidified chamber for overnight at 4° C.
 - 4.3. Decant the mixture solution and wash the cells three times in PBS, 5 mins each wash.
 - 4.4. Incubate cells with **the mixture of two secondary antibodies** which are raised in different species in 1% BSA for 1 hr at room temperature in dark.
 - 4.5. Decant the mixture of the secondary antibody solution and wash three times with PBS for 5 min each in dark.

5. Counter staining

- 5
 - 5.1. Rinse with PBS.
 - 5.2 Incubate tissues on 0.1-1 µg/ml DAPI for 1 min.

6. Mounting

- 6
 - 6.1. Seal coverslip with nail polish to prevent drying and movement under microscope.



6.2. Store in dark at -20°C or 4°C.