



Dec 15, 2023

DQ-BSA and BMV109 microscopy

 Molecular Neurodegeneration

DOI

<https://dx.doi.org/10.17504/protocols.io.261gedk5ov47/v1>

Rebecca Wallings¹

¹University of Florida



Rebecca Wallings

University of Florida

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261gedk5ov47/v1>

External link: <https://doi.org/10.1186/s13024-025-00829-w>

Protocol Citation: Rebecca Wallings 2023. DQ-BSA and BMV109 microscopy. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.261gedk5ov47/v1>

Manuscript citation:

Wallings RL, Gillett DA, Staley HA, Mahn S, Mark J, Neighbarger N, Kordasiewicz H, Hirst WD, Tansey MG (2025) ASO-mediated knock-down of GPNMB in mutant-GRN and in Grn-deficient peripheral myeloid cells disrupts lysosomal function and immune responses. *Molecular Neurodegeneration* 20(). doi: [10.1186/s13024-025-00829-w](https://doi.org/10.1186/s13024-025-00829-w)



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 14, 2023

Last Modified: December 15, 2023

Protocol Integer ID: 92330

Keywords: peritoneal macrophage, bmv109 microscopy dq, bmv109 microscopy, microscopy dq, microscopy, bmv109

Funders Acknowledgements:

ASAP

Grant ID: ASAP-020527

Abstract

DQ-BSA and BMV109 microscopy for peritoneal macrophages



- 1 BMV109 Pan Cathepsin probe (Vergent Bioscience) and DQ red BSA (Invitrogen) were added to each well at a final concentration of 1 μ M and 10 μ G/mL, respectively, and cells were incubated at 37°C for 1-hour.

Cells were washed 3 x DPBS +/- and fixed for 10 minutes at room temperature in Invitrogen eBioscience Intracellular Fixation buffer.

Cells were washed 3 x DPBS +/- and incubated in 1 μ g/ml DAPI (Life Technologies) for 10 minutes at room temperature in DPBS +/-.

Cells were imaged using an EVOS M7000 (Invitrogen) at 20 x magnification. Image analysis was performed using Cellprofiler 4.2.5.