



Feb 11, 2020

# Dual Optical Mapping of Action Potentials and Calcium Transients in the Mouse Heart during Optogenetic Stimulation of the ICNS

DOI

[dx.doi.org/10.17504/protocols.io.bcdtis6n](https://dx.doi.org/10.17504/protocols.io.bcdtis6n)

Pradeep Rajendran<sup>1</sup>, Guy Salama<sup>2</sup>, Ching Zhu<sup>1</sup>, Peter Hanna<sup>1</sup>

<sup>1</sup>UCLA; <sup>2</sup>University of Pittsburgh

SPARC

Tech. support email: [info@neuinfo.org](mailto:info@neuinfo.org)



Ching Zhu

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bcdtis6n>

**Protocol Citation:** Pradeep Rajendran, Guy Salama, Ching Zhu, Peter Hanna 2020. Dual Optical Mapping of Action Potentials and Calcium Transients in the Mouse Heart during Optogenetic Stimulation of the ICNS. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bcdtis6n>

**Manuscript citation:**

Salama G, Hwang SM. Simultaneous optical mapping of intracellular free calcium and action potentials from langendorff perfused hearts. Curr Protoc Cytom. 2009;Chapter 12(Unit 12):17.

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** In development

**We are still developing and optimizing this protocol.**

**Created:** February 11, 2020

**Last Modified:** February 11, 2020

**Protocol Integer ID:** 32915

**Keywords:** optogenetic stimulation, calcium transients in the mouse heart, mouse heart, dual optical mapping, calcium transient, action potential

## Troubleshooting

- 1 The animal is euthanized by a method approved by IACUC (Institutional Animal Care and Use Committee), the chest opened, the heart and lungs removed, and the aorta quickly cannulated on a Langendorff apparatus to perfuse the coronary vessels. The cannulation step should occur in ice-cold Tyrode's solution (112mM NaCl, 5mM KCl, 25mM NaHCO<sub>3</sub>, 1mM KH<sub>2</sub>PO<sub>4</sub>, 1.2mM MgSO<sub>4</sub>, 1.8mM CaCl<sub>2</sub>, 50mM glucose) and completed within 10 minutes of euthanasia.
- 2 With a Langendorff apparatus, hearts can be perfused at either constant coronary pressure or constant flow rate. Hearts were immobilized in a custom-built chamber to reduce motion artifact.
- 3 Hearts were stained with bolus injections of voltage-sensitive RH237 (30uL of 1mg/mL in DMSO, Thermo Fisher Scientific, S1109) and Ca indicator Rhod-2 AM (30uL of 1mg/ml in DMSO, Biotium, 50024) into the coronary perfusate.
- 4 Excitation-contraction uncoupling agents such as blebbistatin were not used.
- 5 Light from a 100-W tungsten lamp was collimated, passed through 530 ± 30 nm interference filters, split by a 560 nm dichroic mirror, and focused on the dorsal epicardial surface of the heart for excitation. Emitted fluorescence was collected with tandem camera lenses (50 mm f/1.2 mm, Nikon and 50 mm f/0.95, Navitar) and split with a 600 nm dichroic mirror.
- 6 The longer wavelength moiety, containing the V<sub>m</sub> signal, was filtered between 610-750 nm and the shorter wavelength moiety, containing the intracellular Ca<sup>2+</sup> signal, was filtered between 570-595 nm. The emitted fluorescence signals were recorded using 2 CMOS cameras (SciMedia, MiCAM ULTIMA) with a sampling rate of 2 kHz and 100 × 100 pixels with a 5 × 5 mm field of view. Pixel resolution of the images was 150 × 150 μm<sup>2</sup>. All stimulations were performed on the dorsal heart at 10 Hz, 10 ms, and 221 mW for 10 s. Data were acquired in 40 s intervals with 15 s collected before and after the stimulation.