

Dec 19, 2019 Version 3

# Dye-terminator DNA sequencing V.3

DOI

[dx.doi.org/10.17504/protocols.io.baq7idzn](https://dx.doi.org/10.17504/protocols.io.baq7idzn)

Diep R Ganguly<sup>1</sup>

<sup>1</sup>University of Pennsylvania

Pogson Group



Diep R Ganguly

University of Pennsylvania, The Australian National Universi...

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.baq7idzn](https://dx.doi.org/10.17504/protocols.io.baq7idzn)

**Protocol Citation:** Diep R Ganguly 2019. Dye-terminator DNA sequencing. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.baq7idzn>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** December 19, 2019

**Last Modified:** December 19, 2019

**Protocol Integer ID:** 31231

**Keywords:** Sanger sequencing, Dye-terminator sequencing

## Abstract

This protocol (based on the BigDye® Terminator v3.1 Cycle Sequencing Kit) is for performing terminator cycling sequencing reactions for Sanger sequencing of amplified PCR products or plasmid DNA on the 3130X genetic analyser (Applied Biosystems).

## Attachments



BigDye Terminator v3...

350KB

## Materials

### MATERIALS

- ⊗ Antarctic Phosphatase - 1,000 units **New England Biolabs Catalog #M0289S**
- ⊗ 96 well PCR Plate Non-skirted **Phenix Research Catalog #MPS-499**
- ⊗ Nuclease-free water (e.g. MilliQ or HPLC grade water)
- ⊗ primers
- ⊗ EDTA
- ⊗ 10 mM dNTPs **Life Technologies Catalog #10297-018**
- ⊗ Ethanol **Merck Millipore (EMD Millipore) Catalog #100983**
- ⊗ BigDye™ Terminator v3.1 Cycle Sequencing Kit **Thermo Fisher Catalog #4337454**
- ⊗ Exonuclease I (E. coli) **NEB Catalog #M0293S**
- ⊗ Hi-Di™ Formamide **Thermo Fisher Scientific Catalog #4311320**

## Before start

Optimize PCR cycling (if sequencing amplified PCR products) to ensure your reaction produces a single product. Perform gel excision or PCR clean-up with the potential inclusion of incubating with Antarctic phosphatase and Exonuclease 1 to dephosphorylate and degrade unincorporated dNTPs in PCR reactions to prepare templates for DNA sequencing.

## Terminator cycling reaction

- 1 Perform sequencing reaction in PCR tubes (or 96-well plate) with BigDye Terminator cycling kit and forward or reverse primers.

30m

| Component                                 | Volume (μl)                            |
|---|--|
| v3.1 Ready reaction mix                   | 1                                      |
| 5X Sequencing buffer                      | 1.5                                    |
| 20 μM F/R Primer                          | 0.5                                    |
| Template (plasmid or cleaned PCR product) | 50-150 ng DNA (plasmid or PCR product) |
| Nuclease-free water                       | to 10 μl                               |

BigDye Terminator Cycling reaction

5x reaction buffer: 400 mM TRIS, 10 mM MgCl<sub>2</sub>

- 2 Run the following thermal cycling protocol:
  1. 1 min at 96 °C
  2. 30-40 cycles: 96 °C for 10 seconds, 50 °C for 5 seconds, and 60 °C for 4 min.
  3. Hold at 4-12 °C.

4h

## Purification

1h 30m



- 3 Transfer PCR reaction to eppendorf tube. To the reaction, add 2.5  $\mu$ L of 125 mM EDTA (make sure it touches bottom of tube).
- 4 Add 30  $\mu$ L of 100% ethanol, mix well (inversion).
- 5 Incubate at room temperature for 15 minutes.
- 6 Centrifuge at 4 °C at max speed for 30 minutes.
- 7 Discard supernatant and add 50  $\mu$ L of ice-cold 70% ethanol.
- 8 Centrifuge at 4 °C at max speed for 5 minutes.
- 9 Discard supernatant and allow to air-dry in the dark for >15 minutes.

## Prepare for sequencing

- 10 Resuspend the pellet (likely transparent) in 7.5  $\mu$ L HiDi Formamide (add to any empty wells). Incubate at RT for 5 minutes then transfer to plate. Spin down briefly.
- 11 Incubate plate at 95 °C for 3 minutes (denature) then place immediately on ice. Spin down briefly.
- 12 Submit for sequencing on 3130X genetic analyser (Applied Biosystems). Keep samples on ice.