

Mar 05, 2019

E. coli K12 DNA Extraction

 Forked from [E. coli K12 DNA Extraction](#)

DOI

dx.doi.org/10.17504/protocols.io.yujfwun

Kenneth Schackart¹, Kattika Kaarj¹


¹University of Arizona

481b Laboratory



Kenneth Schackart

University of Arizona



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yujfwun>

Protocol Citation: Kenneth Schackart, Kattika Kaarj 2019. *E. coli* K12 DNA Extraction. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.yujfwun>



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 05, 2019

Last Modified: March 05, 2019

Protocol Integer ID: 21099

Keywords: k12 dna extraction, dna extraction, extraction, dna

Abstract

How to extract DNA from *E. coli* K12 using Wizard® Genomic DNA Purification Kit by Promega®.

I do not claim any credit for the development of this protocol. It has been adapted from the protocol detailed in:



Wizard Genomic DNA Purification....

Guidelines

Lab coat and gloves must be worn at all times.

Materials

MATERIALS

Wizard(R) Genomic DNA Purification Kit **Promega Catalog #A1620**

Reagents in kit:

- Nuclei Lysis Solution
- RNase Solution
- Protein Precipitation Solution
- DNA Rehydration Solution

Additional Reagents:

- Isopropanol
- 70% Ethanol

Additional Materials:

- 2 mL centrifuge tubes
- Pipettes and tips
- Heating block
- Centrifuge



Troubleshooting

Before start



Spray work area with 70% EtOH solution.



Culture bacteria

- 1 Culture *E. coli* K12 in BHI broth overnight.
 2 mg lyophilized *E. coli* K12 in  10 mL BHI broth.

Pellet the cells



- 2 Add  1 mL cell suspension to 2 mL microcentrifuge tube.
- 3 Label centrifuge tube with your group number.
- 4 Centrifuge at $13,000\text{--}16,000 \times g$ for  00:02:00 .
- 5 Remove supernatant.

Lyse nuclei

- 6 Add  600 μL of Nuclei Lysis Solution.


Note

Nuclei Lysis Solution is marked "NL"

- 7 Gently pipet until the cells are resuspended.
- 8 Incubate at  80 °C on heating block for  00:05:00 to lyse the cells.
- 9 Cool to room temperature.






Degrade RNA

10 Add  3 μL RNase Solution to the cell lysate.

Note

RNase solution is marked "R"

11 Invert 2-5 times to mix.

12 Incubate at  37 °C for  00:15:00 to  01:00:00 .


13 Cool to room temperature.


Precipitate proteins


14 Add  200 μL of Protein Precipitation Solution to the RNase-treated cell lysate.

Note

Protein Precipitation solution is marked "P"


15 Vortex vigorously at high speed for  00:00:20 .

16 Incubate on ice for  00:05:00

17 Centrifuge at 13,000-16,000 $\times g$ for  00:03:00 .



Harvest DNA

- 18 Transfer the supernatant containing the DNA to a clean 1.5 mL microcentrifuge tube containing  600 μL isopropanol.

Note



Some supernatant may remain in the original tube containing the protein pellet. Leave this residual to avoid contaminating the DNA solution with the precipitated protein.

Note

Isopropanol is marked "IPA"

- 19 Label centrifuge tube with your group number.
- 20 Gently mix by inversion until the thread-like strands of DNA form a visible mass.

Wash and dry DNA

- 21 Centrifuge at 13,000–16,000 $\times g$ for  00:02:00 .
- 22 Carefully pour off the supernatant and drain the tube on clean absorbent paper.
- 23 Add  600 μL of 70% ethanol and gently invert the tube several times to wash the DNA pellet.

Note

70% ethanol solution is marked "Et"



- 24 Centrifuge at 13,000-16,000 × *g* for 00:02:00 .
- 25 Carefully pour off the ethanol.
- 26 Drain the tube on clean absorbent paper and allow to air-dry for 10-15 minutes.

Rehydrate DNA

- 27 Add 100 µL of DNA rehydration solution to the tube.

Note

DNA Rehydration solution is marked "DR"

- 28 Rehydrate by incubating the solution overnight at room temperature or 4 °C .
- 29 Store DNA at 2 °C to 8 °C .