

Oct 14, 2019

## E. coli transformation

DOI

[dx.doi.org/10.17504/protocols.io.76jhrcn](https://dx.doi.org/10.17504/protocols.io.76jhrcn)

Igem Dusseldorf<sup>1</sup>

<sup>1</sup>Heinrich-Heine Universität Düsseldorf



Igem Dusseldorf

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.76jhrcn](https://dx.doi.org/10.17504/protocols.io.76jhrcn)

**Protocol Citation:** Igem Dusseldorf 2019. E. coli transformation. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.76jhrcn>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 14, 2019

**Last Modified:** October 14, 2019

**Protocol Integer ID:** 28587



## Abstract

- Take tubes with 50  $\mu$ L of TOP10F' competent cells from -80°C freezer and place on ice.
- Add 5  $\mu$ L of each ligation reaction directly to competent cells and mix by tapping gently. **Do not mix cells by pipetting up and down.** Store any remaining ligation at -20°C.
- Incubate cells on ice for 30 min. [while waiting, turn on 42°C water bath]
- Incubate for exactly 30 sec in 42°C water bath. Do not mix or shake.
- Remove vials and place quickly on ice. Store for 2 min on ice.
- Add 500  $\mu$ L of LB medium (pH. 7.5) to each tube.
- Incubate tubes (taped horizontally to platform) for 60 min at 37°C and 225 rpm.
- Plate 550  $\mu$ L (variable) of each transformation on LB + AB plates.
- Incubate overnight at 37°C, in the dark.

