

Oct 14, 2019

E. coli transformation

DOI

dx.doi.org/10.17504/protocols.io.76jhrcn

Igem Dusseldorf¹

¹Heinrich-Heine Universität Düsseldorf



Igem Dusseldorf

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.76jhrcn>

Protocol Citation: Igem Dusseldorf 2019. E. coli transformation. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.76jhrcn>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 14, 2019



Last Modified: October 14, 2019

Protocol Integer ID: 28587

Keywords: incubate cells on ice, incubate cell, μ L of each ligation reaction, ligation reaction, tube, remaining ligation, incubate

Abstract

- Take tubes with 50 μ L of TOP10F' competent cells from -80°C freezer and place on ice.
- Add 5 μ L of each ligation reaction directly to competent cells and mix by tapping gently. **Do not mix cells by pipetting up and down.** Store any remaining ligation at -20°C.
- Incubate cells on ice for 30 min. [while waiting, turn on 42°C water bath]
- Incubate for exactly 30 sec in 42°C water bath. Do not mix or shake.
- Remove vials and place quickly on ice. Store for 2 min on ice.
- Add 500 μ L of LB medium (pH. 7.5) to each tube.
- Incubate tubes (taped horizontally to platform) for 60 min at 37°C and 225 rpm.
- Plate 550 μ L (variable) of each transformation on LB + AB plates.
- Incubate overnight at 37°C, in the dark.

Troubleshooting

