



Dec 15, 2023

Ea(52-68) uptake assays

 Molecular Neurodegeneration

DOI

<https://dx.doi.org/10.17504/protocols.io.14egn3r3pl5d/v1>

Rebecca Wallings¹

¹University of Florida



Rebecca Wallings

University of Florida

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn3r3pl5d/v1>

External link: <https://doi.org/10.1186/s13024-025-00829-w>

Protocol Citation: Rebecca Wallings 2023. Ea(52-68) uptake assays. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.14egn3r3pl5d/v1>

Manuscript citation:

Wallings RL, Gillett DA, Staley HA, Mahn S, Mark J, Neighbarger N, Kordasiewicz H, Hirst WD, Tansey MG (2025) ASO-mediated knock-down of GPNMB in mutant-GRN and in Grn-deficient peripheral myeloid cells disrupts lysosomal function and immune responses. *Molecular Neurodegeneration* 20(). doi: [10.1186/s13024-025-00829-w](https://doi.org/10.1186/s13024-025-00829-w)

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 14, 2023

Last Modified: December 15, 2023

Protocol Integer ID: 92327

Keywords: uptake assays on pmac, uptake assay, pmac, uptake

Funders Acknowledgements:

ASAP

Grant ID: ASAP-020527

Abstract

Ea(52-68) uptake assays on pMacs



- 1 Reconstitute MHC-II Ea chain (Ea) (52–68) peptide (Anaspec) in sterile distilled H₂O to a final concentrate of 1mg/mL.

Once peritoneal macrophages are adhered to plates, add 5µg per well of 24-well plate and incubate for 18-hours

Harvest cells as per flow cytometry protocol and stain with eBioY-Ae antibody (ThermoFisher) to assess Ea52-68 peptide bound to I-Ab