

Sep 27, 2019

Ebola virus sequencing protocol

DOI

dx.doi.org/10.17504/protocols.io.7nwhmfe

Josh Quick¹

¹University of Birmingham

ARTIC



Josh Quick

University of Birmingham



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.7nwhmfe>

Protocol Citation: Josh Quick 2019. Ebola virus sequencing protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.7nwhmfe>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 25, 2019

Last Modified: September 27, 2019

Protocol Integer ID: 28086





Keywords: ebola virus, protocol

Troubleshooting



cDNA preparation

- 1 Mix the following components in an 0.2mL 8-strip tube;


Component	Volume
50µM random hexamers	 1 µL
10mM dNTPs mix (10mM each)	 1 µL
Template RNA	 10 µL
Total	 12 µL


Note

Viral RNA input from a clinical sample should be between Ct 18-35. If Ct is between 12-15, then dilute the sample 100-fold in water, if between 15-18 then dilute 10-fold in water. This will reduce the likelihood of PCR-inhibition.






- 2 Gently mix by pipetting and pulse spin the tube to collect liquid at the bottom of the tube.

- 3 Incubate the reaction as follows:

 65 °C for  00:05:00

Place on ice for  00:01:00

- 4 Add the following to the annealed template RNA:

Component	Volume
SSIV Buffer	 4 µL
100mM DTT	 1 µL
RNaseOUT RNase Inhibitor	 1 µL
SSIV Reverse Transcriptase	 1 µL
Total	 20 µL

- 5 Gently mix by pipetting and pulse spin the tube to collect liquid at the bottom of the tube.



6 Incubate the reaction as follows:

42 °C 01:30:00

70 °C 00:10:00

Hold at 5 °C

Primer pool preparation

7 If required resuspend lyophilised primers at a concentration of 100µM each

Note

Ebola V2 primers for this protocol were designed using **Primal Scheme** and generate overlapping 400nt amplicons. Primer names and dilutions are listed in the table below.

8 Generate primer pool stocks by adding 5 µL of each primer pair to a 1.5 mL Eppendorf labelled either "Pool 1 (100µM)" or "Pool 2 (100µM)". Total volume should be 505 µL for Pool 1 (100µM) and 530 µL for Pool 2 (100µM). These are your 100µM stocks of each primer pool.

Note

Primers should be prepped and aliquoted prior to departure in a sterile PCR cabinet. At no stage should primers or PCR reagents be anywhere near the template until cDNA addition.

9 Dilute this primer pool 1:10 in molecular grade water, to generate 10µM primer stocks. It is recommend that multiple aliquots of each primer pool are made to in case of degradation or contamination.

	Name	Sequence	Name	Sequence	Po ol	Stoc k
	Ebov-10-Pan_1_LE FT	TGTGTGCGAATAA CTATGAGGAAGA	Ebov-10-Pan_1_RIGHT	TTTCCAATGTTTTA CCCCAAGCTTT	1	100μM
			Ebov-10-Pan_1_RIGHT_a lt1	TTTCCAATGCTTTA CCCCAAGCTTT	1	100μM
			Ebov-10-Pan_1_RIGHT_a lt2	TTTCCAATGTTTTA CCCCAAGTTTT	1	100μM
	Ebov-10-Pan_2_LE FT	CAAGCAAGATTGA GAATTAACCTTGG T	Ebov-10-Pan_2_RIGHT	ATCTCCCTGGTACG CATGATGA	2	100μM
	Ebov-10-Pan_2_LE FT_a lt1	CAAGCAAGATTGA GAATTAACCTTGA T	Ebov-10-Pan_2_RIGHT_ alt1	ATCTCCTTGGTACG CATGATGA	2	100μM
	Ebov-10-Pan_3_LE FT	GGCCTTTGAAGC AGGTGTTGAT	Ebov-10-Pan_3_RIGHT	TCAGTCCTTGCTCT GCATGTAC	1	100μM
	Ebov-10-Pan_4_LE FT	CCTTTGCAAGTCT ATTCCTTCCGA	Ebov-10-Pan_4_RIGHT	CTGAGTGCAGCCT TAAAGGAGT	2	100μM
	Ebov-10-Pan_4_LE FT_a lt1	CTTTTGCAAGTCT ATTCCTTCCGA			2	100μM
	Ebov-10-Pan_5_LE FT	AGTTCGTCTCCAT CCTCTTGCA	Ebov-10-Pan_5_RIGHT	CTGGAAGCTGATTT CGTTCTTTTTCT	1	100μM
	Ebov-10-Pan_6_LE FT	GAGTCTCGCGAA CTTGACCATC	Ebov-10-Pan_6_RIGHT	TCCTCGTCGTCCT CGTCTAGAT	2	100μM

	Ebov-10-Pan_6_LE FT_a It1	GAATCTCGCGAAC TTGACCATC	Ebov-10-Pan_6_RIGHT_a It1	TCCTCATCGTCCTC GTCTAGAT	2	100μM
	Ebov-10-Pan_7_LE FT	AGCTACGGCGAAT ACCAGAGTT	Ebov-10-Pan_7_RIGHT	GTCCTGTCCTGC TCTTCATCA	1	100μM
			Ebov-10-Pan_7_RIGHT_a It1	GTCCTGTCCTGTT CTTCATCA	1	100μM
			Ebov-10-Pan_7_RIGHT_a It2	GTCCTGTCCTGTT CTTCATCG	1	100μM
	Ebov-10-Pan_8_LE FT	TTAACGAAGAGGC AGACCCACT	Ebov-10-Pan_8_RIGHT	TTCCTCTTCAAGG GAGTCTGGA	2	100μM
	Ebov-10-Pan_8_LE FT_a It1	TCAACGAAGAGG CAGACCCACT	Ebov-10-Pan_8_RIGHT_a It1	TTCCTCTTCAAGG GAGTCCGGA	2	100μM
	Ebov-10-Pan_9_LE FT	GTGACAACACCCA GTCAGAACA	Ebov-10-Pan_9_RIGHT	TCTTCCTGTTTTCG TTCCTTGACT	1	100μM
	Ebov-10-Pan_9_LE FT_a It1	GTGACAACACCCA GCCAGAACA	Ebov-10-Pan_9_RIGHT_a It1	TCTTCCTGTTTGCG TTCCTTGACT	1	100μM
			Ebov-10-Pan_9_RIGHT_a It2	TCTTCCTGTTTGCG TTTCTTGACT	1	100μM
	Ebov-10-Pan_10_L EFT	ACAATGGGATGAT TCAACCGACA	Ebov-10-Pan_10_RIGHT	TCGAGTGCTAGAGA ATTCAATTGACG	2	100μM
	Ebov-10-Pan_10_L	ATAATGGGATGAT TTAACCGACA			2	100μM

	EFT_ alt1					
	Ebov -10- Pan_ 11_L EFT	ACCTACTAGCCTG CCCAACATT	Ebov-10- Pan_11_RIGHT	AATTGGGTCCGTTT GGGTTTGA	1	100μ M
	Ebov -10- Pan_ 11_L EFT_ alt1	ACCTACTAGCCTA CCCAACATT	Ebov-10- Pan_11_RIGHT_ alt1	AATTGGATCCGTTT GGGTTTGA	1	100μ M
	Ebov -10- Pan_ 12_L EFT	CCCAAATGCAACA AACGAAGCC	Ebov-10- Pan_12_RIGHT	TCAATCTTACCCCG AATCGCAC	2	100μ M
	Ebov -10- Pan_ 12_L EFT_ alt1	CCCAAATGCAACA AACAAAGCC	Ebov-10- Pan_12_RIGHT_ alt1	TCAATCTTACCCCG AATTGCAC	2	100μ M
	Ebov -10- Pan_ 13_L EFT	TATTGGGCCGAAC ATGGTCAAC	Ebov-10- Pan_13_RIGHT	TGACAGGTGGAGC AGCATCTTG	1	100μ M
	Ebov -10- Pan_ 13_L EFT_ alt1	TATTGGGCTGAAC ATGGTCAAC			1	100μ M
	Ebov -10- Pan_ 14_L EFT	CATTCATGCTGAG TTCCAGGCC	Ebov-10- Pan_14_RIGHT	GCGAGATATGAACA ATTTTATCTTGGTC G	2	100μ M
			Ebov-10- Pan_14_RIGHT_ alt1	GCGAGATAAGGAC AATTTTATCTTGGT CG	2	100μ M
			Ebov-10- Pan_14_RIGHT_ alt2	GCGAGATAAGAAC AATTTTATCTTGGT CG	2	100μ M
	Ebov -10- Pan_ 15_L EFT	TGAGTATCAGCCC TGGATAATATAAGT CA	Ebov-10- Pan_15_RIGHT	TCGATGGAGTGTCC CCATTGAC	1	100μ M

Ebov-10-Pan_15_L EFT_alt1	TGAGTATCAGCCC TAGATAATATAAGT CA	Ebov-10-Pan_15_RIGHT_alt1	TCGATGGAGTGTCT CCATTGAC	1	100μM
Ebov-10-Pan_16_L EFT	GCAACAGCAATAC AGGCTTCCT	Ebov-10-Pan_16_RIGHT	GAAAGCCTGGTTT CCAATTCGC	2	100μM
Ebov-10-Pan_16_L EFT_alt1	GCAACAACAATAC AGGCTTCCT	Ebov-10-Pan_16_RIGHT_alt1	GAAGGCCTGGTTT CCAATTCGC	2	100μM
Ebov-10-Pan_17_L EFT	CCACTTGTCAGAG TCAATCGGC	Ebov-10-Pan_17_RIGHT	GTTTCTGGCACTTC GATTCCCA	1	100μM
		Ebov-10-Pan_17_RIGHT_alt1	GTTTCTGGCACTTC GATACCCA	1	100μM
Ebov-10-Pan_18_L EFT	AAAATCCAAGCAA TAATGACTTCACT CC	Ebov-10-Pan_18_RIGHT	TTGATCAATTAAAA GTGTCTCCTCTAAT GG	2	100μM
		Ebov-10-Pan_18_RIGHT_alt1	TCGATCAATTTAAA GTATCTCCTCTAAT GG	2	100μM
		Ebov-10-Pan_18_RIGHT_alt2	TTGATCAATTAAAA GTATCTCCTCTAAT AG	2	100μM
Ebov-10-Pan_19_L EFT	AGATCCAGTTTTA TAGAATCTTCTCA GGGA	Ebov-10-Pan_19_RIGHT	AGAAGGGCAATGT CTGTACTTGG	1	100μM
Ebov-10-Pan_19_L EFT_alt1	AGATCCAGTTTTA CAGAATCTTCTCA GGGA	Ebov-10-Pan_19_RIGHT_alt1	AGAAGGGCGATGT CTGTGCTTGG	1	100μM
Ebov-10-Pan_20_L EFT	AGCCAGTGTGACT TGGATTGGA	Ebov-10-Pan_20_RIGHT	AGTTTGTGCGACATC ACTAACCTGT	2	100μM



	20_L EFT					
			Ebov-10- Pan_20_RIGHT _alt1	AGTTTGTGCGACATC ACTAACTTGT	2	100μ M
	Ebov -10- Pan_ 21_L EFT	AGAACATTTTCCA TCCCACTTGGA	Ebov-10- Pan_21_RIGHT	AAGCACCTCTTTA TGGAAGGC	1	100μ M
			Ebov-10- Pan_21_RIGHT_ alt1	AAGCACCTCTTT GTGAAGGC	1	100μ M
	Ebov -10- Pan_ 22_L EFT	TGCCGGTATGTGC ACAAAGTAT	Ebov-10- Pan_22_RIGHT	ATATATTGTCTCATT CAGCTGGAGCA	2	100μ M
	Ebov -10- Pan_ 23_L EFT	CGAGGTTGACAA TTGACCTACGT	Ebov-10- Pan_23_RIGHT	GCAAGGGTTGTTA GATGCGACA	1	100μ M
			Ebov-10- Pan_23_RIGHT _alt1	GCAAGGGTTGTCA GATGCGACA	1	100μ M
	Ebov -10- Pan_ 24_L EFT	TGCAATGGTTCAA GTGCACAGT	Ebov-10- Pan_24_RIGHT	CTGGCACTCTCTT CTCCGGTAT	2	100μ M
	Ebov -10- Pan_ 24_L EFT_ alt1	TGCAATGGTTCAA GTGCACAAT			2	100μ M
	Ebov -10- Pan_ 25_L EFT	ACCACAACAAGTC CCCAAAACC	Ebov-10- Pan_25_RIGHT	TAGCTCAGTTGTGG CTCTCAGG	1	100μ M
			Ebov-10- Pan_25_RIGHT _alt1	TAGCTCGGTTGTG GCTCTCAGG	1	100μ M
	Ebov -10- Pan_ 26_L EFT	ATCTGTGGGTTGA GACAGCTGG	Ebov-10- Pan_26_RIGHT	GCTTTTCCATGAAG CAATCTGAAGA	2	100μ M



Ebov-10-Pan_26_L EFT_alt1	ATCTGTGGATTGA GGCAGCTGG	Ebov-10-Pan_26_RIGHT_alt1	GCTTTGCCATGAA GCAATCTGAAGA	2	100μM
Ebov-10-Pan_26_L EFT_alt2	ATCTGTGGGTTGA GGCAGCTGG			2	100μM
Ebov-10-Pan_27_L EFT	TGGAGTTACAGGC GTTATAATTGCA	Ebov-10-Pan_27_RIGHT	AAAGGCTTCTTTC CCTTGTCACCT	1	100μM
Ebov-10-Pan_28_L EFT	TCATCCTTGATTC TACAATCATGACA GT	Ebov-10-Pan_28_RIGHT	AGGTGCTGGAGGA ACTGTTAATG	2	100μM
Ebov-10-Pan_28_L EFT_alt1	TCATCCTTGATTC TACAATCATAACA GT			2	100μM
Ebov-10-Pan_29_L EFT	GAGTACCGTCAAT CAAGGAGCG	Ebov-10-Pan_29_RIGHT	CACAGCACATAGAG TCAACAATGC	1	100μM
Ebov-10-Pan_30_L EFT	GATCAAGACGGCA GAACACTGG	Ebov-10-Pan_30_RIGHT	ATCAGACCATGAGC ATGTCCCC	2	100μM
Ebov-10-Pan_31_L EFT	CTGCTGTCGTTGT TTCAGGGTT	Ebov-10-Pan_31_RIGHT	ATGGGATGGATCGT TGCTACCT	1	100μM
		Ebov-10-Pan_31_RIGHT_alt1	ATGGGATGGATCGT TGCTGCCT	1	100μM
		Ebov-10-Pan_31_RIGHT_alt2	ATGAGATGGATCGT TGCTACCT	1	100μM
Ebov-10-	GCCAAGCATACCT CTTGACACAA	Ebov-10-Pan_32_RIGHT	TGGACTACCCTGA AATAGTACTTTGC	2	100μM

	Pan_32_L EFT					
	Ebov-10-Pan_33_L EFT	TGCGGAGGTCTGA TAAGAATAAACC	Ebov-10-Pan_33_RIGHT	TTCAACCTTGAAA CCTTGCGCT	1	100μM
			Ebov-10-Pan_33_RIGHT_alt1	TTCAACTTTGAAAC CTTGCGCT	1	100μM
	Ebov-10-Pan_34_L EFT	GCTGAAAAGAAG CTTACCTACAACG	Ebov-10-Pan_34_RIGHT	TCCTTGTCATTGAC CATGCAGG	2	100μM
	Ebov-10-Pan_34_L EFT_alt1	GTTGAAAAAAGG CCTACCTACAACG			2	100μM
	Ebov-10-Pan_34_L EFT_alt2	GCTGAAAAGAAG CCCACCTACAACG			2	100μM
	Ebov-10-Pan_35_L EFT	GTGACTCACAAAG GAATGGCCC	Ebov-10-Pan_35_RIGHT	ACAATCCGTTGTAG TTCACGACA	1	100μM
			Ebov-10-Pan_35_RIGHT_alt1	ACAACCCGTTGTA GTTACGACA	1	100μM
	Ebov-10-Pan_36_L EFT	TGCTGTCGTTGAT TCGATCCAA	Ebov-10-Pan_36_RIGHT	AGCAGAGATGTCAA GATAACTATTGAGT	2	100μM
	Ebov-10-Pan_37_L EFT	ACACGAATGCAAA GTTTGATTCTTGA	Ebov-10-Pan_37_RIGHT	TGAAACCTAACACA TGTGACCTGC	1	100μM
			Ebov-10-Pan_37_RIGHT_alt1	TGAAACCTAACACA CGTGACCTGC	1	100μM



	Ebov-10-Pan_38_L EFT	CCCTCAAACAAGA GATTCCAAGACA	Ebov-10-Pan_38_RIGHT	ACAGTTGCGTAGTT GCGGATTA	2	100μM
	Ebov-10-Pan_38_L EFT_alt1	CCCTCAAATAAGA GATTCCAAGACA			2	100μM
	Ebov-10-Pan_38_L EFT_alt2	TCCTCAAATAAGA GATTCCAAGACA			2	100μM
	Ebov-10-Pan_39_L EFT	ACCTAGTCACTAG AGCTTGCGG	Ebov-10-Pan_39_RIGHT	ACATTTGATGTAAA AATTCATTGCCCTG	1	100μM
	Ebov-10-Pan_40_L EFT	GTGGGTGCTCAAG AAGACTGTG	Ebov-10-Pan_40_RIGHT	TGAGATTAGAGTTG TGTTGAATCGACA	2	100μM
	Ebov-10-Pan_40_L EFT_alt1	GTGGGTGCTCAAG AGGACTGTG	Ebov-10-Pan_40_RIGHT_alt1	TGAGATTAGAGTCG TGTTGAATCGACA	2	100μM
	Ebov-10-Pan_41_L EFT	AAGAAGCGGTTCA AGGGCATAC	Ebov-10-Pan_41_RIGHT	CTATGGAATTCACG GATCTTTTGAGC	1	100μM
	Ebov-10-Pan_41_L EFT_alt1	AAGAAGCAGTTCA AGGGCATAC	Ebov-10-Pan_41_RIGHT_alt1	CTATGGAATTCACG GATCTTTTGATC	1	100μM
	Ebov-10-Pan_42_L EFT	TGCATTTAGCTGT AAATCACACCCT	Ebov-10-Pan_42_RIGHT	AATCATTGGCAACG GAGGGAAT	2	100μM
			Ebov-10-Pan_42_RIGHT_alt1	AATCATTGGCAACG GGGGGAAT	2	100μM

Ebov-10-Pan_43_L EFT	GTCAAGGATCTTG GTACAGTGTTACT	Ebov-10-Pan_43_RIGHT	TGAGAAAGAAAAG TTCCGATATTGTGG T	1	100μ M
Ebov-10-Pan_43_L EFT_alt1	GCCAAGGGTCTT GGTACAGTGTTAC T	Ebov-10-Pan_43_RIGHT_alt1	TGAGAAAGAAAAA TTCCGGTATTGTGG T	1	100μ M
Ebov-10-Pan_43_L EFT_alt2	GTCAAGGGTCTTG GTACAGTGTTACT	Ebov-10-Pan_43_RIGHT_alt2	TGAGAAAGAAAAA TTCCGATATTGTGG T	1	100μ M
Ebov-10-Pan_44_L EFT	TTGAGAATGTTCT TTCCTACGCACA	Ebov-10-Pan_44_RIGHT	ACGGTTGCAATATT CTATAAAAGGTGC	2	100μ M
Ebov-10-Pan_44_L EFT_alt1	TTGAGAATGTTCT TTCCTACGCGCA	Ebov-10-Pan_44_RIGHT_alt1	ACGGTTGCAATATT CGATAAAAGGTGC	2	100μ M
		Ebov-10-Pan_44_RIGHT_alt2	ACGGTTACAATATT CTATAAAAGGTGC	2	100μ M
Ebov-10-Pan_45_L EFT	CCACAGTTAGAGG GAGTAGCTTTG	Ebov-10-Pan_45_RIGHT	GCTCGTCTGCGTC AGTCTCTAA	1	100μ M
Ebov-10-Pan_45_L EFT_alt1	CCACAGTTAGAGG GAGTAGTTTTG			1	100μ M
Ebov-10-Pan_46_L EFT	AAGTTACGCTCAG CTGTGATGG	Ebov-10-Pan_46_RIGHT	ATGGAAAGCTGCG GTTATCCTG	2	100μ M
Ebov-10-Pan_47_L EFT	TAGGCACTGCTTT TGAGCGATC	Ebov-10-Pan_47_RIGHT	CACAAAGTCAATG GCAGTGCAG	1	100μ M



Ebov-10-Pan_47_L EFT_alt1	TAGGCACCGCTTT TGAGCGGTC			1	100μM
Ebov-10-Pan_47_L EFT_alt2	TAGGCACTGCTTT TGAACGATC			1	100μM
Ebov-10-Pan_48_L EFT	TCTCCGAATGATT GAGATGGATGATT	Ebov-10-Pan_48_RIGHT	CTCAGTCTGTCCA AAACCGGTG	2	100μM
Ebov-10-Pan_48_L EFT_alt1	TCTCCGAATGATT GGGATGGATGATT			2	100μM
Ebov-10-Pan_49_L EFT	GATATCTTTTCAC GCACGCCGA	Ebov-10-Pan_49_RIGHT	CCACCTGGTTGCT TTGCATTG	1	100μM
Ebov-10-Pan_49_L EFT_alt1	GATATCTTTTCAC GCACGCCCA	Ebov-10-Pan_49_RIGHT_alt1	CCACCAGGTTGCT TTGCATTG	1	100μM
Ebov-10-Pan_50_L EFT	TCAAAGTGTTTTG GCTGAAACCCT	Ebov-10-Pan_50_RIGHT	TCCTGAGTAATGTG AAGGGGTCA	2	100μM
Ebov-10-Pan_50_L EFT_alt1	TCAAAGTGTTTTG GCTGAAACCCT	Ebov-10-Pan_50_RIGHT_alt1	TCCTGAGTAATGTG AAGGAGTCA	2	100μM
Ebov-10-Pan_51_L EFT	AACAGTGACTTGC TAATAAAACCATT TTTG	Ebov-10-Pan_51_RIGHT	AAATACTGAGCTGG TACTTCCCG	1	100μM



Ebov-10-Pan_51_L EFT_alt1	AACAGTGACTTGC TAATAAAGCCATT TTTG			1	100μM
Ebov-10-Pan_51_L EFT_alt2	AACAGTGATTTGC TAATAAAACCATT TTTG			1	100μM
Ebov-10-Pan_52_L EFT	AATCGTGCTCACC TTCATCTAACT	Ebov-10-Pan_52_RIGHT	CCCAAACTGTAC AGAAGTCCTATCT	2	100μM
Ebov-10-Pan_53_L EFT	ACAGACCCAATTA GCAGTGGAGA	Ebov-10-Pan_53_RIGHT	ACAATTGTTCCGC GATTAATTATCCAT	1	100μM
Ebov-10-Pan_53_L EFT_alt1	ACAGACCCAATTA GCAGCGGAGA	Ebov-10-Pan_53_RIGHT_alt1	ACAATTGTTCCGC GATTAATTATCCAC T	1	100μM
Ebov-10-Pan_54_L EFT	TCTCAGATGCGGC CAGTTATT	Ebov-10-Pan_54_RIGHT	TGACCATCACTGTT GTTTGTGCT	2	100μM
Ebov-10-Pan_54_L EFT_alt1	TCTCAGATGCGGC CAGATTATT			2	100μM
Ebov-10-Pan_55_L EFT	TGGAGGAGCAGA CACAGAAACA	Ebov-10-Pan_55_RIGHT	ATGACGTTAATTGG CGTGTCCC	1	100μM
Ebov-10-Pan_55_L EFT_alt1	TGGAGGAGCAGG CACAGAAACA	Ebov-10-Pan_55_RIGHT_alt1	ATGACGTCAATTGG CGTGTCCC	1	100μM



Ebov-10-Pan_55_L EFT_alt2	TGGAGAAGCAGG CACAGAAACA	Ebov-10-Pan_55_RIGHT_alt2	ATGACGTTAATTGG CGCGTCCC	1	100μM
Ebov-10-Pan_56_L EFT	CTCACACCGTCTA GTCCTACCT	Ebov-10-Pan_56_RIGHT	TTTGACATAACAGG TAGAAGCATCCT	2	100μM
Ebov-10-Pan_56_L EFT_alt1	CTCGCACCGTCTA GTCCTACCT			2	100μM
Ebov-10-Pan_56_L EFT_alt2	CTCACATCGTCTA GTCCTACCT			2	100μM
Ebov-10-Pan_57_L EFT	ACACGCTAGCTAC TGAGTCCAG	Ebov-10-Pan_57_RIGHT	ATTGGCTTAATTAA ATAACCAGTGGCA	1	100μM
Ebov-10-Pan_58_L EFT	TGAAAGCAGTGGT CCTTAAAGTCT	Ebov-10-Pan_58_RIGHT	TGCTCTAAGATGTG CTAAGTGCTG	2	100μM
		Ebov-10-Pan_58_RIGHT_alt1	TGCTCTAAGATGTG CCAAGTGCTG	2	100μM
Ebov-10-Pan_59_L EFT	CGTCGATTCAAAA AGAGGTCCACT	Ebov-10-Pan_59_RIGHT	TCAGAAGCCCTGT CAGCCTTTC	1	100μM
Ebov-10-Pan_60_L EFT	AGATTGCAATTGT GAAGAACGTTTCT	Ebov-10-Pan_60_RIGHT	AGAGTGCAGAGTTT ATTATGTTGCGT	2	100μM
Ebov-10-Pan_61_L EFT	TCACAATGCAGCA TGTGTGACA	Ebov-10-Pan_61_RIGHT	AGGTATTTCTGATT TTACAGTCCTGCC	1	100μM











		Ebov-10-Pan_61_RIGHT_alt1	AGGTATTTATGATT TTACAGTCCTGCC	1	100μM
		Ebov-10-Pan_61_RIGHT_alt2	AGGTATTTCTGATT TTACAGTCATGCC	1	100μM
Ebov-10-Pan_62_L EFT	CCTGTCAGATGGA ATAGTGTTTTGGT	Ebov-10-Pan_62_RIGHT	AATTTTTGTGTGCG ACCATTTTTTC	2	100μM

Note

Primers need to be used at a final concentration of 0.015μM per primer. In this case, Pool 1 has 101 primers in it so the requirement is 3.8μL of Pool 1 (10μM) primers per 25μL reaction. Pool 2 has 106 primers the requirement is 4.0μL of Pool 2 (10μM) primers per 25μL reaction. For other schemes, adjust the volume added appropriately.

Multiplex PCR

- 10 In the mastermix hood set up the multiplex PCR reactions as follows in 0.2mL 8-strip PCR tubes:

Component	Pool 1	Pool 2
NEB Q5 Polymerase 2X MasterMix	 12.5 μL	 12.5 μL
Primer Pool 1 or 2 (10μM)	 3.8 μL	 4.0 μL
Water	 6.2 μL	 6.0 μL
Total	 22.5 μL	 22.5 μL

Note

The **mastermix hood** should be prepared by sterilising with UV and treated with MediPal wipes, DNaseZap and RNaseZap reagents. Wipe down the hood with each sequentially, allowing 5 minutes for drying between each. Pipettes should also be treated in the same way, and UV treated for 30 mins between library preparations.

This step should be carried out in the mastermix hood and template should not be taken anywhere near the mastermix hood at any stage.










- 11 In the template hood add  2.5 µL cDNA to each tube and mix well by pipetting.

Note

The **template hood** should be prepared by sterilising with UV and treated with MediPal wipes, DNaseZap and RNaseZap reagents. Wipe down the hood with each sequentially, allowing 5 minutes for drying between each. Pipettes should also be treated in the same way, and UV treated for 30 mins between library preparations.

- 12 Pulse centrifuge the tubes to collect the contents at the bottom of the tube.


- 13 Set-up the following program on the thermal cycler:

Step	Temperature	Time	Cycles
Heat Activation	 98 °C	 00:00:30	1
Denaturation	 98 °C	 00:00:15	25-35
Annealing	 65 °C	 00:05:00	25-35
Hold	 4 °C	Indefinite	1

Note

Cycle number should be 25 for Ct 18-21 up to a maximum of 35 cycles for Ct 35

PCR clean-up

- 14 Combine the entire contents of "Pool 1" and "Pool 2" PCR reactions for each biological sample into to a single  1.5 mL Eppendorf tube.

- 15 Clean-up the amplicons using the following protocol:



Protocol

NAME

Amplicon clean-up using SPRI beads







CREATED BY

Josh Quick








[Preview](#)

- 15.1 Vortex SPRI beads thoroughly to ensure they are well resuspended, the solution should be a homogenous brown colour.

 Agencourt AMPure XP **Beckman Coulter Catalog #A63880**

- 15.2 Add an equal volume (1:1) of SPRI beads to the sample tube and mix gently by either flicking or pipetting. For example add  50 μ L SPRI beads to a  50 μ L reaction.
- 15.3 Pulse centrifuge to collect all liquid at the bottom of the tube.
- 15.4 Incubate for  00:05:00 at room temperature.
- 15.5 Place on magnetic rack and incubate for  00:02:00 or until the beads have pelleted and the supernatant is completely clear.
- 15.6 Carefully remove and discard the supernatant, being careful not to touch the bead pellet.
- 15.7 Add  200 μ L of room-temperature  70 % volume ethanol to the pellet.
- 15.8 Carefully remove and discard ethanol, being careful not to touch the bead pellet.



- 15.9  [go to step #15.7](#) and repeat ethanol wash.
- 15.10 Pulse centrifuge to collect all liquid at the bottom of the tube and carefully remove as much residual ethanol as possible using a P10 pipette.
- 15.11 With the tube lid open incubate for  00:01:00 or until the pellet loses it's shine (if the pellet dries completely it will crack and become difficult to resuspend).
- 15.12 Resuspend pellet in  30 μL Elution Buffer (EB), mix gently by either flicking or pipetting and incubate for  00:02:00 .
-  Elution Buffer (EB) **Qiagen Catalog #19086**
- 15.13 Place on magnet and transfer sample to a clean 1.5mL Eppendorf tube ensuring no beads are transferred into this tube.
- 15.14 Quantify  1 μL product using the Quantus Fluorometer using the ONE dsDNA assay.
-  QuantiFluor(R) ONE dsDNA System, 100rxn **Promega Catalog #E4871**

Equipment

Quantus

NAME

Fluorometer

TYPE

Promega

BRAND

E6150

SKU

<https://www.promega.co.uk/products/microplate-readers-fluorometers-luminometers/fluorometers/quantus-fluorometer>

LINK



Quantification and normalisation

- 16 Quantify the amplicon pools using the Quantus Fluorometer using the ONE dsDNA assay.

Protocol

NAME

DNA quantification using the Quantus fluorometer


CREATED BY

Josh Quick


Preview


- 16.1 Remove Lambda DNA 400 ng/ μ L standard from the freezer and leave on ice to thaw. Remove ONE dsDNA dye solution from the fridge and allow to come to room temperature.

 QuantiFluor(R) ONE dsDNA System, 500rxn **Promega Catalog #E4870**







- 16.2 Set up two  0.5 mL tubes for the calibration and label them 'Blank' and 'Standard'

- 16.3 Add  200 μ L ONE dsDNA Dye solution to each tube.

- 16.4 Mix the Lambda DNA standard 400 ng/ μ L standard by pipetting then add  1 μ L to one of the standard tube.



- 16.5 Mix each sample vigorously by vortexing for  00:00:05 and pulse centrifuge to collect the liquid.



- 16.6 Allow both tubes to incubate at room temperature for  00:02:00 before proceeding.
- 16.7 Selection 'Calibrate' then 'ONE DNA' then place the blank sample in the reader then select 'Read Blank'. Now place the standard in the reader and select 'Read Std'.
- 16.8 Set up the required number of  0.5 mL tubes for the number of DNA samples to be quantified.
- Note**
- Use only thin-wall, clear, 0.5mL PCR tubes such as Axygen #PCR-05-C
- 16.9 Label the tubes on the lids, avoid marking the sides of the tube as this could interfere with the sample reading.
- 16.10 Add  199 μ L ONE dsDNA dye solution to each tube.
- 16.11 Add  1 μ L of each user sample to the appropriate tube.
- Note**
- Use a P2 pipette for highest accuracy.
- 16.12 Mix each sample vigorously by vortexing for  00:00:05 and pulse centrifuge to collect the liquid.
- 16.13 Allow all tubes to incubate at room temperature for  00:02:00 before proceeding.
- 16.14 On the Home screen of the Quantus Fluorometer, select 'Protocol', then select 'ONE DNA' as the assay type.



Note

If you have already performed a calibration for the selected assay you can continue, there is no need to perform repeat calibrations when using ONE DNA pre diluted dye solution. If you want to use the previous calibration, skip to step 11. Otherwise, continue with step 9.

- 16.15 On the home screen navigate to 'Sample Volume' and set it to  1 μL then 'Units' and set it to ng/ μL .
- 16.16 Load the first sample into the reader and close the lid. The sample concentration is automatically read when you close the lid.
- 16.17 Repeat step 16 until all samples have been read.
- 16.18 The value displayed on the screen is the dsDNA concentration in ng/ μL , carefully record all results in a spreadsheet or laboratory notebook.
- 17 Label a  1.5 mL Eppendorf tube for each sample.

Note

This is a 'one-pot ligation' protocol for native barcoded ligation libraries. We have seen no reduction in performance compared to standard libraries, and is made faster by using the Ultra II® ligation module which is compatible with the Ultra II® end repair/dA-tailing module removing a clean-up step.

- 18 Normalise the input but adding  10 ng amplicon pools to each tube and diluting to  10 μL to have an input concentration of [M] 1 ng/ μL

Note

Quantity of amplicons can vary from 10-50ng depending on the amplicon length, any more than this and the molarity of DNA ends will be too high for efficient barcoding. You need to have 6 samples per native barcoded library to have sufficient material at the end.



Native barcoding

19 Barcode the amplicon pools using native barcodes.

Protocol

NAME

One-pot native barcoding of amplicons

CREATED BY

Josh Quick

Preview

19.1 Set up the following reaction for each sample:

Component

Volume

DNA amplicons

🧪 10 μ L

Ultra II End Prep Reaction Buffer

🧪 1.4 μ L

Ultra II End Prep Enzyme Mix

🧪 0.6 μ L

Total

🧪 12 μ L

19.2 Incubate at room temperature for 🕒 00:10:00

Incubate at 🌡 65 °C for 🕒 00:05:00

Incubate on ice for 🕒 00:01:00

19.3 Add the following directly to the previous reactions:

Component

Volume

NBXX barcode

🧪 2.5 μ L

Ultra II Ligation Master Mix

🧪 14.5 μ L

Ligation Enhancer

🧪 0.5 μ L

Total


🧪 29.5 μ L

**Note**

Use one native barcode from the EXP-NBD104 (1-12) or EXP-NBD114 (13-24) per sample. Use from 6 to 24 barcodes in a library, any fewer and there will be insufficient total material to achieve good yields.

19.4 Incubate at room temperature for  00:20:00

Incubate at  70 °C for  00:10:00

Incubate on ice for  00:01:00

Note

The 70°C incubation is to inactivate the DNA ligase to prevent barcode cross-ligation when reactions are pooled in the next step.

19.5 Pool all barcoded fragments together into a new 1.5 ml Eppendorf tube.

Protocol

NAME

Amplicon clean-up using SPRI beads

CREATED BY

Josh Quick



19.6 Quantify the barcoded amplicon pools using the Quantus Fluorometer using the ONE dsDNA assay.

**Protocol**

NAME






DNA quantification using the Quantus fluorometer

CREATED BY

Josh Quick

[Preview](#)


19.7 Set up the following AMII adapter ligation reaction:

Component	Volume
Barcoded amplicon pools	 30 μL
NEBNext Quick Ligation Reaction Buffer (5X)	 10 μL
AMII adapter mix	 5 μL
Quick T4 DNA Ligase	 5 μL
Total	 50 μL

Note




The input of barcoded amplicon pools will depend on the number of barcoded pools and should be between 50 ng (6 barcodes) and 200 ng (24 barcodes).

19.8 Incubate at room temperature for  00:20:00

19.9 Add  50 μL (1:1) of SPRI beads to the sample tube and mix gently by either flicking or pipetting.

Note

Vortex SPRI beads thoroughly before use to ensure they are well resuspended, the solution should be a homogenous brown colour.

- 19.10 Pulse centrifuge to collect all liquid at the bottom of the tube.
- 19.11 Incubate for  00:05:00 at room temperature.
- 19.12 Place on magnetic rack and incubate for  00:02:00 or until the beads have pelleted and the supernatant is completely clear.
- 19.13 Carefully remove and discard the supernatant, being careful not to touch the bead pellet.
- 19.14 Add  200 μL SFB and resuspend beads completely by pipette mixing.


Note

SFB will remove excess adapter without damaging the adapter-protein complexes. Do not use 70% ethanol as in early clean-ups.

- 19.15 Pulse centrifuge to collect all liquid at the bottom of the tube.
- 19.16 Remove supernatant and discard.
- 19.17 Repeat steps 14-16 to perform a second SFB wash.
- 19.18 Pulse centrifuge and remove any residual SFB.

Note

You do not need to allow to air dry with SFB washes.

- 19.19 Add  15 μL EB and resuspend beads by pipette mixing.



19.20 Incubate at room temperature for  00:02:00 .

19.21 Place on magnetic rack.

19.22 Transfer final library to a new 1.5mL Eppendorf tube.

20 Quantify the final library using the Quantus Fluorometer using the ONE dsDNA assay.

Protocol

NAME

DNA quantification using the Quantus fluorometer

CREATED BY


Josh Quick

Preview







Note

Final library can be now be stored in EB at 4°C for up to a week if needed otherwise proceed directly to MinION sequencing.

20.1 Remove Lambda DNA 400 ng/μL standard from the freezer and leave on ice to thaw. Remove ONE dsDNA dye solution from the fridge and allow to come to room temperature.



 QuantiFluor(R) ONE dsDNA System, 500rxn **Promega Catalog #E4870**



- 20.2 Set up two  0.5 mL tubes for the calibration and label them 'Blank' and 'Standard'
- 20.3 Add  200 μ L ONE dsDNA Dye solution to each tube.
- 20.4 Mix the Lambda DNA standard 400 ng/ μ L standard by pipetting then add  1 μ L to one of the standard tube.
- 20.5 Mix each sample vigorously by vortexing for  00:00:05 and pulse centrifuge to collect the liquid.
- 20.6 Allow both tubes to incubate at room temperature for  00:02:00 before proceeding.
- 20.7 Selection 'Calibrate' then 'ONE DNA' then place the blank sample in the reader then select 'Read Blank'. Now place the standard in the reader and select 'Read Std'.
- 20.8 Set up the required number of  0.5 mL tubes for the number of DNA samples to be quantified.


Note

Use only thin-wall, clear, 0.5mL PCR tubes such as Axygen #PCR-05-C

- 20.9 Label the tubes on the lids, avoid marking the sides of the tube as this could interfere with the sample reading.
- 20.10 Add  199 μ L ONE dsDNA dye solution to each tube.
- 20.11 Add  1 μ L of each user sample to the appropriate tube.

Note

Use a P2 pipette for highest accuracy.


20.12 Mix each sample vigorously by vortexing for  00:00:05 and pulse centrifuge to collect the liquid.

20.13 Allow all tubes to incubate at room temperature for  00:02:00 before proceeding.

20.14 On the Home screen of the Quantus Fluorometer, select `Protocol`, then select `ONE DNA` as the assay type.

Note

If you have already performed a calibration for the selected assay you can continue, there is no need to perform repeat calibrations when using ONE DNA pre diluted dye solution. If you want to use the previous calibration, skip to step 11. Otherwise, continue with step 9.


20.15 On the home screen navigate to 'Sample Volume' and set it to  1 μ L then 'Units' and set it to ng/ μ L.

20.16 Load the first sample into the reader and close the lid. The sample concentration is automatically read when you close the lid.

20.17 Repeat step 16 until all samples have been read.

20.18 The value displayed on the screen is the dsDNA concentration in ng/ μ L, carefully record all results in a spreadsheet or laboratory notebook.

MinION sequencing

21 Prime the flowcell and load  30 ng sequencing library onto the flowcell.



Protocol

NAME

Priming and loading a MinION flowcell

CREATED BY

Josh Quick

Preview

Note

From experience we know 30 ng is optimum loading input for short amplicons.

- 21.1 Thaw the following reagents at room temperature before placing on ice:

Sequencing buffer (SQB)

Loading beads (LB)


Flush buffer (FLB)

Flush tether (FLT)

- 21.2 Add  30 μL FLT to the FLB tube and mix well by vortexing.


- 21.3 If required place a new MinION flowcell onto the MinION by flipping open the lip and pushing one end of the flowcell under the clip and pushing down gently.

- 21.4 Rotate the inlet port cover clockwise by 90° so that the priming port is visible.

- 21.5 Take a P1000 pipette and tip and set the volume to  800 μL . Place the tip in the inlet port and holding perpendicularly to the plane of the flowcell remove any air from the inlet port by turning the volume dial anti-clockwise.


Note

Be careful not to remove so much volume that air is introduced onto the rectangular array via the outlet.




21.6 Load  800 µL of FLB (plus FLT) into the flow cell via the inlet port, dispense slowly and smoothly trying to avoid the introduction of any air bubbles.

21.7 Wait for  00:05:00 .

21.8 Gently lift the SpotON cover to open the SpotON port.

21.9 Load another  200 µL of FLB (plus FLT) into the flow cell via the inlet port, this will initiate a siphon at the SpotON port to allow you to load the library dilution.

21.10 In a new tube prepare the library dilution for sequencing:


Component	Volume
SQB	 37.5 µL
LB	 25.5 µL
Final library	 12 µL
Total	 75 µL

Note

Mix LB immediately before use as they settle quickly.

Dilute library in EB if required.

21.11 Mix the prepared library gently by pipetting up and down just prior to loading.

21.12 Add the  75 µL library dilution to the flow cell via the SpotON sample port in a dropwise fashion. Ensure each drop siphons into the port before adding the next.

21.13 Gently replace the SpotON sample port cover, making sure the bung enters the SpotON port, close the inlet port and close the MinION lid.



22 Start the sequencing run using MinKNOW.

Protocol

NAME

Starting a MinION sequencing run using MinKNOW

CREATED BY

Josh Quick

Preview

22.1 If required plug the MinION into the computer and wait for the MinION and flowcell to be detected.

22.2 Choose flow cell 'FLO-MIN106' from the drop-down menu.

22.3 Then select the flowcell so a tick appears.

22.4 Click the 'New Experiment' button in the bottom left of the screen.

22.5 On the New experiment popup screen, select the running parameters for your experiment from the individual tabs:

Experiment: Name the run in the experiment field, leave the sample field blank.

Kit: Selection: Select LSK109 as there is no option for native barcoding (NBD104).

Run Options: Set the run length to 6 hours (you can stop the run once sufficient data has been collected as determined using RAMPART).

Basecalling: Leave basecalling turned but select 'fast basecalling'.



Output: The number of files that MinKNOW will write to a single folder. By default this is set to 4000 but can be reduced to make RAMPART update more frequently.

Click 'Start run'.

22.6 Monitor the progress of the run using the MinKNOW interface.