

Sep 27, 2019

# Ebola virus sequencing protocol

DOI

[dx.doi.org/10.17504/protocols.io.7nwhmfe](https://dx.doi.org/10.17504/protocols.io.7nwhmfe)



Josh Quick<sup>1</sup>

<sup>1</sup>University of Birmingham

ARTIC



Josh Quick

University of Birmingham

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.7nwhmfe](https://dx.doi.org/10.17504/protocols.io.7nwhmfe)

**Protocol Citation:** Josh Quick 2019. Ebola virus sequencing protocol. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.7nwhmfe>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** September 25, 2019





**Last Modified:** September 27, 2019

**Protocol Integer ID:** 28086



## cDNA preparation

- 1 Mix the following components in an 0.2mL 8-strip tube;



Component	Volume
50μM random hexamers	 1 μL
10mM dNTPs mix (10mM each)	 1 μL
Template RNA	 10 μL
<b>Total</b>	 12 μL


### Note

Viral RNA input from a clinical sample should be between Ct 18-35. If Ct is between 12-15, then dilute the sample 100-fold in water, if between 15-18 then dilute 10-fold in water. This will reduce the likelihood of PCR-inhibition.






- 2 Gently mix by pipetting and pulse spin the tube to collect liquid at the bottom of the tube.

- 3 Incubate the reaction as follows:

 65 °C for  00:05:00

Place on ice for  00:01:00

- 4 Add the following to the annealed template RNA:

Component	Volume
SSIV Buffer	 4 μL
100mM DTT	 1 μL
RNaseOUT RNase Inhibitor	 1 μL
SSIV Reverse Transcriptase	 1 μL
<b>Total</b>	 20 μL

- 5 Gently mix by pipetting and pulse spin the tube to collect liquid at the bottom of the tube.



6 Incubate the reaction as follows:

42 °C 01:30:00

70 °C 00:10:00

Hold at 5 °C

## Primer pool preparation

7 If required resuspend lyophilised primers at a concentration of 100µM each

### Note

**Ebola V2** primers for this protocol were designed using **Primal Scheme** and generate overlapping 400nt amplicons. Primer names and dilutions are listed in the table below.

8 Generate primer pool stocks by adding 5 µL of each primer pair to a 1.5 mL Eppendorf labelled either "Pool 1 (100µM)" or "Pool 2 (100µM)". Total volume should be 505 µL for Pool 1 (100µM) and 530 µL for Pool 2 (100µM). These are your 100µM stocks of each primer pool.

### Note

Primers should be prepped and aliquoted prior to departure in a sterile PCR cabinet. At no stage should primers or PCR reagents be anywhere near the template until cDNA addition.

9 Dilute this primer pool 1:10 in molecular grade water, to generate 10µM primer stocks. It is recommend that multiple aliquots of each primer pool are made to in case of degradation or contamination.

<b>Name</b>	<b>Sequence</b>	<b>Name</b>	<b>Sequence</b>	<b>Pool</b>	<b>Stock</b>
Ebov-10-Pan_1_LEF_T	TGTGTGCGAATAAC TATGAGGAAGA	Ebov-10-Pan_1_RIGHT	TTTCCAATGTTTTAC CCCAAGCTTT	1	100μM
		Ebov-10-Pan_1_RIGHT_alt1	TTTCCAATGCTTTAC CCCAAGCTTT	1	100μM
		Ebov-10-Pan_1_RIGHT_alt2	TTTCCAATGTTTTAC CCCAAGTTTT	1	100μM
Ebov-10-Pan_2_LEF_T	CAAGCAAGATTGA GAATTAACCTTGGT	Ebov-10-Pan_2_RIGHT	ATCTCCCTGGTACG CATGATGA	2	100μM
Ebov-10-Pan_2_LEF_T_alt1	CAAGCAAGATTGA GAATTAACCTTGAT	Ebov-10-Pan_2_RIGHT_alt1	ATCTCCTTGGTACG CATGATGA	2	100μM
Ebov-10-Pan_3_LEF_T	GGCCTTTGAAGCA GGTGTGAT	Ebov-10-Pan_3_RIGHT	TCAGTCCTTGCTCT GCATGTAC	1	100μM
Ebov-10-Pan_4_LEF_T	CCTTTGCAAGTCTA TTCCTTCCGA	Ebov-10-Pan_4_RIGHT	CTGAGTGCAGCCTT AAAGGAGT	2	100μM
Ebov-10-Pan_4_LEF_T_alt1	CTTTTGCAAGTCTA TTCCTTCCGA			2	100μM
Ebov-10-Pan_5_LEF_T	AGTTCGTCTCCATC CTCTTGCA	Ebov-10-Pan_5_RIGHT	CTGGAAGCTGATTT CGTTCTTTTTCT	1	100μM
Ebov-10-Pan_6_LEF_T	GAGTCTCGCGAAC TTGACCATC	Ebov-10-Pan_6_RIGHT	TCCTCGTCGTCCTC GTCTAGAT	2	100μM
Ebov-10-Pan_6_LEF_T_alt1	GAATCTCGCGAAC TTGACCATC	Ebov-10-Pan_6_RIGHT_alt1	TCCTCATCGTCCTC GTCTAGAT	2	100μM
Ebov-10-	AGCTACGGCGAAT ACCAGAGTT	Ebov-10-Pan_7_RIGHT	GTCCCTGTCCTGCT CTTCATCA	1	100μM



Pan_7 _LEF T					
		Ebov-10- Pan_7_RIGHT_alt 1	GTCCCTGTCCTGTT CTTCATCA	1	100μ M
		Ebov-10- Pan_7_RIGHT_alt 2	GTCCCTGTCCTGTT CTTCATCG	1	100μ M
Ebov- 10- Pan_8 _LEF T	TTAACGAAGAGGC AGACCCACT	Ebov-10- Pan_8_RIGHT	TTCCTCTTCAAGGG AGTCTGGA	2	100μ M
Ebov- 10- Pan_8 _LEF T_alt1	TCAACGAAGAGGC AGACCCACT	Ebov-10- Pan_8_RIGHT_al t1	TTCCTCTTCAAGGG AGTCCGGA	2	100μ M
Ebov- 10- Pan_9 _LEF T	GTGACAACACCCA GTCAGAACAA	Ebov-10- Pan_9_RIGHT	TCTTCCTGTTTTCGT TCCTTGACT	1	100μ M
Ebov- 10- Pan_9 _LEF T_alt1	GTGACAACACCCA GCCAGAACAA	Ebov-10- Pan_9_RIGHT_al t1	TCTTCCTGTTTTGCG TTCCTTGACT	1	100μ M
		Ebov-10- Pan_9_RIGHT_al t2	TCTTCCTGTTTTGCG TTTCTTGACT	1	100μ M
Ebov- 10- Pan_1 0_LEF T	ACAATGGGATGATT CAACCGACA	Ebov-10- Pan_10_RIGHT	TCGAGTGCTAGAGA ATTCAATTGACG	2	100μ M
Ebov- 10- Pan_1 0_LEF T_alt1	ATAATGGGATGATT TAACCGACA			2	100μ M
Ebov- 10- Pan_1 1_LEF T	ACCTACTAGCCTG CCCAACATT	Ebov-10- Pan_11_RIGHT	AATTGGGTCCGTTT GGGTTTGA	1	100μ M
Ebov- 10- Pan_1 1_LEF T_alt1	ACCTACTAGCCTAC CCAACATT	Ebov-10- Pan_11_RIGHT_al t1	AATTGGATCCGTTT GGGTTTGA	1	100μ M
Ebov- 10- Pan_1	CCCAAATGCAACA AACGAAGCC	Ebov-10- Pan_12_RIGHT	TCAATCTTACCCCG AATCGCAC	2	100μ M



2_LEF T					
Ebov-10- Pan_1 2_LEF T_alt1	CCCAAATGCAACA AACAAAGCC	Ebov-10- Pan_12_RIGHT_a It1	TCAATCTTACCCCG AATTGCAC	2	100μ M
Ebov-10- Pan_1 3_LE FT	TATTGGGCCGAAC ATGGTCAAC	Ebov-10- Pan_13_RIGHT	TGACAGGTGGAGCA GCATCTTG	1	100μ M
Ebov-10- Pan_1 3_LE FT_alt 1	TATTGGGCTGAACA TGGTCAAC			1	100μ M
Ebov-10- Pan_1 4_LE FT	CATTCATGCTGAGT TCCAGGCC	Ebov-10- Pan_14_RIGHT	GCGAGATATGAACA ATTTTATCTTGGTCG	2	100μ M
		Ebov-10- Pan_14_RIGHT_a It1	GCGAGATAAGGACA ATTTTATCTTGGTCG	2	100μ M
		Ebov-10- Pan_14_RIGHT_a It2	GCGAGATAAGAACA ATTTTATCTTGGTCG	2	100μ M
Ebov-10- Pan_1 5_LEF T	TGAGTATCAGCCCT GGATAATATAAGTC A	Ebov-10- Pan_15_RIGHT	TCGATGGAGTGTCC CCATTGAC	1	100μ M
Ebov-10- Pan_1 5_LEF T_alt1	TGAGTATCAGCCCT AGATAATATAAGTC A	Ebov-10- Pan_15_RIGHT_a It1	TCGATGGAGTGTCT CCATTGAC	1	100μ M
Ebov-10- Pan_1 6_LEF T	GCAACAGCAATACA GGCTTCCT	Ebov-10- Pan_16_RIGHT	GAAAGCCTGGTTTC CAATTCGC	2	100μ M
Ebov-10- Pan_1 6_LEF T_alt1	GCAACAACAATACA GGCTTCCT	Ebov-10- Pan_16_RIGHT_a It1	GAAGGCCTGGTTTC CAATTCGC	2	100μ M
Ebov-10- Pan_1 7_LEF T	CCACTTGTCAGAGT CAATCGGC	Ebov-10- Pan_17_RIGHT	GTTTCTGGCACTTC GATTCCCA	1	100μ M



		Ebov-10-Pan_17_RIGHT_alt1	GTTTCTGGCACTTC GATACCCA	1	100μM
Ebov-10-Pan_18_LEFT	AAAATCCAAGCAAT AATGACTTCACTCC	Ebov-10-Pan_18_RIGHT	TTGATCAATTAAAAG TGTCTCCTCTAATG G	2	100μM
		Ebov-10-Pan_18_RIGHT_alt1	TCGATCAATTTAAAG TATCTCCTCTAATGG	2	100μM
		Ebov-10-Pan_18_RIGHT_alt2	TTGATCAATTAAAAG TATCTCCTCTAATAG	2	100μM
Ebov-10-Pan_19_LEFT	AGATCCAGTTTTAT AGAATCTTCTCAGG GA	Ebov-10-Pan_19_RIGHT	AGAAGGGCAATGTC TGTACTTGG	1	100μM
Ebov-10-Pan_19_LEFT_alt1	AGATCCAGTTTTAC AGAATCTTCTCAGG GA	Ebov-10-Pan_19_RIGHT_alt1	AGAAGGGCGATGTC TGTGCTTGG	1	100μM
Ebov-10-Pan_20_LEFT	AGCCAGTGTGACTT GGATTGGA	Ebov-10-Pan_20_RIGHT	AGTTTGTGCGACATCA CTAACCTGT	2	100μM
		Ebov-10-Pan_20_RIGHT_alt1	AGTTTGTGCGACATCA CTAACTTGT	2	100μM
Ebov-10-Pan_21_LEFT	AGAACATTTTCCAT CCCACTTGGA	Ebov-10-Pan_21_RIGHT	AAGCACCTCTTTA TGGAAGGC	1	100μM
		Ebov-10-Pan_21_RIGHT_alt1	AAGCACCTCTTTG TGGAAGGC	1	100μM
Ebov-10-Pan_22_LEFT	TGCCGGTATGTGCA CAAAGTAT	Ebov-10-Pan_22_RIGHT	ATATATTGTCTCATT CAGCTGGAGCA	2	100μM
Ebov-10-Pan_23_LEFT	CGAGGTTGACAATT TGACCTACGT	Ebov-10-Pan_23_RIGHT	GCAAGGGTTGTTAG ATGCGACA	1	100μM
		Ebov-10-Pan_23_RIGHT_alt1	GCAAGGGTTGTCAG ATGCGACA	1	100μM



Ebov-10-Pan_24_LE FT	TGCAATGGTTCAAG TGCACAGT	Ebov-10-Pan_24_RIGHT	CTGGCACTCTCTTC TCCGGTAT	2	100μM
Ebov-10-Pan_24_LE FT_alt1	TGCAATGGTTCAAG TGCACAAT			2	100μM
Ebov-10-Pan_25_LE FT	ACCACAACAAGTC CCCCAAAACC	Ebov-10-Pan_25_RIGHT	TAGCTCAGTTGTGG CTCTCAGG	1	100μM
		Ebov-10-Pan_25_RIGHT_alt1	TAGCTCGGTTGTGG CTCTCAGG	1	100μM
Ebov-10-Pan_26_LE FT	ATCTGTGGGTTGAG ACAGCTGG	Ebov-10-Pan_26_RIGHT	GCTTTTCCATGAAG CAATCTGAAGA	2	100μM
Ebov-10-Pan_26_LE FT_alt1	ATCTGTGGATTGAG GCAGCTGG	Ebov-10-Pan_26_RIGHT_alt1	GCTTTGCCATGAAG CAATCTGAAGA	2	100μM
Ebov-10-Pan_26_LE FT_alt2	ATCTGTGGGTTGAG GCAGCTGG			2	100μM
Ebov-10-Pan_27_LE FT	TGGAGTTACAGGC GTTATAATTGCA	Ebov-10-Pan_27_RIGHT	AAAGGCTTCTTTCC CTTGTCAT	1	100μM
Ebov-10-Pan_28_LE FT	TCATCCTTGATTCT ACAATCATGACAGT	Ebov-10-Pan_28_RIGHT	AGGTGCTGGAGGAA CTGTTAATG	2	100μM
Ebov-10-Pan_28_LE FT_alt1	TCATCCTTGATTCT ACAATCATAACAGT			2	100μM
Ebov-10-Pan_29_LE FT	GAGTACCGTCAATC AAGGAGCG	Ebov-10-Pan_29_RIGHT	CACAGCACATAGAG TCAACAATGC	1	100μM



Ebov-10-Pan_30_LEF T	GATCAAGACGGCA GAACACTGG	Ebov-10-Pan_30_RIGHT	ATCAGACCATGAGC ATGTCCCC	2	100μM
Ebov-10-Pan_31_LEF T	CTGCTGTCGTTGTT TCAGGGTT	Ebov-10-Pan_31_RIGHT	ATGGGATGGATCGT TGCTACCT	1	100μM
		Ebov-10-Pan_31_RIGHT_a lt1	ATGGGATGGATCGT TGCTGCCT	1	100μM
		Ebov-10-Pan_31_RIGHT_a lt2	ATGAGATGGATCGTT GCTACCT	1	100μM
Ebov-10-Pan_32_LEF T	GCCAAGCATACCT CTTGACAA	Ebov-10-Pan_32_RIGHT	TGGACTIONCTGAA ATAGTACTTTGC	2	100μM
Ebov-10-Pan_33_LE FT	TGCGGAGGTCTGA TAAGAATAAACC	Ebov-10-Pan_33_RIGHT	TTCAACCTTGAAAC CTTGCGCT	1	100μM
		Ebov-10-Pan_33_RIGHT_ alt1	TTCAACTTTGAAAC CTTGCGCT	1	100μM
Ebov-10-Pan_34_LE FT	GCTGAAAAGAAGC TTACCTACAACG	Ebov-10-Pan_34_RIGHT	TCCTTGTCATTGAC CATGCAGG	2	100μM
Ebov-10-Pan_34_LE FT_alt 1	GTTGAAAAAAGGC CTACCTACAACG			2	100μM
Ebov-10-Pan_34_LE FT_alt 2	GCTGAAAAGAAGC CCACCTACAACG			2	100μM
Ebov-10-Pan_35_LEF T	GTGACTCACAAAG GAATGGCCC	Ebov-10-Pan_35_RIGHT	ACAATCCGTTGTAGT TCACGACA	1	100μM
		Ebov-10-Pan_35_RIGHT_ alt1	ACAACCCGTTGTAG TTCACGACA	1	100μM



Ebov-10-Pan_36_LEF T	TGCTGTCGTTGATT CGATCCAA	Ebov-10-Pan_36_RIGHT	AGCAGAGATGTCAA GATAACTATTGAGT	2	100μ M
Ebov-10-Pan_37_LEF T	ACACGAATGCAAA GTTTGATTCTTGA	Ebov-10-Pan_37_RIGHT	TGAAACCTAACACA TGTGACCTGC	1	100μ M
		Ebov-10-Pan_37_RIGHT_a It1	TGAAACCTAACACA CGTGACCTGC	1	100μ M
Ebov-10-Pan_38_LE FT	CCCTCAAACAAGA GATTCCAAGACA	Ebov-10-Pan_38_RIGHT	ACAGTTGCGTAGTT GCGGATTA	2	100μ M
Ebov-10-Pan_38_LE FT_alt 1	CCCTCAAATAAGA GATTCCAAGACA			2	100μ M
Ebov-10-Pan_38_LE FT_alt 2	TCCTCAAATAAGAG ATTCCAAGACA			2	100μ M
Ebov-10-Pan_39_LEF T	ACCTAGTCACTAGA GCTTGCGG	Ebov-10-Pan_39_RIGHT	ACATTTGATGTAAAA ATTCATTGCCCTG	1	100μ M
Ebov-10-Pan_40_LEF T	GTGGGTGCTCAAG AAGACTGTG	Ebov-10-Pan_40_RIGHT	TGAGATTAGAGTTGT GTTGAATCGACA	2	100μ M
Ebov-10-Pan_40_LEF T_alt1	GTGGGTGCTCAAG AGGACTGTG	Ebov-10-Pan_40_RIGHT_ alt1	TGAGATTAGAGTCGT GTTGAATCGACA	2	100μ M
Ebov-10-Pan_41_LEF T	AAGAAGCGGTTCA AGGGCATAC	Ebov-10-Pan_41_RIGHT	CTATGGAATTCACG GATCTTTTGAGC	1	100μ M
Ebov-10-Pan_41_LEF T_alt1	AAGAAGCAGTTCA AGGGCATAC	Ebov-10-Pan_41_RIGHT_a It1	CTATGGAATTCACG GATCTTTTGATC	1	100μ M

Ebov-10-Pan_42_LEF T	TGCATTTAGCTGTA AATCACACCCT	Ebov-10-Pan_42_RIGHT	AATCATTGGCAACG GAGGGAAT	2	100μM
		Ebov-10-Pan_42_RIGHT_alt1	AATCATTGGCAACG GGGGAAT	2	100μM
Ebov-10-Pan_43_LE FT	GTCAAGGATCTTG GTACAGTGTTACT	Ebov-10-Pan_43_RIGHT	TGAGAAAGAAAAGT TCCGATATTGTGGT	1	100μM
Ebov-10-Pan_43_LE FT_alt 1	GCCAAGGGTCTTG GTACAGTGTTACT	Ebov-10-Pan_43_RIGHT_alt1	TGAGAAAGAAAAAT TCCGGTATTGTGGT	1	100μM
Ebov-10-Pan_43_LE FT_alt 2	GTCAAGGGTCTTG GTACAGTGTTACT	Ebov-10-Pan_43_RIGHT_alt2	TGAGAAAGAAAAAT TCCGATATTGTGGT	1	100μM
Ebov-10-Pan_44_LE FT	TTGAGAATGTTCTT TCCTACGCACA	Ebov-10-Pan_44_RIGHT	ACGGTTGCAATATTC TATAAAAGGTGC	2	100μM
Ebov-10-Pan_44_LE FT_alt 1	TTGAGAATGTTCTT TCCTACGCGCA	Ebov-10-Pan_44_RIGHT_alt1	ACGGTTGCAATATTC GATAAAAGGTGC	2	100μM
		Ebov-10-Pan_44_RIGHT_alt2	ACGGTTACAATATTC TATAAAAGGTGC	2	100μM
Ebov-10-Pan_45_LEF T	CCACAGTTAGAGG GAGTAGCTTTG	Ebov-10-Pan_45_RIGHT	GCTCGTCTGCGTCA GTCTCTAA	1	100μM
Ebov-10-Pan_45_LEF T_alt1	CCACAGTTAGAGG GAGTAGTTTTG			1	100μM
Ebov-10-Pan_46_LEF T	AAGTTACGCTCAG CTGTGATGG	Ebov-10-Pan_46_RIGHT	ATGGAAAGCTGCGG TTATCCTG	2	100μM
Ebov-10-	TAGGCACTGCTTTT GAGCGATC	Ebov-10-Pan_47_RIGHT	CACAAAGTCAATGG CAGTGCAG	1	100μM



Pan_4 7_LEF T					
Ebov-10- Pan_4 7_LEF T_alt1	TAGGCACCGCTTT TGAGCGGTC			1	100μ M
Ebov-10- Pan_4 7_LEF T_alt2	TAGGCACTGCTTTT GAACGATC			1	100μ M
Ebov-10- Pan_4 8_LE FT	TCTCCGAATGATTG AGATGGATGATT	Ebov-10- Pan_48_RIGHT	CTCAGTCTGTCCAA AACCGGTG	2	100μ M
Ebov-10- Pan_4 8_LE FT_alt 1	TCTCCGAATGATTG GGATGGATGATT			2	100μ M
Ebov-10- Pan_4 9_LEF T	GATATCTTTTCACG CACGCCGA	Ebov-10- Pan_49_RIGHT	CCACCTGGTTGCTT TGCATTG	1	100μ M
Ebov-10- Pan_4 9_LEF T_alt1	GATATCTTTTCACG CACGCCCA	Ebov-10- Pan_49_RIGHT_ alt1	CCACCAGGTTGCTT TGCATTG	1	100μ M
Ebov-10- Pan_5 0_LEF T	TCAAAGTGTTTTGG CTGAAACCCT	Ebov-10- Pan_50_RIGHT	TCCTGAGTAATGTGA AGGGGTCA	2	100μ M
Ebov-10- Pan_5 0_LEF T_alt1	TCAAAGTGTTTTG GCTGAAACCCT	Ebov-10- Pan_50_RIGHT_ alt1	TCCTGAGTAATGTGA AGGAGTCA	2	100μ M
Ebov-10- Pan_5 1_LEF T	AACAGTGA CTTGCT AATAAAACCATTTT TG	Ebov-10- Pan_51_RIGHT	AAATACTGAGCTGG TACTTCCCG	1	100μ M
Ebov-10- Pan_5 1_LEF T_alt1	AACAGTGA CTTGCT AATAAGCCATTTT TG			1	100μ M
Ebov-10-	AACAGTGATTTGCT AATAAAACCATTTT			1	100μ M



Pan_5 1_LEF T_alt2	TG				
Ebov-10- Pan_5 2_LEF T	AATCGTGCTCACCT TCATCTAACT	Ebov-10- Pan_52_RIGHT	CCCCAAACTGTACA GAAGTCCTATCT	2	100μ M
Ebov-10- Pan_5 3_LE FT	ACAGACCCAATTAG CAGTGGAGA	Ebov-10- Pan_53_RIGHT	ACAATTGTTCCGCG ATTAATTATCCAT	1	100μ M
Ebov-10- Pan_5 3_LE FT_alt 1	ACAGACCCAATTAG CAGCGGAGA	Ebov-10- Pan_53_RIGHT_ alt1	ACAATTGTTCCGCG ATTAATTATCCACT	1	100μ M
Ebov-10- Pan_5 4_LE FT	TCTCAGATGCGGC CAGGTTATT	Ebov-10- Pan_54_RIGHT	TGACCATCACTGTT GTTTGTGCT	2	100μ M
Ebov-10- Pan_5 4_LE FT_alt 1	TCTCAGATGCGGC CAGATTATT			2	100μ M
Ebov-10- Pan_5 5_LEF T	TGGAGGAGCAGAC ACAGAAACA	Ebov-10- Pan_55_RIGHT	ATGACGTTAATTGGC GTGTCCC	1	100μ M
Ebov-10- Pan_5 5_LEF T_alt1	TGGAGGAGCAGGC ACAGAAACA	Ebov-10- Pan_55_RIGHT_ alt1	ATGACGTCAATTGG CGTGTCCC	1	100μ M
Ebov-10- Pan_5 5_LEF T_alt2	TGGAGAAGCAGGC ACAGAAACA	Ebov-10- Pan_55_RIGHT_ alt2	ATGACGTTAATTGGC GCGTCCC	1	100μ M
Ebov-10- Pan_5 6_LEF T	CTCACACCGTCTAG TCCTACCT	Ebov-10- Pan_56_RIGHT	TTTGACATAACAGGT AGAAGCATCCT	2	100μ M
Ebov-10- Pan_5 6_LEF T_alt1	CTCGCACCGTCTA GTCCTACCT			2	100μ M



Ebov-10-Pan_56_LEF T_alt2	CTCACATCGTCTAG TCCTACCT			2	100μM
Ebov-10-Pan_57_LEF T	ACACGCTAGCTACT GAGTCCAG	Ebov-10-Pan_57_RIGHT	ATTGGCTTAATTAAA TAACCAGTGGCA	1	100μM
Ebov-10-Pan_58_LEF T	TGAAAGCAGTGGT CCTTAAAGTCT	Ebov-10-Pan_58_RIGHT	TGCTCTAAGATGTG CTAAGTGCTG	2	100μM
		Ebov-10-Pan_58_RIGHT_alt1	TGCTCTAAGATGTG CCAAGTGCTG	2	100μM
Ebov-10-Pan_59_LEF T	CGTCGATTCAAAA AGAGGTCCACT	Ebov-10-Pan_59_RIGHT	TCAGAAGCCCTGTC AGCCTTTC	1	100μM
Ebov-10-Pan_60_LEF T	AGATTGCAATTGTG AAGAACGTTTCT	Ebov-10-Pan_60_RIGHT	AGAGTGCAGAGTTT ATTATGTTGCGT	2	100μM
Ebov-10-Pan_61_LEF T	TCACAATGCAGCAT GTGTGACA	Ebov-10-Pan_61_RIGHT	AGGTATTTCTGATTT TACAGTCCTGCC	1	100μM
		Ebov-10-Pan_61_RIGHT_alt1	AGGTATTTATGATTT TACAGTCCTGCC	1	100μM
		Ebov-10-Pan_61_RIGHT_alt2	AGGTATTTCTGATTT TACAGTCATGCC	1	100μM
Ebov-10-Pan_62_LEF T	CCTGTCAGATGGA ATAGTGT TTTTGGT	Ebov-10-Pan_62_RIGHT	AATTTTTGTGTGCGA CCATTTTTCC	2	100μM









**Note**

Primers need to be used at a final concentration of 0.015μM per primer. In this case, Pool 1 has 101 primers in it so the requirement is 3.8μL of Pool 1 (10μM) primers per 25μL reaction. Pool 2 has 106 primers the requirement is 4.0μL of Pool 2 (10μM) primers per 25μL reaction. For other schemes, adjust the volume added appropriately.



## Multiplex PCR

- 10 In the mastermix hood set up the multiplex PCR reactions as follows in 0.2mL 8-strip PCR tubes:

Component	Pool 1	Pool 2
NEB Q5 Polymerase 2X MasterMix	 12.5 $\mu$ L	 12.5 $\mu$ L
Primer Pool 1 or 2 (10 $\mu$ M)	 3.8 $\mu$ L	 4.0 $\mu$ L
Water	 6.2 $\mu$ L	 6.0 $\mu$ L
<b>Total</b>	 22.5 $\mu$ L	 22.5 $\mu$ L

### Note

The **mastermix hood** should be prepared by sterilising with UV and treated with MediPal wipes, DNaseZap and RNaseZap reagents. Wipe down the hood with each sequentially, allowing 5 minutes for drying between each. Pipettes should also be treated in the same way, and UV treated for 30 mins between library preparations.

This step should be carried out in the mastermix hood and template should not be taken anywhere near the mastermix hood at any stage.

- 11 In the template hood add  2.5  $\mu$ L cDNA to each tube and mix well by pipetting.

### Note

The **template hood** should be prepared by sterilising with UV and treated with MediPal wipes, DNaseZap and RNaseZap reagents. Wipe down the hood with each sequentially, allowing 5 minutes for drying between each. Pipettes should also be treated in the same way, and UV treated for 30 mins between library preparations.

- 12 Pulse centrifuge the tubes to collect the contents at the bottom of the tube.

- 13 Set-up the following program on the thermal cycler:

Step	Temperature	Time	Cycles
------	-------------	------	--------

Heat Activation	98 °C	00:00:30	1
Denaturation	98 °C	00:00:15	25-35
Annealing	65 °C	00:05:00	25-35
Hold	4 °C	Indefinite	1

#### Note

Cycle number should be 25 for Ct 18-21 up to a maximum of 35 cycles for Ct 35

## PCR clean-up

14 Combine the entire contents of "Pool 1" and "Pool 2" PCR reactions for each biological sample into to a single 1.5 mL Eppendorf tube.

15 Clean-up the amplicons using the following protocol:

#### Protocol



NAME

**Amplicon clean-up using SPRI beads**

CREATED BY

Josh Quick

**PREVIEW**

15.1 Vortex SPRI beads thoroughly to ensure they are well resuspended, the solution should be a homogenous brown colour.













Agencourt AMPure XP **Beckman Coulter Catalog #A63880**


15.2 Add an equal volume (1:1) of SPRI beads to the sample tube and mix gently by either flicking or pipetting. For example add 50 µL SPRI beads to a 50 µL reaction.

15.3 Pulse centrifuge to collect all liquid at the bottom of the tube.





- 15.4 Incubate for  00:05:00 at room temperature.
- 15.5 Place on magnetic rack and incubate for  00:02:00 or until the beads have pelleted and the supernatant is completely clear.
- 15.6 Carefully remove and discard the supernatant, being careful not to touch the bead pellet.
- 15.7 Add  200  $\mu\text{L}$  of room-temperature  70 % volume ethanol to the pellet.
- 15.8 Carefully remove and discard ethanol, being careful not to touch the bead pellet.
- 15.9  [go to step #15.7](#) and repeat ethanol wash.
- 15.10 Pulse centrifuge to collect all liquid at the bottom of the tube and carefully remove as much residual ethanol as possible using a P10 pipette.
- 15.11 With the tube lid open incubate for  00:01:00 or until the pellet loses its shine (if the pellet dries completely it will crack and become difficult to resuspend).
- 15.12 Resuspend pellet in  30  $\mu\text{L}$  Elution Buffer (EB), mix gently by either flicking or pipetting and incubate for  00:02:00 .  
 [Elution Buffer \(EB\) Qiagen Catalog #19086](#)
- 15.13 Place on magnet and transfer sample to a clean 1.5mL Eppendorf tube ensuring no beads are transferred into this tube.
- 15.14 Quantify  1  $\mu\text{L}$  product using the Quantus Fluorometer using the ONE dsDNA assay.

 QuantiFluor(R) ONE dsDNA System, 100rxn **Promega Catalog #E4871**

### Equipment

Quantus	NAME
Fluorometer	TYPE
Promega	BRAND
E6150	SKU
<a href="https://www.promega.co.uk/products/microplate-readers-fluorometers-luminometers/fluorometers/quantus-fluorometer">https://www.promega.co.uk/products/microplate-readers-fluorometers-luminometers/fluorometers/quantus-fluorometer</a>	LINK

## Quantification and normalisation

16 Quantify the amplicon pools using the Quantus Fluorometer using the ONE dsDNA assay.

### Protocol




NAME







**DNA quantification using the Quantus fluorometer**

CREATED BY  
Josh Quick

**PREVIEW**



16.1 Remove Lambda DNA 400 ng/μL standard from the freezer and leave on ice to thaw. Remove ONE dsDNA dye solution from the fridge and allow to come to room temperature.

 QuantiFluor(R) ONE dsDNA System, 500rxn **Promega Catalog #E4870**

- 16.2 Set up two  0.5 mL tubes for the calibration and label them 'Blank' and 'Standard'
- 16.3 Add  200  $\mu$ L ONE dsDNA Dye solution to each tube.
- 16.4 Mix the Lambda DNA standard 400 ng/ $\mu$ L standard by pipetting then add  1  $\mu$ L to one of the standard tube.
- 16.5 Mix each sample vigorously by vortexing for  00:00:05 and pulse centrifuge to collect the liquid.
- 16.6 Allow both tubes to incubate at room temperature for  00:02:00 before proceeding.
- 16.7 Selection 'Calibrate' then 'ONE DNA' then place the blank sample in the reader then select 'Read Blank'. Now place the standard in the reader and select 'Read Std'.
- 16.8 Set up the required number of  0.5 mL tubes for the number of DNA samples to be quantified.


#### Note

Use only thin-wall, clear, 0.5mL PCR tubes such as Axygen #PCR-05-C

- 16.9 Label the tubes on the lids, avoid marking the sides of the tube as this could interfere with the sample reading.
- 16.10 Add  199  $\mu$ L ONE dsDNA dye solution to each tube.
- 16.11 Add  1  $\mu$ L of each user sample to the appropriate tube.

**Note**

Use a P2 pipette for highest accuracy.


16.12 Mix each sample vigorously by vortexing for  00:00:05 and pulse centrifuge to collect the liquid.

16.13 Allow all tubes to incubate at room temperature for  00:02:00 before proceeding.

16.14 On the Home screen of the Quantus Fluorometer, select 'Protocol', then select 'ONE DNA' as the assay type.

**Note**

If you have already performed a calibration for the selected assay you can continue, there is no need to perform repeat calibrations when using ONE DNA pre diluted dye solution. If you want to use the previous calibration, skip to step 11. Otherwise, continue with step 9.

16.15 On the home screen navigate to 'Sample Volume' and set it to  1  $\mu$ L then 'Units' and set it to ng/ $\mu$ L.

16.16 Load the first sample into the reader and close the lid. The sample concentration is automatically read when you close the lid.



16.17 Repeat step 16 until all samples have been read.

16.18 The value displayed on the screen is the dsDNA concentration in ng/ $\mu$ L, carefully record all results in a spreadsheet or laboratory notebook.

17 Label a  1.5 mL Eppendorf tube for each sample.

**Note**

This is a 'one-pot ligation' protocol for native barcoded ligation libraries. We have seen no reduction in performance compared to standard libraries, and is made faster by using the Ultra II® ligation module which is compatible with the Ultra II® end repair/dA-tailing module removing a clean-up step.

- 18 Normalise the input but adding  10 ng amplicon pools to each tube and diluting to  10 µL to have an input concentration of [M] 1 ng/µL

**Note**

Quantity of amplicons can vary from 10-50ng depending on the amplicon length, any more than this and the molarity of DNA ends will be too high for efficient barcoding. You need to have 6 samples per native barcoded library to have sufficient material at the end.

## Native barcoding

- 19 Barcode the amplicon pools using native barcodes.

**Protocol**

NAME

**One-pot native barcoding of amplicons**

CREATED BY

Josh Quick

**PREVIEW**


- 19.1 Set up the following reaction for each sample:

**Component****Volume**

DNA amplicons

 10 µL

Ultra II End Prep Reaction Buffer

 1.4 µL

Ultra II End Prep Enzyme Mix

 0.6 µL

**Total**

12 µL

19.2 Incubate at room temperature for 00:10:00

Incubate at 65 °C for 00:05:00

Incubate on ice for 00:01:00

19.3 Add the following directly to the previous reactions:

**Component****Volume**

NBXX barcode

2.5 µL

Ultra II Ligation Master Mix

14.5 µL

Ligation Enhancer

0.5 µL

**Total**

29.5 µL

**Note**

Use one native barcode from the EXP-NBD104 (1-12) or EXP-NBD114 (13-24) per sample. Use from 6 to 24 barcodes in a library, any fewer and there will be insufficient total material to achieve good yields.

19.4 Incubate at room temperature for 00:20:00

Incubate at 70 °C for 00:10:00

Incubate on ice for 00:01:00

**Note**

The 70°C incubation is to inactivate the DNA ligase to prevent barcode cross-ligation when reactions are pooled in the next step.

19.5 Pool all barcoded fragments together into a new 1.5 ml Eppendorf tube.

## Protocol



NAME

### Amplicon clean-up using SPRI beads

CREATED BY

Josh Quick

**PREVIEW**

- 19.6 Quantify the barcoded amplicon pools using the Quantus Fluorometer using the ONE dsDNA assay.

## Protocol



NAME

### DNA quantification using the Quantus fluorometer

CREATED BY

Josh Quick

**PREVIEW**

- 19.7 Set up the following AMII adapter ligation reaction:

#### Component

#### Volume

Barcoded amplicon pools

30 µL

NEBNext Quick Ligation Reaction Buffer (5X)

10 µL

AMII adapter mix

5 µL

Quick T4 DNA Ligase

5 µL


**Total**

50 µL

**Note**

The input of barcoded amplicon pools will depend on the number of barcoded pools and should be between 50 ng (6 barcodes) and 200 ng (24 barcodes).


19.8 Incubate at room temperature for  00:20:00


19.9 Add  50  $\mu\text{L}$  (1:1) of SPRI beads to the sample tube and mix gently by either flicking or pipetting.

**Note**

Vortex SPRI beads thoroughly before use to ensure they are well resuspended, the solution should be a homogenous brown colour.

19.10 Pulse centrifuge to collect all liquid at the bottom of the tube.

19.11 Incubate for  00:05:00 at room temperature.

19.12 Place on magnetic rack and incubate for  00:02:00 or until the beads have pelleted and the supernatant is completely clear.

19.13 Carefully remove and discard the supernatant, being careful not to touch the bead pellet.

19.14 Add  200  $\mu\text{L}$  SFB and resuspend beads completely by pipette mixing.

**Note**

SFB will remove excess adapter without damaging the adapter-protein complexes. Do not use 70% ethanol as in early clean-ups.

19.15 Pulse centrifuge to collect all liquid at the bottom of the tube.




19.16 Remove supernatant and discard.

19.17 Repeat steps 14–16 to perform a second SFB wash.

19.18 Pulse centrifuge and remove any residual SFB.

#### Note

You do not need to allow to air dry with SFB washes.

19.19 Add  15  $\mu$ L EB and resuspend beads by pipette mixing.

19.20 Incubate at room temperature for  00:02:00 .

19.21 Place on magnetic rack.

19.22 Transfer final library to a new 1.5mL Eppendorf tube.

20 Quantify the final library using the Quantus Fluorometer using the ONE dsDNA assay.

#### Protocol



NAME

**DNA quantification using the Quantus fluorometer**

CREATED BY


Josh Quick


**PREVIEW**


**Note**


Final library can now be stored in EB at 4°C for up to a week if needed otherwise proceed directly to MinION sequencing.


- 20.1 Remove Lambda DNA 400 ng/μL standard from the freezer and leave on ice to thaw. Remove ONE dsDNA dye solution from the fridge and allow to come to room temperature.

 QuantiFluor(R) ONE dsDNA System, 500rxn **Promega Catalog #E4870**

- 20.2 Set up two  0.5 mL tubes for the calibration and label them 'Blank' and 'Standard'


- 20.3 Add  200 μL ONE dsDNA Dye solution to each tube.

- 20.4 Mix the Lambda DNA standard 400 ng/μL standard by pipetting then add  1 μL to one of the standard tube.

- 20.5 Mix each sample vigorously by vortexing for  00:00:05 and pulse centrifuge to collect the liquid.

- 20.6 Allow both tubes to incubate at room temperature for  00:02:00 before proceeding.

- 20.7 Selection 'Calibrate' then 'ONE DNA' then place the blank sample in the reader then select 'Read Blank'. Now place the standard in the reader and select 'Read Std'.


- 20.8 Set up the required number of  0.5 mL tubes for the number of DNA samples to be quantified.

**Note**

Use only thin-wall, clear, 0.5mL PCR tubes such as Axygen #PCR-05-C


20.9 Label the tubes on the lids, avoid marking the sides of the tube as this could interfere with the sample reading.

20.10 Add  199  $\mu\text{L}$  ONE dsDNA dye solution to each tube.

20.11 Add  1  $\mu\text{L}$  of each user sample to the appropriate tube.

#### Note

Use a P2 pipette for highest accuracy.


20.12 Mix each sample vigorously by vortexing for  00:00:05 and pulse centrifuge to collect the liquid.

20.13 Allow all tubes to incubate at room temperature for  00:02:00 before proceeding.

20.14 On the Home screen of the Quantus Fluorometer, select 'Protocol', then select 'ONE DNA' as the assay type.

#### Note

If you have already performed a calibration for the selected assay you can continue, there is no need to perform repeat calibrations when using ONE DNA pre diluted dye solution. If you want to use the previous calibration, skip to step 11. Otherwise, continue with step 9.

20.15 On the home screen navigate to 'Sample Volume' and set it to  1  $\mu\text{L}$  then 'Units' and set it to ng/ $\mu\text{L}$ .


20.16 Load the first sample into the reader and close the lid. The sample concentration is automatically read when you close the lid.

20.17 Repeat step 16 until all samples have been read.



- 20.18 The value displayed on the screen is the dsDNA concentration in ng/μL, carefully record all results in a spreadsheet or laboratory notebook.

## MinION sequencing

- 21 Prime the flowcell and load  30 ng sequencing library onto the flowcell.

### Protocol



NAME

### Priming and loading a MinION flowcell

CREATED BY

Josh Quick

**PREVIEW**

### Note

From experience we know 30 ng is optimum loading input for short amplicons.

- 21.1 Thaw the following reagents at room temperature before placing on ice:

Sequencing buffer (SQB)

Loading beads (LB)

Flush buffer (FLB)


Flush tether (FLT)

- 21.2 Add  30 μL FLT to the FLB tube and mix well by vortexing.

- 21.3 If required place a new MinION flowcell onto the MinION by flipping open the lip and pushing one end of the flowcell under the clip and pushing down gently.


- 21.4 Rotate the inlet port cover clockwise by 90° so that the priming port is visible.




- 21.5 Take a P1000 pipette and tip and set the volume to  800  $\mu\text{L}$  . Place the tip in the inlet port and holding perpendicularly to the plane of the flowcell remove any air from the inlet port by turning the volume dial anti-clockwise.


**Note**

Be careful not to remove so much volume that air is introduced onto the rectangular array via the outlet.





- 21.6 Load  800  $\mu\text{L}$  of FLB (plus FLT) into the flow cell via the inlet port, dispense slowly and smoothly trying to avoid the introduction of any air bubbles.

- 21.7 Wait for  00:05:00 .

- 21.8 Gently lift the SpotON cover to open the SpotON port.

- 21.9 Load another  200  $\mu\text{L}$  of FLB (plus FLT) into the flow cell via the inlet port, this will initiate a siphon at the SpotON port to allow you to load the library dilution.


- 21.10 In a new tube prepare the library dilution for sequencing:

Component	Volume
SQB	 37.5 $\mu\text{L}$
LB	 25.5 $\mu\text{L}$
Final library	 12 $\mu\text{L}$
<b>Total</b>	 75 $\mu\text{L}$

**Note**

Mix LB immediately before use as they settle quickly.

Dilute library in EB if required.

- 21.11 Mix the prepared library gently by pipetting up and down just prior to loading.
- 21.12 Add the  75 µL library dilution to the flow cell via the SpotON sample port in a dropwise fashion. Ensure each drop siphons into the port before adding the next.
- 21.13 Gently replace the SpotON sample port cover, making sure the bung enters the SpotON port, close the inlet port and close the MinION lid.
- 22 Start the sequencing run using MinKNOW.

#### Protocol



NAME

**Starting a MinION sequencing run using MinKNOW**

CREATED BY

Josh Quick

**PREVIEW**

- 22.1 If required plug the MinION into the computer and wait for the MinION and flowcell to be detected.
- 22.2 Choose flow cell 'FLO-MIN106' from the drop-down menu.
- 22.3 Then select the flowcell so a tick appears.
- 22.4 Click the 'New Experiment' button in the bottom left of the screen.
- 22.5 On the New experiment popup screen, select the running parameters for your experiment from the individual tabs:



**Experiment:** Name the run in the experiment field, leave the sample field blank.

**Kit: Selection:** Select LSK109 as there is no option for native barcoding (NBD104).

**Run Options:** Set the run length to 6 hours (you can stop the run once sufficient data has been collected as determined using RAMPART).

**Basecalling:** Leave basecalling turned but select 'fast basecalling'.

**Output:** The number of files that MinKNOW will write to a single folder. By default this is set to 4000 but can be reduced to make RAMPART update more frequently.

Click 'Start run'.

22.6 Monitor the progress of the run using the MinKNOW interface.