

May 01, 2025

Eco-Friendly Efficient DNA Staining Using Modified Methylene Blue Solution: A Frugal and Safer Alternative to Ethidium Bromide in Low-Resource Biotechnology Laboratories

DOI

<https://dx.doi.org/10.17504/protocols.io.dm6gp27p8lzp/v1>

Joseph shenekji¹, Kamar Shayah¹, Sara Shanta¹, Ladi Serawan¹, Elaf Bobaky¹, Ahmad Mouzayek¹

¹Department of Biotechnology Engineering, Faculty of Technical Engineering, University of Aleppo

Joseph shenekji: PhD in Biotechnology Engineering - josephshenekji94@gmail.com

Kamar Shayah: PhD in Biotechnology Engineering

Sara Shanta: Graduate student in Biotechnology Engineering

Ladi Serawan: Graduate student in Biotechnology Engineering

Elaf Bobaky: Graduate student in Biotechnology Engineering

Ahmad Mouzayek: Graduate student in Biotechnology Engineering

Reclone.org (The Reagent Collaboration Network)


Tech. support email: protocols@recode.org

[Click here to message tech. support](#)



Joseph shenekji

department of biotechnology engineering, faculty of technica...



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

[Create free account](#)

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gp27p8lzp/v1>

Protocol Citation: Joseph shenekji, Kamar Shayah, Sara Shanta, Ladi Serawan, Elaf Bobaky, Ahmad Mouzayek 2025. Eco-Friendly Efficient DNA Staining Using Modified Methylene Blue Solution: A Frugal and Safer Alternative to Ethidium Bromide in Low-Resource Biotechnology Laboratories. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.dm6gp27p8lzp/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 01, 2025

Last Modified: May 01, 2025

Protocol Integer ID: 206713

Keywords: DNA, staining, methylene blue, low-cost, eco-friendly, friendly efficient dna staining, need for carcinogenic dye, safe dyes like sybr, safe dye, carcinogenic dye, using modified methylene blue solution, visualizing dna, modified methylene blue solution, safer alternative to ethidium bromide, ethidium bromide, resource biotechnology laboratories this protocol, research biotechnology laboratory, resource biotechnology laboratory, biotechnology, agarose gel electrophoresis, using methylene, ethidium, laboratory, dna

Abstract

This protocol presents a validated, low-cost, and environmentally sustainable method for visualizing DNA using methylene blue staining after agarose gel electrophoresis. It eliminates the need for carcinogenic dyes like ethidium bromide and does not require UV transilluminators, knowing that safe dyes like SYBR safe are considered expensive for low-budget labs, making this method especially suitable for educational and research biotechnology laboratories in developing regions, in Aleppo in our case.

Materials

Materials and Reagents

	A	B
	Reagent	Purpose
	Methylene blue (powder)	Intercalating dye that binds to nucleic acids
	Sodium acetate (CH ₃ COONa)	Enhances DNA-dye retention via ionic stabilization and pH control
	Glacial acetic acid	Adjusts pH for optimal dye-DNA interaction (acidic range)
	Sodium chloride (NaCl)	Increases ionic strength, promoting stable dye binding to DNA
	Ethanol (99%) or Acetic acid 5%	Fixation of stained DNA to prevent fading
	Distilled water (dH ₂ O)	Solvent and diluent
	Loading dye	Adds density and tracking dyes to DNA samples for gel loading
	Agarose	Medium for DNA electrophoresis
	TAE or TBE buffer	Electrophoresis running buffer
	DNA ladder/marker	Reference for molecular weight visualization

Equipment

- Electrophoresis apparatus with power supply
- Gel casting tray and combs
- Light box or white-light transilluminator
- pH meter or pH paper
- Balance (precision: ±0.001 g)
- Pipettes and microcentrifuge tubes
- Beakers, measuring cylinders
- Heating source (optional for increasing stain uptake)

Introduction

- 1 Ethidium bromide (EtBr), though widely used for DNA staining, poses significant health and environmental risks due to its mutagenicity and the need for UV light exposure. Methylene blue offers a non-toxic, UV-free alternative. This protocol describes an optimized approach for using methylene blue to visualize DNA ladder bands, focusing on the molecular roles of each component for improved clarity and reproducibility.

Preparation of Enhanced Methylene Blue Staining Solution

2 Prepare 0.5 M Sodium Acetate (pH 5.5):

To prepare a robust, low-toxicity staining solution that maximizes DNA interaction.

- Weigh  4.1 g sodium acetate and dissolve in ~  80 mL distilled water.
- Adjust pH to  5.5 using glacial acetic acid (kitchen vinegar could do just fine).
- Bring volume to  100 mL .

Molecular Role: Sodium acetate maintains an acidic environment that enhances the electrostatic binding between methylene blue and DNA phosphate backbones.

3 0.25 gr methylene blue Preparation of Methylene Blue Stock Solution:

- Weigh  0.25 g  0.25 g methylene blue and dissolve in  1 L  1 L distilled water.



0.25 gr methylene blue



dissolving methylene blue in 1 L of distilled water

- Add 100 mL of 0.5 Mass Percent sodium acetate (final concentration ~ 0.05 Mass Percent).
- **Optional:** Add 0.1 Mass Percent NaCl (0.58 g /100 mL) to stabilize staining.

Molecular Role: Methylene blue intercalates into the minor groove of DNA, especially under acidic conditions where the positive charge on the dye is stabilized.

4 Optional Fixative Preparation:

- 10% Ethanol: Mix 10 mL ethanol + 90 mL water.
- 5% Acetic Acid: Mix 5 mL glacial acetic acid + 95 mL water.

DNA Electrophoresis

1h 57m

5 Prepare 1% Agarose Gel:

- Dissolve 1 g agarose in 100 mL TAE buffer.
- Heat until fully melted.
- Cool slightly, then pour into casting tray with comb.

6 loading ladder markerLoad DNA Ladder:

5m

- Mix 5 μ L loading dye with DNA ladder. (1kb Gene Ruler Fermentas)
- Load samples into wells and run gel at:
- 50 V for 00:05:00, then
- 100 V for 00:30:00.



loading ladder marker


Function: Electrophoresis separates DNA based on molecular weight.

7

Staining the Gel

1h 30m

- Post-Electrophoresis Staining:
- Immerse gel in prepared methylene blue staining solution.
- Stain at room temperature for  00:20:00 , or for  01:00:00 for enhanced sensitivity in lower temperatures in the winter.



staining solution

Optional Fixation (Recommended):

Soak stained gel in  10 % (v/v) ethanol or  5 % (v/v) acetic acid for

 00:10:00 to fix bands, you could put  10 mL of ethanol or  5 mL acetic acid

to the same solution of staining containing methylene blue.

8

Destaining:

20m

- Rinse gel in dH₂O or 1X TAE for  00:15:00 .

- Change wash solution every  00:05:00 until background is sufficiently clear.
- use  30 °C water to have faster results.

9 Visualization and Interpretation

2m

- Visualize DNA bands under white light or a light box. (we used a white plexiglass and a fixed a light under it).
- Bands appear as dark blue lines against a lighter background.
- Faint bands may require extended destaining or longer staining time.
- in case you rinsed too much and the bands disappeared, you can simply put the gel back into staining solution for 5-10 minutes and rinse for  00:02:00 , and the bands will reappear.
- you could use smartphone image editing apps to enhance contrast.

10 Troubleshooting and Optimization

A	B	C
Issue	Cause	Solution
Faint bands	Insufficient staining time or low DNA amount	Increase staining time or DNA concentration
High background	Incomplete rinsing or too high dye concentration	Increase washing steps, optimize dye concentration
Bands fade quickly	No fixation step applied	Use 10% ethanol or 5% acetic acid to fix dye
Uneven staining	Poor mixing or uneven gel immersion	Ensure uniform staining solution coverage

Troubleshooting several problems

11 Comparison Table: Methylene Blue vs. Ethidium Bromide


A	B	C
Feature	Methylene Blue	Ethidium Bromide
Cost	Low	High
Safety	Non-toxic	Carcinogenic
UV Light Required	No (white light suffices)	Yes

	A	B	C
	Environmental Impact	Eco-friendly	Hazardous waste
	Staining Time	20 - 30 minutes	15–30 min
	Band Persistence	Long-lasting (if fixed)	Fades quickly under UV
	Sensitivity	Moderate	High
	Ease of Disposal	Simple	Requires hazardous disposal

comparison between EtBr and Methylene Blue

Results

12 Interpretation of Gel Electrophoresis Results Using Ethidium Bromide and Methylene Blue Compared to a Commercial DNA Ladder



comparing DNA staining using EtBr, Methylene blue, and Commercial standard.

This image presents a comparative analysis of DNA band visualization using two common staining methods-ethidium bromide (EtBr) and methylene blue-on a 1% agarose

gel in TAE buffer, alongside a commercial 1 kb DNA ladder (GeneRuler, Fermentas) as a molecular weight reference.

Band Patterns and Sensitivity

The **EtBr**-stained lane displays **14 distinct bands**, visualized as bright orange-yellow under UV illumination. Ethidium bromide intercalates between DNA bases, providing high sensitivity and clear band definition, especially for fragments of medium and lower molecular weight.

The **methylene blue**-stained lane reveals **12 visible bands**, appearing as dark blue bands on a lighter blue background. Methylene blue enables DNA visualization under normal white light, offering a safer alternative to EtBr, though with slightly reduced sensitivity, particularly for smaller DNA fragments.

The **reference DNA ladder** lane contains **15 bands**, corresponding to fragment sizes ranging from 75 to 20,000 base pairs, with key marker bands (5000, 1500, 500 bp) highlighted for orientation. Red arrows in the image indicate the correspondence between visible bands in the methylene blue lane and their respective sizes in the ladder.

Comparative Analysis

Ethidium bromide demonstrates higher sensitivity, allowing detection of more bands and lower DNA quantities compared to methylene blue. This is evident from the 14 bands visualized with EtBr versus 12 with methylene blue, while the original ladder contains 15 bands.

Methylene blue, while less sensitive, still provides adequate resolution for most DNA fragments, and its bands align accurately with the expected sizes in the commercial ladder, as indicated by the red arrows.

Both staining methods allow for reliable estimation of DNA fragment sizes when compared to the reference ladder, but EtBr is superior for detecting low-abundance or small fragments.

Practical Implications

Ethidium bromide is recommended when maximum sensitivity is required, such as for detecting low-concentration DNA samples. However, due to its mutagenic properties, appropriate safety measures must be implemented.

Methylene blue offers a safer, non-mutagenic alternative suitable for routine applications and for downstream recovery of DNA, albeit with lower sensitivity for small or faint bands.



The choice of stain should be guided by experimental requirements: EtBr for maximal sensitivity, methylene blue for user safety and DNA integrity.

13 **Conclusion**

The comparison demonstrates that while both stains are effective for DNA visualization in agarose gels, ethidium bromide provides greater sensitivity and band resolution, closely approaching the full band pattern of the original commercial ladder. Methylene blue, though less sensitive, remains a practical and safer option for many standard Biotechnology laboratory applications, with low-cost and high availability, it is the best candidate for students labs, in cases where safe DNA stains are not available or expensive.

Protocol references

- Hu, Z., & Tong, C. (2007). Synchronous fluorescence determination of DNA based on the interaction between methylene blue and DNA.. *Analytica chimica acta*, 587 2, 187-93 . <https://doi.org/10.1016/J.ACA.2007.01.050>.
- Vardevanyan PO, Antonyan AP, Parsadanyan MA, Shahinyan MA, Hambardzumyan LA. Mechanisms for binding between methylene blue and DNA. *Journal of Applied Spectroscopy*. 2013 Sep;80:595-9.
- Merdekawati F, Gustira Rahayu I, Khoirul Abror Y. Check Similarity: METHYLENE BLUE AS ALTERNATIVE DNA STAINING IN ELECTROPHORESIS. In *Proceeding of The 5th International Conference on Interprofessional Health Collaboration and Community Empowerment 2022* (Vol. 4, No. 1, pp. 87-92). PPM Poltekkes Kemenkes Bandung M Poltekkes Kemenkes Bandung.