



May 01, 2025

EdU Click-iT Administration

 The Journal of Cell Biology

DOI

<https://dx.doi.org/10.17504/protocols.io.8epv5m5bd1b/v1>

Pia Rodriguez¹

¹Duke University



Pia Rodriguez

Duke University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv5m5bd1b/v1>

External link: <https://doi.org/10.1083/jcb.202410130>

Protocol Citation: Pia Rodriguez 2025. EdU Click-iT Administration. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.8epv5m5bd1b/v1>

Manuscript citation:

Salazar MPR, Kolanukuduru S, Ramirez V, Lyu B, Manigault G, Sejourne G, Sesaki H, Yu G, Eroglu C (2025) Mitochondrial fission controls astrocyte morphogenesis and organization in the cortex. The Journal of Cell Biology 224(10). doi: [10.1083/jcb.202410130](https://doi.org/10.1083/jcb.202410130)



License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 30, 2025

Last Modified: May 01, 2025

Protocol Integer ID: 203629

Keywords: EdU, imaging kit, thermo fisher click, imaging, it administration this protocol, c10640, fluorophore, protocol

Funders Acknowledgements:

Paul and Daisy Soros Fellowship

HHMI Gilliam Fellowship

Grant ID: GT15076

Chan Zuckerberg Initiative, Neurodegeneration Challenge Network Collaborative

Grant ID: DAF2018-191999

Chan Zuckerberg Initiative, Neurodegeneration Challenge Network Collaborative

Grant ID: DAF2021-237435

BRAIN Initiative

Grant ID: U19-NS123719

Aligning Science Across Parkinson's initiative

Grant ID: ASAP-020607

Howard Hughes Medical Institute

Abstract

This protocol uses the Thermo Fisher Click-iT EdU Imaging Kit (I used the 647 fluorophore, C10640, <https://www.thermofisher.com/order/catalog/product/C10640>).

Materials

Thermo Fisher Click-iT EdU Imaging Kit

Prepare EdU Stock Solution

- 1 Add 2 ml of DPBS to 5 mg EdU stock vial (component A).
- 2 Aliquot into 100 μ l aliquots and store at -80°C .
- 3 This makes a 10 μM solution of EdU (2.5 $\mu\text{g}/\mu\text{l}$).

EdU Injection into Pups - IP Injection

- 4 Must inject 10-15 μg of EdU/gram of animal weight, so 4-6 μl of working 10 μM stock per gram of animal weight.
- 5 Thaw out appropriate number of aliquots depending on weight and number of pups/mice you wish to inject. Combine into one vial and take to the mouse house. For CD1/bl6 pups, P3 pups weigh ~ 3 g for example.
- 6 Pull all volume into syringe to inject pups as you weigh them. I inject 6 $\mu\text{l}/\text{g}$ with a 2 ml syringe (this syringe uses ~ 70 μl of void volume, so add that to your total volume).
- 7 Weigh one pup at a time and inject IP appropriate volume (I round up for injection, i.e., 3 g pup would get 18 μl of 10 μM stock, so 20 μl injection).

Tissue Collection and Staining

- 8 Collect tissue as normal, section, block with NGS.
- 9 After blocking, follow Click-iT EdU Thermo Fisher staining protocol (Starting with "Detect EdU section, step 4.1).
- 10 EdU staining takes ~ 1 hr.



- 11 After EdU staining, can continue with regular primary antibody incubation/staining protocol as needed.