

Jan 11, 2019

EIT data acquisition in rat sciatic nerve using stimulation of tibial and peroneal branches

In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.ww7ffhn

Kirill Aristovich¹, Matteo Donega², David Holder¹

¹University College London, University of London; ²Galvani Bioelectronics


SPARC

Tech. support email: info@neuinfo.org



Kirill Aristovich

University College London



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ww7ffhn>

External link: <http://iopscience.iop.org/article/10.1088/1741-2552/aad78e/pdf>



Protocol Citation: Kirill Aristovich, Matteo Donega, David Holder 2019. EIT data acquisition in rat sciatic nerve using stimulation of tibial and peroneal branches. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.ww7ffhn>

Manuscript citation:

Kirill Aristovich et al 2018 J. Neural Eng. 15 056025

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: January 09, 2019

Last Modified: January 11, 2019

Protocol Integer ID: 19135

Keywords: imaging, neural traffic, EIT, electroseuticals, cuff electrode, eit imaging dataset, eit data acquisition in rat, sciatic nerve, eit data acquisition, imaging

Abstract


High-level protocol to obtain the EIT imaging dataset and reconstruct images for the paper.

Guidelines

Please consult UCL EIT team before attempting the protocol, the complexity of the required hardware and software is such, that the collaboration will be much more effective. It is possible to create customised easy-to-operate systems and software packages which are tailored for the specific paradigm

Troubleshooting

Safety warnings

 Chemicals used during the electrode interface construction are toxic.

Before start

consult UCL EIT team before attempting the protocol

Hardware

1

Prepare nessesary hardware:

Equipment

UCL EIT system

NAME

EIT system

TYPE

Custom

BRAND

n/a

SKU

<https://github.com/EIT-team/ScouseTom>

LINK

Number of channels: up to 256

SPECIFICATIONS

Current frequency range: DC - 20kHz

Sampling frequency: 100kHz

Onboard stimulator



Neural Interface

2


Prepare neural interface:

Follow the instructions in teh paper <http://iopscience.iop.org/article/10.1088/1741-2552/aae868>

Safety information

Some of the chemical components are toxic
Chemical hood cabinet is required

Expected result



Software

3

Software

UCL EIT system

NAME

UCL EIT group

DEVELOPER

<https://github.com/EIT-team/ScouseTom/tree/master/src>

SOURCE LINK

Animal Preparation

- 4 Procedures were performed on a vibration isolated table (Thorlabs Inc., USA) and, throughout experiments, the core body temperature of the rat was controlled with a homeothermic heating unit (Harvard Apparatus, Kent, UK) and maintained at 37°C

Animal Preparation

5 Sprague-Dawley adult male rats weighing 300 to 450g.

6 Anaesthesia is induced with a 5% isoflurane in 100% O₂.

7 ntubated using a small animal laryngoscope

Equipment

Small animal laryngoscope

NAME

Laryngoscope

TYPE

Custom 3D printed

BRAND

n/a

SKU

<https://www.thingiverse.com/thing:148315>

LINK



8 18G cannula and mechanical ventilation is provided, using a Harvard Apparatus Inspira Ventilator (Harvard Apparatus, Ltd, UK), with a 50/50% gas mixture of O₂ and air.

9 Arterial and venous access was established through cannulation (BD Insyte/Vialon, Becton, Dickinson U.K. Ltd.) of the right femoral vessels.

10 The arterial blood pressure should is monitored (Cardiicap 5, Datex Ohmeda) and the mean arterial pressure (MAP) kept between 90 and 110mmHg using labetalol and adrenaline as necessary.

11 Once intravenous access had been established, a constant infusion of propofol is initiated. Monitor the depth with pinch reflex.

Surgery

- 12 Access sciatic trunc, common peroneal nerve and tibial nerve following the guide:
<https://onlinelibrary.wiley.com/doi/full/10.1002/mus.21652>


Note

Access to the common peroneal nerve was established through a 2mm lateral incision in biceps femoris near the knee joint.

Electrode placement

- 13 Implant the stimulation electrodes on common peroneal and tibial branches (bipolar Ag hook electrodes)
- 14 Implant the EIT array on main sciatic trunc

Expected result




Experimental design and setup. (a) The EIT array was placed on the main sciatic nerve running in the posterior compartment of the thigh approximately 20 mm distal to the greater sciatic foramen. (b-c) Two stimulating electrode pairs were placed on the tibial and common peroneal nerves. EIT images were recorded during the repeated activity in the sciatic nerve evoked with 5Hz stimulation of each branch at a time.



Data aquisition

- 15 For all studies, a square biphasic (positive first) constant current temporal waveform was delivered using a balanced current source (Keithley, UK model no 6221) with:
1-3 mA amplitude,
50 μ s pulse width,
5Hz frequency (200 ms inter-stimulus time)
- 16 Program the EIT system with the following paerameters:

 EIT_parameters.mat
- 17 Connect stimulator to the tibial stimulation electrode. Collect the data
- 18 Connect stimulator to the common peroneal stimulation electrode. Collect the data

Data Processing and Image reconstruction

- 19 Load and process the data using the following software:

Software

UCL LoadData

NAME

UCL EIT team

DEVELOPER

https://github.com/EIT-team/Load_data

SOURCE LINK

- 20 Forward solution



Software

PEITS

NAME

UCL EIT team

DEVELOPER

<https://github.com/EIT-team/PEITS>

SOURCE LINK

21 Inverse solution and images generation

Software

Image Reconstruction

NAME

UCL EIT team

DEVELOPER

<https://github.com/EIT-team/Reconstruction>

SOURCE LINK