



Jan 25, 2025

Electroporation & Electrocompetent Cell Preparation Protocol for DAP auxotroph E.coli WM3064



Forked from [Electroporation Protocol](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.e6nvw1rxwlmk/v1>

Atharva Karde^{1,2}, Hans Christopher Bernstein^{1,2,3}

¹The Norwegian College of Fishery Science, Faculty of Biosciences, Fisheries and Economics, UiT - The Arctic University of Norway, Tromsø, Norway;

²Microalgae & Microbiomes Research Group (M2RG), UiT - The Arctic University of Norway, Tromsø, Norway;

³The Arctic Centre for Sustainable Energy (ARC), UiT - The Arctic University of Norway, Tromsø, Norway



Atharva Karde

UiT The Arctic University of Norway

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvw1rxwlmk/v1>

Protocol Citation: Atharva Karde, Hans Christopher Bernstein 2025. Electroporation & Electrocompetent Cell Preparation Protocol for DAP auxotroph E.coli WM3064. [protocols.io](#) <https://dx.doi.org/10.17504/protocols.io.e6nvw1rxwlmk/v1>



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We regularly use this protocol.

Created: October 02, 2024

Last Modified: January 25, 2025

Protocol Integer ID: 108862

Keywords: Electroporation, Competent cells, DAP auxotroph, WM3064, Transformation, Electrocompetent, Ecoli, E.coli, electrocompetent cell preparation protocol, dna through electroporation, electrocompetent cell, aliquots of electrocompetent cell, prep of electroporation, electroporation, cell, dna, protocol, preparation,

Funders Acknowledgements:

Kometen – Nord-Norges forskningsfond

Abstract

This is a protocol for making electrocompetent E.coli WM3064 cells and transforming them with DNA through electroporation. This protocol produces  30 µL aliquots of electrocompetent cells, each enough for one prep of electroporation.

Guidelines

Electroporation Protocol

The electroporation protocol will vary depending on the strain so this protocol may need to be optimized.

Transformation efficiency calculation

If the culture was diluted 1000-fold when plated, the total cfu per ml is 1000 times the number of colonies counted. The cfu is divided by the amount of plasmid (e.g. 10 pg per ml)

Transformation efficiency = cfu/ µg = (colonies counted*1000) / (0.00001 µg plasmid)

Materials

MATERIALS

☒ Magnesium sulfate heptahydrate Merck MilliporeSigma (Sigma-Aldrich) Catalog #M2773

☒ NaCl Merck MilliporeSigma (Sigma-Aldrich) Catalog #53014

☒ Tryptone Fisher Scientific Catalog #BP1421-500

☒ Glucose Merck MilliporeSigma (Sigma-Aldrich) Catalog #G8270

☒ Magnesium chloride hexahydrate Merck MilliporeSigma (Sigma-Aldrich) Catalog #M2670

☒ Potassium chloride Merck MilliporeSigma (Sigma-Aldrich) Catalog #P9333

☒ Glycerol Thermo Fisher Scientific Catalog #17904

☒ Yeast Extract Thermo Fisher Catalog #211930

☒ 2,6-Diaminopimelic acid Thermo Scientific

Media

SOB + DAP:

2% tryptone

0.5% yeast extract

10 mM NaCl

2.5 mM KCl

10 mM MgCl₂

10 mM MgSO₄

0.3 mM DAP

SOC + DAP :

SOB + 20 mM glucose + 0.3 mM DAP

LB agar + DAP + antibiotic

LB agar + 0.3 mM DAP + appropriate antibiotic

Appropriate Antibiotics for Your Application

Antibiotics for Plasmid selection

	A	B
	Antibiotic	Working Concentration



A	B
Ampicillin	100 μg/ ml
Carbenicillin	100 μg/ ml
Chloramphenicol	33 μg/ ml
Kanamycin	50 μg/ ml
Streptomycin	25 μg/ ml
Tetracycline	15 μg/ ml

Protocol materials

⊗ 2,6-Diaminopimelic acid **Thermo Scientific**

⊗ Glycerol **Thermo Fisher Scientific Catalog #17904**

⊗ Potassium chloride **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P9333**

⊗ Yeast Extract **Thermo Fisher Catalog #211930**

⊗ Magnesium sulfate heptahydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M2773**

⊗ Tryptone **Fisher Scientific Catalog #BP1421-500**

⊗ Glucose **Merck MilliporeSigma (Sigma-Aldrich) Catalog #G8270**

⊗ Magnesium chloride hexahydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M2670**

⊗ NaCl **Merck MilliporeSigma (Sigma-Aldrich) Catalog #53014**



Safety warnings

⚠ Please refer to the Safety Data Sheets (SDS) for health and environmental hazards.

Before start

Quite a bit of glycerol gets left behind in the tube during the washing steps while discarding the supernatant. Note that this dilutes the cells a bit. Prepare aliquots based on the cell quantity required specifically for your use.

Dilute the plasmid to have 1 - 3 μL of DNA (<100 ng) per electroporation prep. Calculate the amount of DNA used per prep to determine the transformation efficiency later.

The electroporation mixture ratio will vary depending on the strain and DNA size so this protocol may need to be optimized.

Make sure to include a negative control which also serves as a "survival control". This mixture should contain only electrocompetent cells and no plasmid. Plate this mixture on LB+DAP agar to assess whether the cells survive the electroporation pulse.

This protocol was developed as part of a master's thesis affiliated with the Microalgae & Microbiomes Research Group at NFH, UiT, Norway.

M2RG Website : <https://uit.no/research/micro>



Growing overnight culture

- 1 Add a loopful of bacterial culture to  5 mL LB broth.
- 2 Incubate  Overnight at  37 °C with shaking.



Making electrocompetent cells

- 3 Add  10 mL SOB medium containing  0.3 millimolar (mM)  2,6-Diaminopimelic acid **Thermo Scientific** in a 50 mL conical flask.
- 4 Inoculate each flask with overnight culture to OD600 = 0.05
- 5 Grow cells until mid-exponential phase (OD600 = 0.6) or  02:00:00 
- 6 Transfer the cells to a 15 mL Falcon tube
- 7 Chill the culture  On ice for  00:10:00 .
- 8 Centrifuge the cells at  6000 rpm, 4°C, 00:05:00 . 
- 9 Discard the supernatant.
- 10 Wash the pellet by adding  10 mL chilled  10 % volume  Glycerol **Thermo Fisher Scientific Catalog #17904** and vortex vigorously.
- 11 Centrifuge  6000 rpm, 4°C, 00:03:30 and discard the supernatant. 



12 go to step #10 Repeat for a total of 4 wash cycles .

13 Resuspend the washed pellet in 100 μ L 10 % volume

Glycerol Thermo Fisher Scientific Catalog #17904

14 Make 30 μ L aliquots in sterile microcentrifuge tubes. Proceed to electroporation or store at -80 °C . One tube allows for one prep of electroporation.

Thawing and preparation

15 Turn on electroporator and set to **1.5 kv, 200 ohms** and **25 μ F**.

16 Pre-warm recovery SOC+DAP medium and LB-antibiotic plates at 37 °C .

17 Thaw microcentrifuge tubes with electrocompetent cells On ice or use freshly made cells.

18 Place appropriate number of 1 mm-electroporation cuvettes On ice .

Electroporation and recovery

19 Add 1 μ L DNA solution (<100 ng) to 30 μ L electrocompetent cells in microcentrifuge tubes.

Note

This mixture can be optimized for your strain and plasmid. Trying a few combination ratios is recommended in order to determine which one yields the highest transformation efficiency. A lower amount of DNA usually yields higher efficiency.

20 Transfer the DNA-cell mixture to the cold cuvette, tap on countertop 2X, wipe water from exterior of cuvette and place in the electroporation module and press pulse (**don't hold**



the button down).

- 21 Immediately add  970 μL 37°C SOC + DAP , mix by pipetting up and down once and transfer to a microcentrifuge tube. 

Note

Try to add the recovery medium as soon as possible. Transformation efficiency sharply declines with increase in time between pulse and addition of recovery medium.

- 22 Incubate at  37 °C for  01:00:00 . 

Dilution plating and efficiency calculation

- 23 Make appropriate dilutions and plate  100 μL on LB agar + DAP + antibiotic plates. 

- 24 Incubate  Overnight at  37 °C . 

Note

Transformation efficiency calculation:

If the culture was diluted 1000-fold when plated, the total cfu per ml is 1000 times the number of colonies counted. The cfu is divided by the amount of plasmid (e.g. 10 pg per ml)

Transformation efficiency = cfu/ μg = (colonies counted*1000) / (0.00001 μg plasmid)

Protocol references

1. Short Protocols in Molecular Biology, Chapter 1.
2. Datsenko, K.A., and Wanner, B.L. 2000. One-step inactivation of chromosomal genes in Escherichia coli K-12 using PCR products. Proc. Natl. Acad. Sci. U. S. A. 97: 6640-6645.
3. <https://barricklab.org/twiki/bin/view/Lab/ProtocolsElectrocompetentCells>



Acknowledgements

This protocol was developed as part of a master's thesis affiliated with the Microalgae & Microbiomes Research Group at NFH, UiT, Norway.

M2RG Website : <https://uit.no/research/micro>

The E.coli strain used in this project was generously provided by the Lowe-Power Lab at UC Davis (<https://lowepowerlab.ucdavis.edu/>).