

Oct 05, 2025

Electroporation

DOI

<https://dx.doi.org/10.17504/protocols.io.14egnrneyl5d/v1>

igemnthuth¹

¹iGEM_NTHU_Taiwan

iGEM NTHU



igemnthuth

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egnrneyl5d/v1>

Protocol Citation: igemnthuth 2025. Electroporation. protocols.io <https://dx.doi.org/10.17504/protocols.io.14egnrneyl5d/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: October 03, 2025

Last Modified: October 05, 2025

Protocol Integer ID: 228975

Keywords: plasmid into dh10b, designed plasmid, plasmid, electroporation, dh10b, supporting orthogonal replication, orthogonal replication

Abstract

Introduce the designed plasmid into DH10B to establish *E. coli* supporting orthogonal replication.

Materials

- Plasmid
- Pre-prepared DH10B competent cells
- SOB medium pre-warmed to 37 °C
- Sterile electroporation cuvette (2 mm gap)
- Antibiotic agar plates (according to the resistance gene carried by the plasmid)
- Electroporator
- Shaking incubator at 37 °C (220 rpm)
- Biosafety cabinet
- Ice bucket with ice



Preparation of the Mixture

- 1 Take out pre-prepared DH10B competent cells and thaw on ice (~10 minutes).
- 2 In a sterile workspace, add 1.0 μL of plasmid DNA (dissolved in ultrapure water, ~50–100 ng/ μL) to the thawed competent cells. Do not stir, pipette harshly, or create bubbles.
- 3 Carefully transfer the mixture into a pre-chilled 2 mm electroporation cuvette (kept on ice ≥ 10 minutes). Close the lid and leave on ice for 2 minutes.

Electroporation

- 4 Place the cuvette into the electroporator and set the following conditions:
 - **Voltage:** 2500 V (standard for 2 mm gap)
 - **Time constant:** 4.5–5.5 ms (ideal range)
 - **Pulse number:** 1 pulse, duration should not exceed 6 ms
- 5 A faint spark sound should be heard during electroporation. If the instrument displays “Arcing” (discharge failure), replace the cuvette immediately and check whether the sample contains salt.

Immediate Recovery

- 6 Immediately add 1 mL of pre-warmed (37 °C) LB medium into the cuvette and gently pipette to mix.
- 7 Transfer the cuvette contents into a sterile 15 mL culture tube.
- 8 Incubate at 37 °C in a shaking incubator at 220 rpm for 60 minutes (recovery culture).

Antibiotic Selection (see detailed steps in next protocol)

- 9 After 1 hour of recovery, plate 100 μL of the culture onto an LB agar plate containing the appropriate antibiotic.
- 10 Incubate the plate at 37 °C in an inverted position (agar side up) overnight (12–16 hours).

