

Oct 25, 2019

## Electroporation of *Thalassiosira pseudonana*

DOI

[dx.doi.org/10.17504/protocols.io.yvqfw5w](https://dx.doi.org/10.17504/protocols.io.yvqfw5w)

Joshua Bugge<sup>1</sup>, Deborah Robertson<sup>1</sup>

<sup>1</sup>Clark University

Protist Research to Opti...

Robertson Lab



Deborah Robertson

Clark University

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.yvqfw5w](https://dx.doi.org/10.17504/protocols.io.yvqfw5w)

**Protocol Citation:** Joshua Bugge, Deborah Robertson 2019. Electroporation of *Thalassiosira pseudonana*. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.yvqfw5w>

**Manuscript citation:**

Bugge, J. 2015. Electroporation-mediated Transformation and post-transcriptional gene regulation of nitrate reductase in the marine diatom *Thalassiosira pseudonana*. M.S., Clark University, 89 pp.

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** In development

**We are still developing and optimizing this protocol**

**Created:** March 06, 2019

**Last Modified:** October 25, 2019

**Protocol Integer ID:** 21136

**Keywords:** electroporation, diatoms, *Thalassiosira*



## Abstract

An electroporation-mediated genetic transformation of the marine diatom *Thalassiosira pseudonana* was developed. Using a sorbitol-based buffer, *T. pseudonana* cells were successfully transformed with the *Tpfc*p/nat plasmid at an efficiency of 2820 per  $10^8$  cells. This represents a six-fold improvement compared with previously published methods.

## Materials

### MATERIALS

 f/2 medium **NCMA Catalog #MKF250L**

375 mM sorbitol (sterile filtered)

2 mM electroporation cuvettes (chilled on ice)

50 mL conical centrifugation tubes

1.5 mL centrifuge tubes

Selection agent

NEPC plates

- 1 Grow cultures of *T. pseudonana* on f/2 supplemented sterile seawater to a density of approximately  $1.1 \times 10^6$  cells mL<sup>-1</sup>

## Cell collection option 1

- 2 All procedures are carried out at 4 C and cells were kept on ice. Collect cells from 500 mL of culture (see step 1) by centrifugation for 10 min at 3000 x g in 10 x 50 mL conical centrifuge tubes. Discard supernatant and resuspend pellet in 1 mL of 375 mM sorbitol (filter sterilized)
- 3 Combine resuspended cells into one 50 mL conical centrifuge tube and centrifuge at 3000 x g for 10 min. Discard supernatant and resuspend pellet in 800 µL of 375 mM sorbitol (filter sterilized). Transfer to 1.5 mL tube and store on ice.

## Cell collection option 2

- 4 All procedures are carried out at 4 C and cells were kept on ice. Concentrate cells by vacuum filtration to approximately 300 mL. Collect cells by centrifugation in 50 mL conical centrifugation tubes at 3000 x g for 10 min. Discard supernatant and resuspend pellet in 2 mL of 375 mM sorbitol (filter sterilized).
- 5 Combine cells into new tube and centrifuge at 3000 x g for 10 min. Discard supernatant and resuspend in 1 mL of 375 mM sorbitol.

## Electroporation

- 6 On ice, add linearized plasmid (in 250 mM sorbitol) to a final concentration of 0.15 µg mL<sup>-1</sup> to the concentrated cells from above.
- 7 Transfer 100 µL of cells + plasmid to a pre-chilled 2 mm gap electroporation cuvette.
- 8 Wipe cuvette dry with kimwipe and place into Biorad Gene Pulser electroporator set to 25 µF and 400 ohm resistance. In our hands, the highest transformation efficiencies were obtained using a field strength of 2.50 kV/cm.
- 9 Transfer electroporated cells to 10 mL of sterile, f/2 supplemented artificial seawater and incubate overnight at 17 C under constant light (170 µE m<sup>-2</sup> s<sup>-1</sup>).



- 10 The following day, plate ca.  $5 \times 10^8$  cells on selective NEPC agar plates (in our experiment we used  $100 \mu\text{g mL}^{-1}$  NAT). Controls should include no-plasmid and no electroporation treatments.
- 11 Grow plates under constant light ( $170 \mu\text{E m}^{-2} \text{s}^{-1}$ ). Colonies appear approximately 10 days after plating.
- 12 Transfer colonies to 300  $\mu\text{L}$  of f/2 supplemented artificial seawater in 96 well plates. Optical density (600 or 680 nm) can be used as a proxy for growth using a plate reader.