



Apr 27, 2018

Electroporation of *Vibrio natriegens* (Weinstock *et al.* 2016, modified .)

DOI

dx.doi.org/10.17504/protocols.io.ptmdnk6

Benjamin Daniel¹

¹iGEM Marburg



Benjamin BBJD Daniel

iGEM Marburg

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ptmdnk6>

Protocol Citation: Benjamin Daniel 2018. Electroporation of *Vibrio natriegens* (Weinstock et al. 2016, modified .). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.ptmdnk6>

Manuscript citation:

***Vibrio natriegens* as a fast-growing host for molecular biology Weinstock M, Hesek E, Wilson C, Gibson D. Nature Methods 2016 vol: 13 (10) pp: 849-851**

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: April 27, 2018

Last Modified: April 27, 2018

Protocol Integer ID: 11853

Keywords: Electroporation , electroporation of vibrio natriegen, competent cells for vibrio natriegen, chilled electroporation cuvette electroporation, vibrio natrigen, v2 salts overnight culture inoculation, electroporation cuvette electroporation, vibrio natriegen, preparation of electrocompetent cell, electroporation, ml electroporation, electrocompetent cell, freeze in liquid nitrogen store, aliquot with electrocompetent cell, liquid nitrogen store, preparation, centrifuge, brain haert infusion, v2 salts plate, haert infusion



Abstract

A protocol outlining the Preparation and transformation of electro-competent cells for *Vibrio natriegens*

Preparation of electrocompetent cells (*Vibrio natriegens*)

1. 10 mL BHI (Brain heart infusion) + v2 salts Overnight culture
2. Inoculation of a new BHI + v2 media with 1% of the overnight culture as inoculum.
3. Grow at 37°C shaking to an OD of 0.5
4. Chill the culture on ice for 15 min
5. Use a chilled (4°C) centrifuge at 4,500 r.p.m. for 20 min
6. Decant the supernatant diligently and carefully
7. Resuspend gently in 5-10 mL Electroporation-Buffer (680 mM sucrose, 7 mM K₂HPO₄, pH 7), then fill the falcon tube to the top.
8. To wash, spin again at 4.500 r.p.m. for 15 min at 4 °C, then resuspend gently in 5 mL Electroporation-Buffer. Repeat washing twice for a total of three times.
9. Spin again at 4.500 r.p.m. for 15 min at 4 °C, decant supernatant and resuspend carefully in the residual buffer.
10. Adjust volume for a OD of 16
11. Aliquot in chilled Micro-Reaction-Tubes, then snap freeze in liquid nitrogen
12. Store at -80°C


Electroporation


1. remove one aliquot with electrocompetent cells from storage
2. keep on ice until thawed
3. add plasmid and mix gently
4. Transfer to chilled electroporation cuvette
5. Electroporation is done with the following parameters: **700**(-900) V, 25 µF, 200 Ω in a 1mm cuvette
6. Recover for one hour at 37°C in brain heart infusion with v2 salts
7. Plate on LB2 agar plates.



Materials

MATERIALS

 sucrose

 K₂HPO₄

 v2 salt

 Brain haert infusion

Troubleshooting



1 **Preparation of electrocompetent cells**

- 2 10 mL BHI (Brain haert infusion) + v2 salts Overnight culture
- 3 Inoculation of a new BHI + v2 media with 1% of the overnight culture as inoculum.
- 4 Grow at 37°C shaking to an OD of 0.5
- 5 Chill the culture on ice for 15 min
- 6 Use a chilled (4°C) centrifuge at 4,500 r.p.m. for 20 min
- 7 Decant the supernatant diligently and carefully
- 8 Resuspend gently in 5-10 mL Electroporation-Buffer (680 mM sucrose, 7 mM K₂HPO₄, pH 7), then fill the falcon tube to the top.
- 9 To wash, spin again at 4.500 r.p.m. for 15 min at 4 °C, then resuspend gently in 5 mL Electroporation-Buffer. Repeat washing twice for a total of three times.
- 10 Spin again at 4.500 r.p.m. for 15 min at 4 °C, decant supernatant and resuspend carefully in the residual buffer.
- 11 Adjust volume for a OD of 16
- 12 Aliquot in chilled Micro-Reaction-Tubes, then snap freeze in liquid nitrogen
- 13 Store at -80°C



Electroporation

- 14 remove one aliquot with electrocompetent cells from storage
- 15 keep on ice until thawed
- 16 add plasmid and mix gently
- 17 Transfer to chilled electroporation cuvette
- 18 Electroporation is done with the following parameters: **700**(–900) V, 25 μ F, 200 Ω in a 1mm cuvette
- 19 Recover for one hour at 37°C in brain heart infusion with v2 salts
- 20 Plate on LB2 agar plates.