

Nov 10, 2023

Elution of nasal lining fluid collected via nasosorption

DOI

<https://dx.doi.org/10.17504/protocols.io.14egn3nn6l5d/v1>

Samuel Montgomery^{1,2}

¹Telethon Kids Institute; ²University of Western Australia



Samuel Montgomery

Telethon Kids Institute

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn3nn6l5d/v1>

Protocol Citation: Samuel Montgomery 2023. Elution of nasal lining fluid collected via nasosorption. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.14egn3nn6l5d/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 09, 2023



Last Modified: November 10, 2023

Protocol Integer ID: 90751

Keywords: nasosorption, nasal lining fluid, elution of nasal lining fluid, eluting nasal lining fluid, nasosorption this procedure, nasosorption device, somascan assay, mucosal diagnostics fx, dna extraction, dna extraction for 16s rrna, proteomic, metagenomic, extraction

Abstract

This procedure outlines the materials and processes for eluting nasal lining fluid collected using a Mucosal Diagnostics FX-i nasosorption device. This protocol has been successfully used to elute nasal lining fluid for both DNA extraction for 16s rRNA sequencing & metagenomics, and for proteomics using the SomaScan assay.

Guidelines

This protocol should be performed in a sterile environment to reduce cross contamination between nasosorption devices, to enable metagenomic DNA extraction and analyses

Materials

- Phosphate Buffered Saline (PBS) (#10010023)
- Tween 20 (Sigma, #P9416-50ML)
- Corning SpinX columns (Sigma, #CLS9301)
- 2mL Corning microfuge tubes (Sigma, #CLS3213)
- Sterile forceps (two pairs per swab)
- 1.6mL sterile microfuge tubes



Making elution buffer

- 1 Add **[M] 0.05 % volume** of Tween20 to PBS in a sterile environment to reduce cross-contamination
- 2 Store at **4 °C** until required

Elution of nasal lining fluid

30m 30s

- 3 Thaw nasosorption device in cryotube on ice
- 4 Add **300 µL** of elution buffer to a labelled 2mL microcentrifuge tube and place on ice
- 5 Remove the synthetic absorptive matrix (SAM) from the nasosorption device using sterile forceps by tearing it from the handle and place in the microcentrifuge tube containing elution buffer
- 6 Vortex the microcentrifuge tube for **00:00:30** 30s
- 7 Using a new pair of sterile forceps, remove the SAM from the elution buffer, and place into a spinX column on a new 2mL microcentrifuge tube
- 8 Add the elution buffer from the first tube onto the SAM in the spinX column
- 9 Centrifuge the SAM and spinX column at **16000 x g** for **00:30:00** at **4 °C** 30m
- 10 Discard the spinX column and SAM



- 11 Aliquot the eluate into labelled microtubes, and measure total protein concentration using the assay of your choice (Pierce BCA assay, Qubit protein assay etc)

Protocol references

<https://www.jove.com/t/56413/absorption-nasal-bronchial-fluids-precision-sampling-human>