Jun 03, 2019

Version 1

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Heather Etchevers

Aix-Marseille Université, INSERM

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

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**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol in our group and it is working / Eh oui, ça marche !**

**Created:** May 27, 2019

**Last Modified:** June 03, 2019

**Protocol Integer ID:** 23708

**Keywords:** cryostat, coupes, sectioning, histology, histologie, freezing tissue in oct, freezing tissue, cooled isopentane, using isopentane, tissus ou embryons dan, des coupes cryostat, downstream cryo, oct, embryons dan, ice, tissue

## Abstract

# Introduction

This protocol describes how to freeze tissues in OCT using dry-ice cooled isopentane (for downstream cryo-sectioning).

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|                                                                                                         |
|---------------------------------------------------------------------------------------------------------|
| <a href="#">Method</a>                                                                                  |
| Prior to commencing work decontaminate the work area.                                                   |
| Work with tissues from other species (e.g. mouse) should be performed in a fume hood or downflow table. |

## Safety warnings

### Chemical risks

- This protocol should be performed under a fume hood or under a ducted MSC.
- Always wear correct PPE (which includes eye protection, nitrile gloves, thermal gloves for handling Liquid Nitrogen / dry-ice cooled iso-pentane and appropriate labcoat) when handling any chemical.
- For more chemical information see the **COSHH** forms or **MSDS** for each chemical.

**Any chemicals which have specific risks and handling instructions will be outlined in the appropriate SOP method section**

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## Before start

Some fixed tissues freeze better having been impregnated in 15%, then 30% sucrose in neutral buffered PBS. Particularly for embryos or brain tissues. Let equilibrate by falling to bottom of well/tube.

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- 1 Cool 1 pair of forceps in the dry ice. Keep the other pair at room temperature. If used, place labelled cryotube on dry ice.

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- 2 On a downflow or fume hood: Place metal container on dry ice and half fill with isopentane.(use enough iso-pentane to completely immerse the tissue pieces).Iso-pentane is extremely flammable and may be fatal if swallowed or enters airways. To be used in a fume-hood or downflow table wearing correct PPE.

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- 3 Put lumps of dry ice into the iso-pentane and monitor with the -100°C thermometer until it measures -70°C.

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- 4 Label a cryomold and foil with the sample id

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
- 5 Trim the tissue into a piece small enough to fit into the cryomold, if required.

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- 6 Half fill the cryomold with OCT and position the tissue in the cryomold. Note: for spatial transcriptomics use small cryomolds (maximum 10 ×10 mm). One can prepare multiple samples in this way on wet (not dry) ice, and reposition just before freezing.

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Ensure that the tissue is oriented correctly to achieve the desired cryosection. *E.g.* Place tissue with region of interest to section, pushed flat at the bottom of the mold. Ensure there is a thin layer of OCT at the bottom of the mold so the tissue is fully covered. Dotted line in the image above shows how cross-sections would be cut from the tissue. Sections would be generated from the bottom of the block upwards.

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- 7 Add more OCT if needed to ensure the tissue is fully submerged and covered in OCT.

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- 8 Using forceps, submerge the cryomold in the container of isopentane, chilled to -70 °C on dry ice. Only samples without OCT should be submerged entirely, otherwise float the cryomold to immerse 5 of the 6 sides. It is possible to add powdered dry ice to the top of the cryomold.

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- 9 When completely frozen, wrap sample in foil (still in cryomold), and keep on dry ice. Alternatively, remove sample from cryomold and place in labelled cryotube on dry ice.

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- 10 Transfer samples on dry ice, to labelled box in  -80 °C freezer.

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- 11 If freezing more samples monitor the isopentane temperature and adjust as necessary.

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- 12 Wearing cold-resistant gloves or the wooden forceps, carefully remove the metal container of iso-pentane from the dry ice. The iso-pentane can allowed to warm to -10C, then decanted to a suitable bottle for re-use. The bottle must be labelled and kept in a fume hood/ downflow table, with the lid loosened until the iso-pentane returns to room temperature.

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