Jun 03, 2019 Version 2

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License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol in our group and it is working / Eh oui, ça marche toujours !

Created: June 03, 2019

Last Modified: June 03, 2019

Protocol Integer ID: 23978

Keywords: cryostat, coupes, sectioning, histology, histologie

Abstract

Introduction

This protocol describes how to freeze tissues in OCT using dry-ice cooled isopentane (for downstream cryo-sectioning). In this new version we add that it is also possible to use liquid nitrogen to cool the isopentane, but it's more complicated.

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Guidelines

Safety information

Liquid nitrogen is very dangerous! Appropriate handling training is necessary.

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Method
Prior to commencing work decontaminate the work area.
Work with tissues from other species (e.g. mouse) should be performed in a fume hood or downflow table.

Safety warnings

Chemical risks

- This protocol should be performed under a fume hood or under a ducted MSC.
- Always wear correct PPE (which includes eye protection, nitrile gloves, thermal gloves for handling Liquid Nitrogen / dry-ice cooled iso-pentane and appropriate labcoat) when handling any chemical.
- For more chemical information see the **COSHH** forms or **MSDS** for each chemical.

Any chemicals which have specific risks and handling instructions will be outlined in the appropriate SOP method section

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Before start

Some fixed tissues freeze better having been impregnated in 15%, then 30% sucrose in neutral buffered PBS. Particularly for embryos or brain tissues. Let equilibrate by falling to bottom of well/tube.

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- 1 Cool 1 pair of forceps in the dry ice. Keep the other pair at room temperature. If used, place labelled cryotube on dry ice. This step and dry ice remains useful, even when cooling the isopentane in liquid nitrogen.

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- 2 On a downflow or fume hood: Place metal container on dry ice and half fill with isopentane (use enough iso-pentane to completely immerse the tissue pieces).

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Safety information

Iso-pentane is extremely flammable and may be fatal if swallowed or enters airways. To be used in a fume-hood or downflow table wearing correct PPE.

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- 3 Put lumps of dry ice into the iso-pentane and monitor with the -100°C thermometer until it measures -70°C. OR : immerge

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- 4 Label a cryomold and foil with the sample id

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
- 5 Trim the tissue into a piece small enough to fit into the cryomold, if required.

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- 6 Half fill the cryomold with OCT and position the tissue in the cryomold. Note: for spatial transcriptomics use small cryomolds (maximum 10 ×10 mm). One can prepare multiple samples in this way on wet (not dry) ice, and reposition just before freezing.

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Ensure that the tissue is oriented correctly to achieve the desired cryosection. *E.g.* Place tissue with region of interest to section, pushed flat at the bottom of the mold. Ensure there is a thin layer of OCT at the bottom of the mold so the tissue is fully covered. Dotted line in the image above shows how cross-sections would be cut from the tissue. Sections would be generated from the bottom of the block upwards.

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- 7 Add more OCT if needed to ensure the tissue is fully submerged and covered in OCT.

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- 8 Using forceps, submerge the cryomold in the container of isopentane, chilled to  -70 °C on dry ice. Only samples without OCT should be submerged entirely, otherwise float the cryomold to immerse 5 of the 6 sides. It is possible to add powdered dry ice to the top of the cryomold.

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- 9 When completely frozen, wrap sample in foil (still in cryomold), and keep on dry ice. Alternatively, remove sample from cryomold and place in labelled cryotube on dry ice.

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- 10 Transfer samples on dry ice, to labelled box in  -80 °C freezer.

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- 11 If freezing more samples monitor the isopentane temperature and adjust as necessary.

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- 12 Wearing cold-resistant gloves or the wooden forceps, carefully remove the metal container of iso-pentane from the dry ice. The iso-pentane can allowed to warm to -10C, then decanted to a suitable bottle for re-use. The bottle must be labelled and kept in a fume hood/ downflow table, with the lid loosened until the iso-pentane returns to room temperature.

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