

Oct 21, 2019

Encasing LAMP mixture in paraffin

DOI

<https://dx.doi.org/10.17504/protocols.io.8jbhuin>

Niels Appelman¹

¹Wageningen University

iGEM Wageningen 2019



Niels Appelman

Wageningen University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8jbhuin>

Protocol Citation: Niels Appelman 2019. Encasing LAMP mixture in paraffin. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.8jbhuin>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: October 21, 2019

Last Modified: October 21, 2019

Protocol Integer ID: 28995

Keywords: lamp mixture in paraffin, primers in paraffin, encasing lamp mixture, paraffin, lamp mastermix, primer, lamp, protocol

Abstract

This protocol can be used to encase LAMP mastermix and primers in paraffin.

- 1 Prepare a 10X primermix according to the following scheme

	Prim er	Con cent ratio n
	FIP	16 μM
	BIP	16 μM
	F3	2 μM
	B3	2 μM
	LOO P F	8 μM
	LOO P B	8 μM

- 2 Add the following components to a micro centrifuge tube

	Component	Volume added
	Mastermix	5 ul
	Lamp primers (10X)	1 ul

- 3 Open the centrifuge tube and let dry at 37 C until only solid is left.

Note: The longer this step takes, the more activity the enzyme loses. It is recommended to dry in a storage area with a low humidity, as to speed up this step. To minimize the loss of activity, lyophilisation is recommended.

- 4 Add 4 mg of paraffin wax



- 5 Heat at 60 C for 5 minutes to let the wax melt and encapsulate the dried LAMP mixture
- 6 After the wax dries, the continued activity of the enzyme can be tested. To do this, add the following components

	Component	Volume added
	Template	x uL
	MilliQ	10 - x uL

- 7 Incubate at 65 C for 30 minutes
- 8 Record the color of all samples. Yellowing of a sample indicates a positive result.
- 9 In case of doubt regarding the result of a sample, run it on an agarose gel. The presence of many bands, often not individually recognizable, indicates a positive result.