



Mar 14, 2020

## End-preparation of amplicon pools

DOI

[dx.doi.org/10.17504/protocols.io.bdp3i5qn](https://dx.doi.org/10.17504/protocols.io.bdp3i5qn)

Josh Quick<sup>1</sup>

<sup>1</sup>University of Birmingham



Josh Quick

University of Birmingham

---

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.bdp3i5qn](https://dx.doi.org/10.17504/protocols.io.bdp3i5qn)

External link: <http://lab.loman.net/protocols/>

**Protocol Citation:** Josh Quick 2020. End-preparation of amplicon pools. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bdp3i5qn>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** March 14, 2020

**Last Modified:** March 14, 2020

**Protocol Integer ID:** 34267

## Abstract

This is a subprotocol for generating a library from a single amplicon pool


## Attachments



One-pot native barco...






64KB

## Safety warnings

 See SDS (Safety Data Sheet) for safety warnings and hazards.




- 1 Set up the following reaction for each sample:

Component	Volume
DNA amplicons (5ng/ul)	 10 µL
Nuclease-free water	 2.5 µL
Ultra II End Prep Reaction Buffer	 1.75 µL
Ultra II End Prep Enzyme Mix	 0.75 µL
<b>Total</b>	 15 µL

- 2 Incubate at room temperature for  00:05:00

Incubate at  65 °C for  00:05:00

Incubate on ice for  00:01:00